



Table of volume 18 · issue 1 Contents

M.A.R.C.H. Products

Products designed to address Massive Hemorrhage, Airway Management, Respiratory Support, Circulation/Vascular Access and Hypothermia Prevention

4

Combat Casualty Response Kits · CCRK®

Kits designed to meet the needs/challenges of our customers and to address the leading causes of potentially preventable death



DoD Bleeding Control Response Kits · D-BCRK

Individual kits, packs and stations designed to equip bystanders and first responders with the essential medical gear needed to help save lives



TAC-EVAC

The most effective and highest quality evacuation platforms and accessories for use in the tactical environment



Medical/Surgical Products

Medical and surgical equipment for use in austere environments, during prolonged field care situations and in forward operating field hospitals



Rescue Equipment

Includes an array of rescue products ranging from lightweight, versatile Pegasus® Ropes to tactical lights and accessories



Treatment Support

From personal protection products to splinting and immobilization products to help from the point of injury throughout the entire continuum of care



Training Products & Services

NAR Training Course Overview • Products designed to train users on TCCC protocols to reduce preventable death at the point of wounding

97

Product Index

Product offerings indexed by name, NSN number, and Item number for quick reference and information

111





MASSIVE HEMORRHAGE | Tourniquets

C-A-T[®] Combat Application Tourniquet - Gen 7

• Scientifically proven 100% effective in occluding blood flow in both upper and lower extremities by the USAISR and recommended by CoTCCC

• Affordable, lightweight, rugged design for true circumferential pressure featuring patented free-moving internal band, proprietary Red Tip Technology® elliptical tip, and new single routing buckle for simplified training with single protocol for all applications

 Reinforced windlass and optimized windlass securing system with writeable time stamp to document application time

SPECIFICATIONS

● (€ **⊕** Packaged: H 6.5 in. x W 2.4 in. x D 1.5 in. Device: L 37.5 in.

Weight: 2.7 oz



online video available

< scan to view

COLOR	STYLE	ITEM#	NSN#
Black	Tactical	30-0001	6515-01-521-7976
Orange	Rescue	30-0023	6515-01-612-7129
Blue	Trainer	30-0033	6910-01-560-2972

CAT TOURNIQUET

Abdominal Aortic and Junctional Tourniquet -STABILIZED • Easy to use and quick to apply - Lightest, most compact

- Stabilizes the pelvis in addition to stopping bleeding
- Large area of effective pressure for greater reliability
- Provides stable and complete occlusion of blood flow to the lower extremities, upper extremities, or pelvis
- Can be applied at junctional sites (inquinal, axillary) or mid-abdomen
- Demonstrated ability to remain in place and effective during patient movement
- Similar hemostatic efficacy to Zone 3 REBOA¹ ¹Rall, Jason M., et al. ^{*}Comparison of Zone 3 Resuscitative Endovascular Balloon Occlusion of the Aorta and the Abdominal Aortic and Junctional Tourniquet in a Model of Junctional Hemorrhage in Swine." Journal of Surgical Research, vol. 226, 2018, pp. 31–39
- Lower training burden compared to REBOA procedure

junctional tourniquet on the market

AAJT-S

ITEM#	NSN#
30-0248	6515-01-695-4457

SPECIFICATIONS

Packaged: L 9 in. x W 8 in. x D 2 in. Stored: L 7.5 in. x W 6.5 in. x D 2 in. Weight: 17 oz

-® Junctional **Emergency Treatment Tool**

- Immediately treat high extremity wounds in the thigh/groin area that is not amenable to TQ application
- Unilateral or bilateral treatment solution
- Simultaneous occlusion of blood flow to both lower limbs that will not impede respiration

JETT

ITEM#	NSN#
30-0088	6515-01-616-5841



online video available

< scan to view

SPECIFICATIONS

Packaged: H 5.5 in. x W 7 in. x D 3.75 in. Overall Device: L 56.75 in. Buttocks Pad: L 19 in. x W 6.5 in. Compression Pad: L 4.85 in. x W 2.5 in. x D 1.1 in. T-Handle: H 4.62 in. x W 4.04 in. x D 0.63 in. Weight: 1 lb 9.6 oz



JETT Pack: Item# 80-0355 & 80-0356 a nylon pouch designed specifically for carrying the JETT

TQ Accessories | MASSIVE HEMORRHAGE



RIGID C-A-T® TQ CASE with Shirt Shield

- Allows carriage of a Gen 7 C-A-T® Tourniquet on your belt while protecting your shirt from undue wear caused by friction from the windlass clip
- Blade-Tech TekLok Belt Attachment fits up to a 2.25 in. belt
- Can be mounted either vertically or horizontally
- Injection Molded Nylon Case stands up to rugged usage
- · Cross medical identifier patch included

RIGID CAT TO CASE

COLOR	ITEM#
Black	80-0944
Coyote	80-0945
Ranger Green	80-0946

SPECIFICATIONS

H 5.88 in. x W 2.5 in. x D 2 in. Weight: 5.36 oz

Explore other TQ carrier options available at NARescue.com





- PALS/MOLLE-style connectors provide versatile attachment capabilities (mounts horizontally or vertically)
- Includes easy-open pull elastic tab and tourniquet "TQ" ID patch
- Compatible with both Gen 6 & Gen 7 C-A-T® Tourniquets (sold separately)

SPECIFICATIONS



L 6.25 in. x W 1.75 in. x D 1.5 in. Weight: 2.17 oz

NYLON CAT HOLDER

NATO NSN*

COLOR	ITEM#	NSN#
Black	30-0057	
Coyote	30-0044	
MultiCam	30-0059	
Ranger Gr.	30-0208	
OD Green	30-0053	6530-01-605-9432*


Cummerbund C-A-T® HOLDER

- Allow attachment to tactical armor such as assault and plate carriers
- Utilizes 4 in. elastic coupled with hook and loop attachment to external cummerbund flap
- · Rapid, easy access with ambidextrous one handed TQ removal
- Compatible with both Gen 6 & Gen 7 C-A-T® Tourniquets (sold separately)

SPECIFICATIONS

L 5 in. x W 4 in. x D 0.25 in. Weight: 0.6 oz

CUMMERBUND CAT HOLDER



Elastic C-A-T® HOLDER

- PALS/MOLLE-style strap with tuck finish eliminating the need for a button snap
- Versatile attachment to body armor, gear or belt for rapid & easy access
- Easy-open quick release hood & extended elastic body for TQ release outwards under pressure
- Compatible with both Gen 6 & Gen 7 C-A-T® Tourniquets (sold separately)

FLASTIC TO HOLDER

LLASTIC I WITOLDER		
COLOR	ITEM #	
Black	30-0184	
Coyote	30-0185	
Multi-Cam	30-0186	
Ranger Green	30-0191	





MASSIVE HEMORRHAGE | Hemostatics

• Sterile, non-woven Z-fold gauze impregnated with kaolin for temporary external control of bleeding from penetrating trauma

- Effective, simple and safe to use on wounds not amenable to tourniquet application
- Flexible & pliable, contours to all wounds
- Does not generate heat
- Individually wrapped in easy-rip, military-grade pouch
- Features an embedded x-ray detectable strip

SPECIFICATIONS

Standard Gauze: L 12 ft x W 3 in.

Packaged: H 5 in. x W 5.5 in. x D 0.625 in. Weight: 0.8 oz XL Gauze: L 4 yd x W 4 in. (2 ply)

Packaged: H 6.5 in. x W 6.5 in. x D 0.75 in. Weight: 1.75 oz

COMBAT GAUZE

SIZE	ITEM#	NSN#
Standard	30-0039	6510-01-562-3325
XL	30-0070	6510-01-600-4705



COUNCION COMBAT GAUZE

CHITOGAUZE®

SPECIFICATIONS

C E STERLE

Gauze: L 12 ft x W 3 in. Packaged: H 6 in. x W 5 in. x D 0.65 in. Weight: 0.8 oz

CHITOGAUZE

ITEM#	NSN#
30-0049	6510-01-591-7740



CELOX" RAPID

C E STENLE

SPECIFICATIONS

Gauze: W 3 in. x 5 ft Packaged: H 5.25 in. x W 4.25 in. x D 0.5 in. Weight: 1.1 oz

CELOX RAPID

ITEM#	NSN#
30-0069	6510-01-623-9903

NUSTAT®

- Patented dual-fiber technology controls bleeding through absorption at the injury site, which results in a physical aggregation of blood cells and clotting factors
- Compact packaging allows for convenient storage and rapid application

NUSTAT

SPECIFICATIONS

Tactical: W 2 in. x L 36 in.

Pads: 4 in. x 4 in. (Pk of 2)

OPTION ITEM# 30-0236 4 in. x 4 in. Tactical 30-0170



Gauze & Bandages | MASSIVE HEMORRHAGE



CE R STERLE



Emergency Trauma Dressing

- Premium cotton construction for max quality and absorption
- High quality elastic wrap with a sterile non-adherent pad
- Durable C-clasp securing device

SPECIFICATIONS 4 in. ETD: L 70 in. x W 4 in.

Pad: H 6 in. x W 4 in

Weight: 2.6 oz

Abd/Stump ETD: L 70 in. x W 6 in.

Weight: 5oz Std Pkg: L 4.5 in. x W 2.5 in. x D 1.75 in. Pad: H 12 in. x W 12 in. Pkg: L 5.25 in. x W 6 in. x D 1.51 in.

Flat Pkg: L4.5 in. x W 3.5 in. x D 1.13 in. *6 in. ETD*: ∠ 70 in. x W 6 in. Pad: H 6 in. x W 6 in. Weight: 4 oz

Std Pkg: L 6 in. x W 2.5 in. x D 2.1 in. Flat Pkg: L6 in. x W 3.5 in. x D 1 in.

8x10 inch ETD:

L 70 in. x W 6 in. Weight: 4oz Pad: H 8 in. x W 10 in. Pkg: L 4.13 in. x W 7 in. x D 1 in.

EMERGENCY TRAUMA DRESSINGS

SIZE	STYLE	ITEM#	NSN#
4 inch	Standard	30-0031	6510-01-558-4108
4 inch	Flat	30-0093	
6 inch	Standard	30-0032	6510-01-558-4114
6 inch	Flat	30-0094	
Abd/Stump		30-0012	6510-01-541-8121
8x10 inch		30-0165	6510-01-654-6801

COMPRESSED GAUZE

- Low-cube, ultra-compact sterile gauze
- Premium rolled 100% cotton gauze for maximum absorption



S-ROLLED GAUZE

- Unique dispenser, allows controlled linear feed application
- Premium 100% cotton



S-ROLLED GAUZE

ITEM#	NSN#
30-0003	6510-01-605-2234

CRAVAT (Triangular Bandage)

Bandage: 37 in. x 37 in. x 52 in. Pkg: 3.5 in. x 2 in. x 0.63 in. Wgt: 1.6 oz

CRAVAT

ITEM#	NSN#
30-0089	6510-00-201-1755



TACTICAL COMPRESSION WRAP



Wrap: L 38 in. x W 5 in. Wgt: 6oz Pkg: 5.4 in. x 2.7 in. x 2 in.

TAC-WRAP

THE WILL	
ITEM#	NSN #
30-0016	6510-01-614-0938



AIRWAY MANAGEMENT

SPECIFICATIONS

Weight: 5.2 oz U.S. Patents: Pending

BAC-PACK ITEM#

10-0052

Packaged: H 8 in. x W 4 in. x D 1.75 in.

· All-in-one cricothyroidotomy kit for use in any environment

- Pre-inserted bougie in tube allows for a more rapid, easy procedure
- Tube meets TCCC recommendations for bougieaided open surgical technique
- #10 Scalpel included for safety and ease of use

1 #10

- Small, durable pouch keeps the items organized and ready to use
- · Lightweight and ruggedized
- Components configured by order of use
- Scalpel, tube, syringe all provided sterile
- Patent pending

INCLUDED IN THE BAC-PACK:



BAC-PACK[™] Bougie-Aided Cricothyroidotomy Pack


















NAR Tracheostomy Kit with

BOUGIE INTRODUCER



TACTICAL CRICKIT® Tactical Surgical Cricothyroidotomy Kit

- Contains all features required by CoTCCC for a surgical airway kit in the tactical setting
- Ideal length & size tracheal tube with securing device
- Configured in the sequence of use, ideal for high-stress, low-light, austere environments
- · Lightweight, ruggedized and minimal in cube size
- Included rolled pouch provides a clean working surface

INCLUDED IN THE TACTICAL CRICKIT:



















SPECIFICATIONS

Packaged: H 6.5 in. x W 2.5 in. x D 1.25 in. Unrolled: H 6.5 in. x W 7.5 in. x D 1.25 in. Weight: 3.7 oz

TACTICAL CRICKIT

ITEM#	NSN#
10-0017	6515-01-540-7568

NAR TRACHEAL HOOK



NAR TRACHEAL HOOK

ITEM#	NSN #	SPECIFICATIONS
10-0020	6515-01-562-9734	L 4.2 in. x W 0.5 in. x D 0.8 in. Weight: 0.5 oz

NAR TACTICAL SUCTION DEVICE

SPECIFICATIONS

L 9.75 in. x W 3.75 in. x D 3 in. | Wgt: 7.3 oz Suction Capacity: 1000 mL Vacuum Force: 100 mmHg

TACTICAL SUCTION	
ITEM#	NSN #
10-0018	6515-01-540-7206

NASOPHARYNGEAL AIRWAY

with lubricant

SPECIFICATIONS

L 6 in. x W 1.6 in. x D 0.5 in. Weight: 0.5 oz



22F-36F available at NARescue.com

TM Supraglottic **KING LTS-D** Airway Kit

 Unique dual channel design, anatomical shaped distal tip/cuff & open second lumen

€ (€

 Tapered ramp design provides additional ventilation outlet

SPECIFICATIONS

Folded: H 8 in. x W 11.43 in. x D 1.92 in. Unfolded: H 8 in. x W 13.36 in. x D 1.22 in. Weight: 5.7 oz

KING LTS-D

CLR	SIZE	ITEM#	NSN#
YEL	3	10-0025	
RED	4	10-0026	6515-01-591-1578
PUR	5	10-0027	6515-01-591-0397

INCLUDED WITH KING LTS-D™











RWAY MANAGEMENT

AMBU RDIC Mark III Resuscitator

- Genuine Ambu Mark III reusable resuscitator equipped with impermeable elastomeric components on the compression bag and patient valve (butyl rubber)
- Specialized inlet filtration system, including C2 canister, for filtering contaminated air

SPECIFICATIONS L 14 in. x W 7.5 in. x H 7 in. Weight: 5.4 lb

AMBU RDIC (
ITEM#	NSN#
10-0047	6515-01-338-6602



CYCLONE[®] *Improved* Pocket BVM

• Unique collapsible design packaged 25% smaller than conventional BVMs

• Requires minimum pressure for compression of the bag to ventilate

SPECS	CYCLON
J. 200	

Diameter: 5.3 in. x H 2.9 in. Weight: 17.4 oz

CYCLONE BVM	
ITEM#	NSN#
10-0019	6515-01-568-0193



Operations LARYNGOSCOPE SET

- Perform direct visual laryngoscopy even in the most severe environments
- Lightweight, compact kit to equip medical care providers in the field with tools to secure an airway

SPECIFICATIONS

Open: L 8 in. x W 5.5 in. x D 2 in. Rolled: L 8 in. x W 2 in. x D 2 in. Weight: 10 oz

• Contains 2 single-use fiber optic blades, textured handle & removable light source all rolled up in a rugged nylon pouch

LARYNGOSCOPE SET

ITEM#	NSN#
10-0001	6515-01-450-9790

INCLUDED IN THE 1 Pouch, LARYNGOSCOPE SET:













NAR Fiber Optic LARYNGOSCOPE BLADES

LARYNGOSCOPE BLADES (
ITEM#	NSN# available on all sizes offered	
See NARescue.com	See NARescue.com for details	
See NARescue.com	See NARescue.com for details	
	ITEM # See NARescue.com	

Nylon



HYFIN® Chest Seals

Twin Pack

- Two individually packaged, vented/valve-type chest seals in one unique package designed for treatment of both entry/exit wounds or multiple penetrating injuries.
- Patented 3 channel pressure relief vents
- Ruggedized dual perforated pouch designed to allow the provider to open only one dressing at a time, as needed
- Advanced adhesive technology designed to provide a superior seal



HYFIN CHEST SEAL TWIN PACK



ITEM#	NSN#
10-0037	6510-01-642-6212

SPECIFICATIONS

Chest Seals | RESPIRATORY SUPP

Device: 6 in. x 6 in. (excludes tab) Wgt: 2.5 oz Packaged: H 7.5 in. x W 4.5 in. x D 0.25 in. Unfolded: H 7.5 in. x W 9 in. x D 0.13 in.



online video available

< scan to view

HyFin® Vent Individual

HYFIN VENT INDIVIDUAL

SPECIFICATIONS

ITEM#	NSN#
10-0029	6510-01-624-0840



Device: 6 in. x 6 in. (excludes tab) Wgt: 1.31 oz Packaged: H 4.63 in. x W 7.5 in. x D 0.1 in.



HyFin® Occlusive Individual

HYFIN INDIVIDUAL

ITEM#	NSN#
10-0015	6510-01-532-8019



Device: 6 in. x 6 in. (excludes tab) Wgt: 1.31 oz Packaged: H 4.63 in. x W 7.5 in. x D 0.1 in.

USMC HyFin® Chest Seal **COMBO PACK**

- Versatile combo pack designed for the treatment of both entry/ exit or multiple penetrating injuries to the chest
- Includes one HyFin® Occlusive & one HyFin® Vented Chest Seal packaged in a clear resealable pouch

SPECIFICATIONS

Packaged: L 9 in. x W 5 in. x D 0.2 in. Each Chest Seal (when deployed): H 6 in. x W 6 in. x D 0.1 in. Total Weight: 2.9 oz US Patent 7,504,549 & Patent(s) Pending

USMC COMBO PACK:







ITEM #	NSN #
85-1045	6510-01-654-6787





RESPIRATORY SUPPORT

Tension Pneumothorax


10 ga x 3.75 inch

SPEAR[™] Simplified Pneumothorax Emergency Air Release

- Latest innovation in Needle Decompression technology
- Designed based on known gaps in the management of tension pneumothorax
- Enhanced needle tip for improved insertion control
- User friendly features like spin-lock connection, depth markings, and detachable low-pressure one-way valve
- Fenestrations along catheter for additional air flow

- Familiar, rugged, pen-style case for ease of storage and deployment
- · Allows for either lateral decompression or anterior needle thoracostomy
- Catheter made of flexible kinkresistant material with a soft tip
- Single use only. Supplied sterile





3.25 inch ARS® NEEDLE DECOMPRESSION KIT

 Penetrate the pleural space in 99% of casualties presenting the signs and symptoms of a tension pneumothorax

 Quality, reliable devices with 3.25 in, length needle/catheter

• Capless flash chamber provides immediate confirmation of correct placement

 Patented design packaged in an easy-to-open, crush-proof protective tube featuring easy-to-identify textured twist top with clip

SPECIFICATIONS

Pkg: H 5.7 in. x W 0.7 in. x D 0.8 in. Weight: 0.6 oz

ARS NEEDLE DECOMPRESSION KIT

GAUGE	ITEM#	NSN #
14 ga	ZZ-0056	6515-01-541-0635
10 ga	ZZ-0298	6515-01-673-1701

U.S. Design Patents: D584,409 S; D595,847S; D584,410S; U.S. Patent 7,874,426B2; @ North American Rescue® LLC., 127820, 127822, 127823, 001002372-0001, 001013940-0001, 001013940-0002





14 ga

online video

NAR SIMPLE THORACOSTOMY KIT



Packaged: L 4.5 in. x W 7.5 in. x D 0.75 in. Weight: 5.28 oz

- Compact kit that contains all necessary equipment to perform simple thoracostomy
- Includes 3-pack skin prep for performing the procedure bilaterally if needed
- HyFin Compact Vented Twin Pack Chest Seal to cover site if required

SIMPLE THORACOSTOMY KIT

ITEM#	UNIT	NSN#
83-0027	Single	













CHEST TUBE INSERTION KIT


 $\overline{\mathbb{X}}$

1 Valve

 Designed to provide the necessary equipment needed to perform emergency chest tube thoracostomy

 Compact and lightweight kit packaged in an easy to open, vacuum sealed pouch

• Can be stored in



SPECIFICATIONS

Packaged: L 10 in. x W 6 in. x D 2.25 in. Weight: 11.2 oz

CHEST TUBE INSERTION KIT

ITEM #	UNIT	NSN#
83-0007	Single	























NAR SALINE LOCK KIT Meets all requirements for US Army Task: 081-831-1011



Essential components for establishing a saline lock for safer, quicker and more effective IV cannulation in the most demanding environments

NAR SALINE LOCK KIT

ITEM#	NSN#
30-0013	6515-01-537-4094

SPECIFICATIONS

Pkg: H 5.5 in. x W 3.25 in. x D 1 in. Weight: 1.3 oz

SALINE LOCK KIT CONTENTS:



















NAR NEEDLELESS SALINE LOCK KIT

IV safety catheter & needle-free valve port reduce or eliminate accidental needle-stick exposures

NAR NEEDLELESS SALINE LOCK

ITEM#	NSN#
30-0046	6515-01-607-1943

SPECIFICATIONS

Pkg: H 5.5 in. x W 3.25 in. x D 1 in. Weight: 1.3 oz



NEEDLELESS SALINE LOCK KIT CONTENTS:





















BOA® IV CONSTRICTING BAND

- Direct peripheral IV access utilizing a revolutionary rolling action for transfusing venous blood to the distal veins while still allowing arterial flow
- Quick Click connector & two-finger Quick Release improves "time on task"
- XL version stretches approximately 50% longer than the standard BOA®

SPECIFICATIONS

Standard: L 7.5 in. to 18 in. | XL: L 10 in. to 27 in. Packaged: L 4 in. x W 3 in. x D 0.9 in. | Weight: 1.5 oz US Patents: 5,304,202 & 5,607,448




online video available

< scan to view







CE

SHARPS SHUTTLE

SPECIFICATIONS

L 6.38 in. x Diameter 1.11 in. Weight: 1.5 oz

SHARPS SHUTTLE

UNIT	ITEM #	NSN#
Single	30-0014	6515-01-518-8551
6 pack	30-00146	6515-01-619-5671

RAPTOR® IV Securing Device

SPECIFICATIONS

Extended: L 19 in. x W 1 in. Weight: 1.9 oz



RAPTOR IV SECURING DEVICE

UNIT	ITEM#	NSN#
6 pack	30-0006	6515-01-523-3317

G3 POWER DRIVER by EZ-10[®]

- Battery-powered solution for providing rapid vascular access to humerus & tibia in less than 10 seconds (Not for sternal use)
- Trigger guard prevents accidental activation
- Molded Hard Side Carrying Case (sold separately)

G3 POWER DRIVER

UNIT	ITEM #	
Driver	30-0075	
Case	30-0077	

SPECIFICATIONS = CE

H 4.5 in. x W 6.5 in. x D 2.5 in. Weight: 15 oz



MANUAL EXTREMITY NEEDLE SET



POWER DRIVER **NEEDLE SET**

for use with the EZ-IO® G3 Power Driver

POWER DRIVER NEEDLE SET

OWER DRIVER REEDEL SE		
COLOR	SIZE	ITEM#
Pink	15mm/PD	30-0081
Blue	25mm/AD	30-0083
Yellow	45mm/LD	30-0085

SPECS (E STERLE Packaged: H 4.75 in. W 4.25 in. D 0.7 in. Weight: 1 oz

POWER DRIVER NEEDLE SET:







T.A.L.O.N.TM MANUAL NEEDLE SET Manual I/O insertion at 7 extremity sites – including the sternum

SPECIFICATIONS



Device: L 6 in. x W 3.25 in. Weight 1.5 oz

T.A.L.O.N. MANUAL NEEDLE SET

SIZE	ITEM#
38.5mm	30-0132















LIFEWARMER[™] IV Fluid & Blood Warming System

• **Lightweight** The Quantum™ has the lightest operational weight of any portable blood & fluid warming system (Total Operational Weight: 22 oz)


• **Powerful** Warms two units of 39°F refrigerated whole blood at 100mL/min+ on a single charge

• Compact Fewest components to carry (three)

• Fast Warm-up time of only 24 seconds

• **Intuitive** The thermal transfusion/infusion set looks, handles, and sets up/primes like any conventional IV set

• **Safe** Tested safe for blood, blood products & IV solutions



• **Complete** Does not require additional standard transfusion or infusion set – all other warmers do

• **Tactical** Only system FDA cleared allowing easy disable of audible alerts to protect operators in combat

QUANTUM COMPLETE SYSTEM

ITEM#	NSN#
35-0006	6515-01-693-7810

OPERATIONAL CAPABILITIES

IV FLUIDS From a single battery charge, capable of warming 1700 mL of 68°F IV fluids to an output temperature of 100.4°F (38°+/- 2°C).

Warm up time: **24 seconds.**Delivery Flow rate: **200mL/min.**



1 x *Quantum* Charger

1 x Quantum Battery

1 x Quantum Controller





2 x TTS-B Thermal Transfusion Set - Blood

(Blood, blood products & IV solutions)

SPECIFICATIONS

Total Operational Weight: 22 oz (Quantum™ Controller, Battery & TTS-B)

Battery: 16 oz Controller: 3 oz Battery Charger: 22 oz Controller Cover: 0.23 oz

Thermal Transfusion Set (TTS-B): 2.7 oz

BLOOD or BLOOD PRODUCTS

From a *single battery charge*, capable of warming *2 units* of *39°F Whole Blood* to an output temperature of *100.4°F (38°+/- 2°C)*.

Warm up time: *24 seconds*.

Delivery Flow rate: *100mL/min*.

REV062121

TTS-B Thermal Transfusion Set

- Compatible with whole blood, component therapy (including FDP), colloids, and isotonics
- Similar size and weight to traditional IV tubing
- Intelligent sensing & heating elements integrated into the walls of the tubing set
- No disruption to the fluid path or laminar flow of the fluid
- Easy bi-directional connection to the Quantum Controller

SPECIFICATIONS

TTS-B package (each): 6.5 in. x 11 in. Tubing Length: 80 in. TTS-B Wat (each): 2.6 oz

Packaged Wgt (2 TTS-B & box): 11 oz



TTS-B	(2 Pack)	
-------	----------	--

ITEM#	NSN#
35-0005	6515-01-695-4444

TIS Thermal Infusion Set

- Similar size and weight to traditional IV tubing
- Intelligent sensing and heating elements integrated into the walls of the tubing set
- No disruption to the fluid path or laminar flow of the fluid
- Easy bi-directional connection to the Quantum Controller
- Compatible with standard IV fluids such as colloids and isotonics

SPECIFICATIONS

TIS packaged (each): 6.5 in. x 11 in. Tubing Length: 80 in.

Weight: TIS (each): 2.1 oz Packaged (2 TIS & box):

10 oz



ITEM#	NSN#
35-0004	



Quantum **BATTERY**

- Lightweight at only 16 oz
- Designed for the tactical environment- water resistant, subdued color, and the option to disable the audible alarm
- 44V Lithium-Polymer Battery
- Intelligent, state-of-the-art battery technology
- Status button to check charge state at any time

Quantum

HARDCASE

SPECS I 10.63 in. x W 9.63 in. x D 6.73 in. Weight: 3.22 lb

- Rugged hard case designed to protect your Quantum from dust, dirt, and moisture • Features pressure relief
- valve for automatic pressure equalization
- Custom laser cut high density foam holds components securely in place
- Flexible storage capacity for multiple batteries, controllers and thermal tubing sets

QUANTUM HARDCASE

ITEM#	NSN#
35-0008	

QUANTUM BATTERY

ITEM#	NSN#
35-0001	

SPECIFICATIONS

H 4.5 in. x W 3.25 in. x D 1.25 in. Weight: 16 oz

Quantum

CONTROLLER

- Moderates the energy flow from the Quantum Battery to the Thermal Tubing
- Designed for ease of connecting and durability
- Colored LED lights constantly show during administration of blood/fluid to show status
- Lightweight & durable
- IP53 rated for dust & water resistance with jacket

• Minimum service life -1000 insertions/removals

SPECIFICATIONS

Replacement cover

is also available

Controller Body: W 1.25 in. x H 3 in. x D 0.75 in. Controller Total Length: 51.25 in. Weight (Controller only): 2.5 oz Packaged Weight: 5 oz

QUANTUM CONTROLLER

OPTION	ITEM#	NSN#
Controller	35-0003	
Cover (Replacement)	35-0010	

Fully charge the Quantum Battery in < 90 minutes

- Smart, safe charging
- Includes 6 ft AC power cord
- Colored LED lights to inform the user of charging and battery status
- IP 22 rated for dust and water resistance
- 100V to 264V AC input; 50Hz, 60Hz
- Minimum service life-2500 insertions/removals

QUANTUM CHARGER

COMMITTED OF THE CONTROLLER	
ITEM#	NSN#
35-0002	

Quantum CHARGER

SPECIFICATIONS

Charger: H 6 in. x W 2.25 in. x D 1.75 in. AC power cord length: 6 feet Charger Wgt: 1 lb 4 oz Wgt with cord: 1 lb 10 oz



RCULATION & Vascular Access Field Blood Transfusion Kits

- Comprehensive kit designed for use with the Quantum Blood & Fluid Warming System allowing rapid, consistent normothermic field blood transfusion in far forward and near point of injury (POI) prehospital casualty management.
- Includes the Quantum Thermal Transfusion Set (Blood) that incorporates an integrated heating system that does not disrupt the fluid path
- Equipped with Donor/Recipient Modules to verify blood types of the donor and recipient
- Also allows collecting of blood samples from the donor for additional surveillance **QUANTUM FIELD** testing

SPECIFICATIONS

L 10 in. x W 7 in. x D 2.5 in. Weight: 1.09 lb

BLOOD TRANSFUSION KIT

ITEM#	NSN#
83-0028	

QUANTUM FIELD BLOOD TRANSFUSION KIT CONTENTS:









1 DD572

Questionnaire







MADE IN

RECIPIENT

DONO

RECIPIEN

*Each Donor Module Kit contains:



1 Permanent



2 Large Black

Nitrile Gloves




1 Beaded

Cord, 11 in.









DESCRIPTION OF REAL PROPERTY.

DONOR















System Clinical Practice Guidelines

**Each Recipient Module Kit contains:





















FBTK Field Blood Transfusion Kits

- Comprehensive kit that provides all the needed supplies to allow the user to rapidly perform a field blood transfusion
- Equipped to allow for verifying the blood types of the donor and recipient and collecting blood samples from the donor for additional testing

FBTK

OPTIONS	ITEM#
FBTK Basic	83-0022
Blood Surveillance	83-0025
FBTK Complete*	83-0000

Recipient Module:

5 x Alcohol Prep Pad

1 x IV Catheter - 18 ga x 1.25 in.

2 x 2-pack Gauze Pad - 4 in. x 4 in.

1 x IV Constricting Band - 1 in.x 18 in.

1 x Clear Adhesive Dressing - 4 in. x 4.75 in.

1 x IV Catheter - 16 ga x 1.25 in.

Additional Supplies:

1 x roll 1 in. Surgical Tape

2 x pair Black Nitrile Gloves, Large

1 x Black Permanent Marker

1 x Kelly Curved Hemostat - 5.5 in.

1 x Single Blood Bag (450ml)

1 x 550 Cord - 9.5 in.

1 x Beaded Security Tie - 11 in.

1 x IV Set Y-Type with Filter 170-260 UM

1x Questionnaire Card - DD572 Modified

Surveillance Module:

Donor Module:

1 x IV Hypodermic Needle - 18 ga x 1.5 in.

1 x IV Syringe, 10cc

1 x Eldon Card

1x Interlink Injection Site Needle

1x Blood Transfer Device

2 x Vacutainer, Gold Top (5ml)

1 x Vacutainer, Purple Top (5ml)

2 x Vacutainer, Pearl Top (5ml)

Recipient Module:

1 x 3 Way Stopcock

1 x Eldon Card

1 x Interlink Injection Site Needle

*FBTK Complete:

Contains all components of FBTK Basic & Blood Surveillance Module



LifeFlow® **FLUID INFUSER**

LifeFlow volume resuscitation. when minutes matter.

Air Check™

air enters the system

makes priming and setup easier

and will stop the flow of fluid if

- Hand-operated, rapid infusion device that allows health care providers to deliver fluids quickly and effectively
- Deliver 500ml of fluid in less than 2 minutes through a 20G IV catheter
- Designed to protect the syringe from bacterial contamination
- Can be operated with one hand

SPECIFICATIONS

12 in. x 7.5 in. x 2.75 in. Weight: 0.95 lb

LIFEFLOW INFUSER

ITEM#	NSN#
30-0286	

Clear Canopy Syringe is protected from contact & bacterial contamination during use

> Fast and Accurate 10ml of fluid with each pull

Reducer™ smooths the pulsatile flow, reducing the amount of effort required by the user

Force

Improve Resuscitation Workflow with LifeFlow®

PRESUSCITATE

500 ml in < 2 minutes Stav with Your Patient

Deliver first bolus quickly to reverse hypotension in patients with shock

- Provide active resuscitation
- Maintain control over delivery

PRESTORE

Improve Perfusion / Reverse Shock

Assess using vitals & bedside markers of perfusion, decide if additional fluid bolus is needed.

Early fluids more rapidly reverse shock & improve outcomes.1-3

>REASSESS

Rapid Decision / Disposition

Each bolus delivered provides an opportunity to evaluate and determine next steps.

Faster resuscitiation may lead to:

- Less mechanical ventilation⁴
- More rapid disposition
- Less donning/doffing of PPE

LifeFlow® PLUS

BLOOD & FLUID INFUSER

- Hand-operated, rapid infusion device that allows health care providers to deliver fluid & blood products quickly and effectively
- Deliver 500ml of fluid in less than 2 minutes through a 20G IV catheter
- Deliver 1 unit of whole blood in 2 minutes
- Goes from the field to the hospital
- Can be operated with one hand
- Handle and tubing packaged together for quick and easy set-up
- Designed to protect the syringe from bacterial contamination
- Average time to set up and prime: less than 3 minutes

SPECIFICATIONS

12 in. x 7.5 in. x 2.75 in. Weight: 0.95 lb

LIFEFLOW PLUS INFUSER

ITEM# NSN# 30-0287



MMINIM

A Leisman, D. B. Wie, M. Doerfler, A. Bianculli, M. F. Ward, M. Akerman, J. K. D'Angela and J. A. Zemmel D'Amore (2016). "Association of Fluid Resuscitation Initiation Within 30 Minutes of Severe Sepsis and Septic Shock Recognition With Reduced Mortality and Length of Stay." Ann Emerg Med 68(3): 298-311.

0.9% Sodium Chloride

bservational study." Eur merg Med 25(2): 97-104.

*Leisman D. et al. (2017).
Patterns and Outcomes
Associated With Timeliness of
Initial Crystalloid Resuscitation
in a Prospective Sepsis and
Septic Shick Cohort. Crit Care
Med. Oct. 45(10):1596-1606.

CE

HYPOTHERMIA Prevention & Management

HPMK® Hypothermia Prevention & Management Kit

• Includes reinforced Heat Reflective Shell (HRS™) of 4-ply, composite fabric with a protected non-conductive reflective layer

• Built-in hood, integrated fluid absorption pad & top-to-bottom taper providing excellent thermal insulation

 Continuous 1.5 in. Velcro® closures provide rapid 360° access to the casualty with minimal exposure to the elements

 Self-heating, 4-cell liner (Ready-Heat™ Blanket) provides 10 hours of continuous dry heat up to 125°F (52°C)

 Vacuum-packaged for a low-cube footprint

SPECIFICATIONS

Shell Pkg: H 6.75 in. x W 10.5 in. x D 5.5 in. Shell Open: L 78 in. x W 43 in. US Patent: 7,766,950 B2

Weight: 3 lb 8 oz

HPMK® CONTENTS:



HPMK

ITEM#	NSN #
80-0027	6515-01-532-8056

READY HEAT" BLANKET

NOTE: Not for sale in Canada

NOTE: Litter is

sold separately

HRS® Heat Reflective Shell

As an individual unit, reinforced 4-ply, composite shell with built-in hood, weight tested to support 250+ lb.

HEAT REFLECTIVE SHELL

ITEM#	NSN#
ZZ-0051	6532-01-627-7888

SPECIFICATIONS

Shell Packaged: L 10 in. x W 6.75 in. x D 2.5 in.

Open: L 78 in. x W 43 in.

US Patent: 7,766,950 B2

Weight: 1 lb 13 oz

the self-heating, 4-cell liner provides 10 hours of continuous dry heat after reaching peak temperature

reaching peak tempera approximately 15-20 minutes after opening

READY HEAT BLANKET

ITEM #

ZZ-0052

As an

individual unit,

SPECIFICATIONS

Open: L 48 in. x W 36 in. Weight: 1 lb 11 oz

EMERGENCY SURVIVAL WRAP

SURVIVAL WRAP

COLOR	ITEM#	
Orange/Silver	ZZ-0834	
OD Green/Silver	ZZ-1390	

Designed to reflect heat back to the body

Reversible Silver / high visibility Orange or OD Green

SPECIFICATIONS

CE

Packaged: L 5 in. x W 3.75 in.

Open Wrap: L 96 in. x W 60 in. | Weight: 3.7 oz











North American Rescue® is *THE market leader* for designing, producing, integrating and assembling tactical field medical sets and kits *supporting the services of the Department of Defense, federal agencies* as well as *state and local law enforcement.*

Over the past several years, our ability to provide rapid response products and services has aided in the delivery of critical battlefield medical care when and where it's needed most with an *end result of saving more lives*.



A comprehensive list of supply chain services affords a unique opportunity to North American Rescue's customers in the procurement, assembly, storage and management of their *complete medical* and material logistics requirements. NAR maintains a 75,000 square foot warehouse facility which consists of 50,000 square feet for receiving and distribution staging and an additional 14,000 square feet solely dedicated to the assembly room where sets, kits and custom medical solutions are built to order.

Highly qualified logistics and warehouse staff ensure that orders, inventory and distribution processes *maintain the superior quality standards* associated with the North American Rescue brand. Our integral partnership with the South Carolina Vocational Rehabilitation *NISH Partner* is a pillar in the successful execution of these processes.

SIRK Supplemental IFAK ReSupply Kits









SIRK Gen 2



- Vacuum-sealed first aid kits ideal for re-stocking an existing IFAK that contains the essential items to treat leading causes of preventable combat death
- Designed for small cube space in rugged, clear package allowing contents to be easily visualized



























SPECIFICATIONS H 7 in. x W 5 in. x D 2 in. Weight: 14 oz

SIRK GEN 2

ARS OPTION	ITEM#
w/o ARS Needle	85-1738
with ARS Needle	85-1899

SIRK Gen 1



1 Hyfin®Chest



2 Nitrile

Gloves, Lg



1 Surgical

Tape, 1 in.









SPECIFICATIONS H 7 in. x W 5 in. x D 2 in. Weight: 14 oz

SIRK GEN 1

HEMOSTATIC	ITEM#	
Combat Gauze	80-0291	
ChitoGauze	80-0225**	

IPOK Individual Patrol Officer Kit

- Individual hemorrhage control kit small enough to be carried on your person
- Compact, vacuum-sealed & packaged to fit into a BDU pocket, vest pouch or individual bag
- Rugged, durable pouch features NAR's signature Red Tip Technology® ensuring quality & performance









IPOK

HEMOSTATIC/GAUZE	ITEM#
S-Rolled Gauze	80-0167**
Combat Gauze, LE	80-0168



SPECIFICATIONS

H 4 in. x W 6.5 in. x D 2.75 in. Weight: 7.4 oz









- Contains the essential point-of-wounding equipment needed for treating penetrating, blast and other traumatic injuries
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access

M-FAK MINI FIRST AID KIT

COLOR	LEVEL	ITEM #
Black	Basic	80-0494
Coyote	Basic	80-0495
Multi-Cam	Basic	80-0496
OD Green	Basic	80-0497
Ranger Green	Basic	80-0996
Black	Combat Gauze LE	80-0498
Coyote	Combat Gauze LE	80-0499
Multi-Cam	Combat Gauze LE	80-0500
OD Green	Combat Gauze LE	80-0501
Ranger Green	Combat Gauze LE	80-0998
Black	Advanced	80-0502
Coyote	Advanced	80-0503
Multi-Cam	Advanced	80-0504
OD Green	Advanced	80-0505
Ranger Green	Advanced	80-0997























M-FAK ADVANCED:















ROO M-FAK Mini First Aid Kit



Rapid one handed access

to your tourniquet without ever having

to open the full kit!







- Contains the essential point-of-wounding equipment needed for treating penetrating, blast and other traumatic injuries
- Unique front pouch that easily opens with the pull of a tab for rapid access to your tourniquet
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access

SPECIFICATIONS

H 6 in. x W 3 in. x D 3 in. | Weight: 13 oz

ROO M-FAK MINI FIRST AID KIT

COLOR	LEVEL	ITEM#
Black	Basic	80-1049
Coyote	Basic	80-1050
OD Green	Basic	80-1051
Ranger Green	Basic	80-1052
Black	Combat Gauze LE	80-1053
Coyote	Combat Gauze LE	80-1054
OD Green	Combat Gauze LE	80-1055
Ranger Green	Combat Gauze LE	80-1056
Black	Advanced	80-1057
Coyote	Advanced	80-1058
OD Green	Advanced	80-1059
Ranger Green	Advanced	80-1060

Bag only option also available at NARescue.com













ROO M-FAK with **COMBAT GAUZE**:





























E-FAK Expeditionary First Aid Kit











• Rugged nylon platform designed to meet the Expeditionary Military Forces First Responder requirements for critical point-of-wounding medical supplies

- Includes an additional TQ Holder and Tourniquet for primary placement separate from the EFAK, meeting the (2) Tourniquet Expeditionary Military requirement
- Allows attachment both vertically and horizontally
- Advanced Treatment pack out

SPECIFICATIONS

L 6 in. x W 3 in. x D 3 in. Weight: 1 lb 3 oz

E-FAK

COLOR	ITEM#	
Coyote	85-2533	
Black	80-1044	











2 C-A-T®

TQ Holder



1 Combat

Gauze Z-fold



KIT CAPABILITIES INCLUDE


S-Rolled



1 ETD[™]4 in.

Dressing





• Equipped with a 3 in. belt loop attachment to accommodate a duty/pistol belt

• Constructed with rip stop nylon, elastic loop organizational storage and easy open zipper

• Purposely designed to look like a non-tactical, accessory bag



SPECIFICATIONS L 7 in. x H 2.25 in.

x D 3.25 in. Weight: 13 oz



PLUS+

BASIC















HyFin®





Nitrile

Gloves, Lg



COLOR	LEVEL	ITEM#
Black/Gray	Basic	80-0374
Blue/Gray	Basic	80-0375
Red/Gray	Basic	80-0915
Black/Gray	Combat Gauze	80-0921
Blue/Gray	Combat Gauze	80-0922
Red/Gray	Combat Gauze	80-0917
Black/Gray	Advanced	80-0372
Blue/Gray	Advanced	80-0373
Red/Gray	Advanced	80-0916

See OCHO Accessory Bag only on page 56



- Can connect to the assault vest/ body armor or to a duty/pistol belt
- Contains the essential equipment needed for self-aid/buddy-aid to address the leading causes of potentially preventable death in a tactical environment
- Quick Pull Tab easily opens the kit with just one hand
- Unique, detachable shock-cord inner panel that can be adjusted to your components
- Inner panel is removable for easy access to your kit

SPECIFICATIONS

H 8 in. x W 4 in. x D 3 in. Weight: 1 lb 2 oz

EAGLE IFAK

COLOR	LEVEL	HEMOSTATIC	ITEM#
Black	Basic	None	85-0177
Coyote	Basic	None	85-0178
Multi-Cam	Basic	None	85-0179
OD Green	Basic	None	85-0180
Black	Advanced	Combat Gauze	85-0189
Coyote	Advanced	Combat Gauze	85-0190
Multi-Cam	Advanced	Combat Gauze	85-0191
OD Green	Advanced	Combat Gauze	85-0192
4 options online	Advanced	ChitoGauze	See NARescue.com

EAGLE IFAK BASIC:

























1 Nitrile


Gloves, Lg



KIT CAPABILITIES INCLUDE:

COMPA





- Protective hood with quick pull latch release helps protect IFAK components and includes sewn in embroidered medical cross identifier
- Detachable (hook & loop backed) shock-cord inner panel adjusts to secure internal components and is easily removable for rapid access

RIG EAGLE IFAK

COLOR	LEVEL	ITEM#
Black	Basic	80-0509
Coyote	Basic	80-0510
Multi-Cam	Basic	80-0511
Ranger Green	Basic	80-0991
Black	Combat Gauze LE	80-0517
Coyote	Combat Gauze LE	80-0518
Multi-Cam	Combat Gauze LE	80-0519
Ranger Green	Combat Gauze LE	80-0993
Black	Advanced	80-0513
Coyote	Advanced	80-0514
Multi-Cam	Advanced	80-0515
Ranger Green	Advanced	80-0992

Bag only option also available at NARescue.com











RIG SERIES EAGLE with COMBAT GAUZE LE:















KIT CAPABILITIES INCLUDE: •

SPECIFICATIONS H 8 in. x W 4 in. x D 3 in. Kit Packed Weight: 1.18 lb

Empty Bag Weight: 7 oz





















Tactical Operator Response Kit T.O.R.K.



• Streamlined, compact IFAK designed to MOLLE to the assault vest or body armor

• Includes essential life-saving equipment to be used at the point of wounding

• For self-aid/buddy-aid injuries caused by penetrating or blast trauma

SPECIFICATIONS

H 7 in. x W 4.5 in. x D 2.75 in. Weight: 1 lb 5 oz

Bag only options also available at NARescue.com

T.O.R.K. TACTICAL OPERATOR RESPONSE KIT

COLOR	LEVEL/HEMOSTATIC	ITEM #
Black	Basic	85-0008
Coyote	Basic	85-0009
Digital Camo	Basic	85-2300
Multi-Cam	Basic	85-0012
OD Green	Basic	85-0054
Red	Basic	85-1805
Ranger Green	Basic	85-2301
Black	Basic w/Combat Gauze	80-1097
Coyote	Basic w/Combat Gauze	80-1098
Digital Camo	Basic w/Combat Gauze	80-1099
Multi-Cam	Basic w/Combat Gauze	80-1100
OD Green	Basic w/Combat Gauze	80-1101
Ranger Green	Basic w/Combat Gauze	80-1102
Red	Basic w/Combat Gauze	80-1108
Black	Advanced w/Combat Gauze	85-0048
Coyote	Advanced w/Combat Gauze	85-0049
Digital Camo	Advanced w/Combat Gauze	85-0050
Multi-Cam	Advanced w/Combat Gauze	85-0052
OD Green	Advanced w/Combat Gauze	85-0053
Ranger Green	Advanced w/Combat Gauze	85-2302
Red	Advanced w/Combat Gauze	85-1806
7 options online	Advanced w/ChitoGauze	NARescue.com



T.O.R.K. BASIC / BASIC with Combat Gauze:







































KIT CAPABILITIES INCLUDE: -

KIT CAPABILITIES INCLUDE:

USCG *Mini* IFAK



• Featuring our modular system that can be worn by the operator, hooked directly to the duty belt, or attached to an assault vest or body armor with MOLLE-style connectors

SPECIFICATIONS

H 6.5 in. x W 3.5 in. x D 3.5 in. | Weight: 1 lb





















Low Profile Basic IFAK designed for the US Coast Guard Boarding Teams

Versatile modular system can be worn by the operator as a drop-down leg rig, duty-belt addon or assault vest/body armor attachment (via MOLLE-style connectors)

USCG IFAK BASIC

ITEM#	NSN#
80-1023	6545-01-682-2596



Bag only option also available at NARescue.com



















Versatile modular system can be worn by the operator as a drop-down leg rig, duty-belt add-on or assault vest/body armor attachment (via MOLLE-style connectors)

USCG IFAK

COLOR	ITEM #	NSN#
Coyote	80-0158	6545-01-587-8513
Black	80-0149	6545-01-572-7248
Resupply*	85-0104	



2 Gauze,



1 HyFin®



2 Nitrile







Eye Shield





Bag only option also available at NARescue.com *Resupply does not include the Nylon IFAK Bag

KIT CAPABILITIES INCLUDE:

Meets specific requirements of the Navy Expeditionary Combat Command

Designed to meet the unique & specific requirements of the Navy Expeditionary Combat Command for a small, efficient IFAK

SPECIFICATIONS

H 7.4 in. x W 6.9 in. x D 5 in Weight: 1 lb 13 oz

MARITIME ASSAULT KIT

COLOR	ITEM #	NSN #
Multi-Cam	80-0245	
Coyote	80-0015	6545-01-529-4187
OD Green	80-0014	6545-01-531-6626
Black	80-0013	6545-01-529-4185
Digital Camo	80-0058	
































DIVIDUAL KIT

Designed with all the tools for the Warfighter/Operator to mitigate extremity bleeding, tension pneumothorax & airway obstruction

SPECIFICATIONS

H 8.5 in. x W 8.75 in. x D 4.5 in. Weight: 1 lb 12 oz

INDIVIDUAL KIT

COLOR	ITEM#	NSN#
Coyote	80-0005	
OD Green	80-0004	
Black	80-0003	6530-01-529-1134
Digital Camo	80-0091	

Bag only option also available at NARescue.com

















KIT CAPABILITIES INCLUDE:



KIT CAPABILITIES INCLUDE:





Designed for Tactical Healthcare Professionals configurable based on individual skill set & mission requirement

SPECIFICATIONS

H 9.5 in. x W 9 in. x D 7.75 in. Weight: 3 lb 9 oz

Bag only option also available at NARescue.com

MEDIC LEG RIG

COLOR	ITEM#	NSN#
Multi-Cam	80-0249	
Coyote	80-0048	6545-01-543-3170
OD Green	80-0047	6545-01-543-3174
Black	80-0046	6545-01-543-3185
Digital Camo	80-0087	



2 Saline

Lock Kit













__ Leash"

























- Contoured profile, large central compartment, multiple attachment points & front trauma-shear storage
- Multiple-load carriage capabilities (fanny pack, over the shoulder, attachments for 3-day assault packs & MOLLE-style webbing

Bag only option also available at NARescue.com































- KIT CAPABILITIES INCLUDE:



Digital Camo

SPECIFICATIONS H 7.75 in. x W 13 in. x D 9 in. Weight: 4 lb 6.5 oz

SQUAD KIT

Multi-Cam

OD Green Black

COLOR

Coyote



ITEM#

80-0251 80-0039

80-0038

80-0037

80-0093







NSN#

6545-01-543-2956

6545-01-543-3159

6545-01-543-5721


6545-01-543-5733

MARINER KIT

ITEM# NSN# 80-0107 6545-01-543-5770

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in. Weight: 4 lb 6.5 oz



MARINE

 Rugged rescue orange nylon bag for high visibility, easy transport & convenient storage

Designed to meet the rigorous requirements of the

Emergency Med Assist Team (EMAT) Response Bag

- Compact, contoured profile with large central compartment, multiple attachment points & front trauma-shear storage
- Designed to be worn as a belly bag or slung over the shoulder

































Bag only option also available at NARescue.com

HEST POUCH Rapid access to critical medical equipment from your chest

KIT CAPABILITIES INCLUDE:

















NAR-4 CHEST POUCH

COLOR	HEMOSTATIC	ITEM#
Multi-Cam	Combat Gauze	80-0247
Coyote	Combat Gauze	80-0175
OD Green	Combat Gauze	80-0289
Black	Combat Gauze	80-0173
Digital Camo	Combat Gauze	80-0176
5 options	Chito Gauze	Online

 Auxiliary NAR-4 companion kit streamlined to provide rapid access to critical basic & advanced equipment from the chest

• Designed to be secured to the front of the body armor using a MOLLE-style attachment system

 Ambidextrous accessibility to trauma shears











SPECIFICATIONS





ARS®

Needle, 14g











Bag only option also available at NARescue.com

KIT CAPABILITIES INCLUDE:









Configurable for individual skill set and mission requirement





- Low-profile, load carriage for Tactical Healthcare Professionals
- Modular system with 5 segregated removable equipment pouches to carry equipment consistent with **TCCC Guidelines**
- Designed to attach directly to the body armor with four rapid attachment straps or carried using shoulder straps

NAR-4 AID KIT

COLOR	HEMOSTATIC	ITEM#
Black	Combat Gauze	80-0181
Coyote	Combat Gauze	80-0183
Digital Camo	Combat Gauze	80-0184
Multi-Cam	Combat Gauze	80-0244
OD Green	Combat Gauze	80-0308
5 options	ChitoGauze	Online

SPECIFICATIONS

H 19 in. x W 15 in. x D 7 in. | Weight: 12 lb 8.3 oz































































2 HyFin®



2 BurnTec

4"x4" Dressing



2 SAM®

Splint II





SPECIFICATIONS H 14 in. x W 11 in. x D 5 in. (Depth expandable to 7.25 In.)

MINI MEDIC KIT

I-III I-IIII I-III				
COLOR	LEVEL	ITEM#		
Black	Basic	80-0979		
Coyote	Basic	80-0980		
OD Green	Basic	80-0981		
Ranger Green	Basic	80-1003		
Black	Combat Gauze	80-0985		
Coyote	Combat Gauze	80-0986		
OD Green	Combat Gauze	80-0987		
Ranger Green	Combat Gauze	80-1005		
Black	Advanced	80-0982		
Coyote	Advanced	80-0983		
OD Green	Advanced	80-0984		
Ranger Green	Advanced	80-1004		



NAR













COMBAT GAUZE also includes:













- Slim profile that can be worn with single bandolier style sling strap and with a detachable waist strap
- Two tear away internal panels designed to hold standard IFAK contents
- Two additional tear away clear accessory pouches for additional equipment
- Internal shock cord area to hold compact, soft litter or additional equipment
- External MOLLE on sides of bag for external tourniquet pouches
- Large panel on exterior of kit for ID patches

RUF KIT

LEVEL	ITEM#
Basic	80-0974
Combat Gauze LE	80-0975
Advanced	80-0976

Every RUF KIT includes these core items:















Weight: 6.3 lb





















ADVANCED also includes:







Can be worn with single bandolier style sling strap as well as with a detachable waist strap

MedicTRAUMA PACK



- 9 removable, color-coded pull bags and 2 removable side compartments provide modular organization
- Over 1,000 cubic inches of gear capacity and MOLLE-style attachment system compatible with removable side pouches

MEDIC TRAUMA PACK

**NATO NSN

COLOR	ITEM#	NSN#
Black	80-0070	6545-01-572-7251
Coyote	80-0072	6545-01-588-1862
Digital Camo	80-0096	
OD Green	80-0071	6545-01-594-5617**

SPECIFICATIONS

H 19 in. \times W 19 in. \times D 10 in Weight: 18 lb 1 oz

MEDIC TRAUMA PACK:















2 28F Naso



King

LT-D Size 4



1 Cyclone

Pocket BVM



1 Tactical

Suction Device



























KIT CAPABILITIES INCLUDE:

















- Official vehicle life-saving kit of the U.S. Army
- For treatment of multiple trauma casualties and the immediate evacuation of one patient
- Includes full-size Talon® 90C quad-fold litter
- Pocket system holds self or buddy-aid medical supplies
- Back pack straps provide easy transport from the vehicle to the point of need
- Includes Hypothermia Prevention & Management Kit[™] to help stabilize patients during transport
- Interfaces with any standard or non-standard TACEVAC asset

WARRIOR AID & LITTER KIT

COLOR	ITEM #	NSN #
Black	80-0097	
Coyote	80-0089	6545-01-587-1199
Digital Camo	80-0059	6545-01-532-4962
OD Green	80-0060	

SPECIFICATIONS

H 22.75 in. x W 16.5 in. x D 11.5 in. Weight: 29 lb 15 oz

WARRIOR AID & LITTER KIT: \(\text{\tint{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\tint{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\tint{\text{\text{\text{\tint{\text{\tint{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\text{\tint{\text{\tint{\text{\tint{\text{\tint{\text{\text{\text{\text{\text{\tint{\tint{\text{\tint{\tint{\tint{\tint{\text{\tint{\text{\tint{\tint{\tint{\tint{\tint{\tint{\tint{\text{\tint{\text{\tinit}}\tint{\text{\tinithtet{\text{\tinit{\text{\tinit{\text{\tinit{\text{\text{\text{\text{\text{\text{\tinit{\tinit{\tinit{\text{\tinithtet{\tinit{\text{\tinit{\text{\tinit{\tinit{\text{\tinit{\tex{\tinit{\tinithtet{\text{\tinit{\text{\tinithtet{\tinit{\text{\ti}\tinit{\text{\tinit{\tinit{\tinit{\texit{\tinit{\tinit{\tinit{\tinit{\tinit{\tinit{\tinit{\tinit{\tert{\tinit{\tinit{\tiit{\tinit{\tiitil\tiint{\tinit{\tinit{\tinit{\tinit{\tinit{\tiitht{\tin











































Contains WALK components for resupply - excluding nylon bag & litter

Designed to stow away perfectly in the Naval First Aid Box



- First responder rapid access and organization of critical life saving equipment for rendering on scene emergency care at sea
- Designed to meet the rigorous requirements of the U. S. Naval Ships, Military Sealift Command Ships and afloat commercial vessels
- Lightweight, durable slimline pouch is securable using anti-pilferage seals and features a clear window to readily view contents
- Pouch easily fits into the standard Naval First Aid Box (not included) and is ideal for restocking the Naval Gun Response Aid Kit as needed

SPECIFICATIONS

H 15 in. x W 6 in. x D 4 in Weight: 1 lb 15 oz

NAVAL FIRST AID BOX RESPONSE KIT

OPTION	ITEM#	NSN#
Kit	80-0419	6545-01-648-2705
Pouch only	80-0357	8465-01-666-7328

















Easily fits into the standard Naval First Aid Box (not included) and is ideal for restocking the Naval Gun Response Aid Kit (see below) as needed

Bag only option also available at NARescue.com











AVAL GUN Response Aid Kit



















esigned for the rigorous requirements of the US Navy



Weight: 1 lb 15 oz

See Above for Contents of the Naval First Aid Box Response Kit

- Designed for the rigorous requirements of the U.S. Navy for first responder emergency care equipment on board U.S. Naval Ships, Military Sealift Command Ships & afloat commercial vessels
- Ergo dynamic, low profile casualty response system features 2 removable equipment pouches that attach via hook-&-loop to a lightweight, durable outer bag
- Two internal pouches are identical & provide rapid access & organization of critical life saving equipment needed to render on-scene care at sea
- Both outer bag & 2 inner pouches are securable using anti-pilferage seals
- Slim-line, lightweight kit stores easily on a wall and can be carried over the shoulder or bandolier-style



NAVAL GUN RESPONSE AID KIT

OPTION	ITEM#	NSN #
Kit	80-0418	6545-01-648-2690
Bags only	80-0417	8465-01-666-7169

Bag only option also available at NARescue.com

Specifically designed for U.S. Navy Search and Rescue Swimmers





KIT CAPABILITIES INCLUDE: -

- The SAR Level A Medical Bag designed for U.S. Navy Search and Rescue (SAR) Swimmers
- Low profile & durable response bag featuring 6 internal pouches that attach via hook-and-loop allowing access to the specific pouch needed to provide on-scene medical care to the patient
- Various sized modules can be self labeled for the specific mission set; 4 smaller pouches for hemorrhage control products and 2 larger ones to accommodate airway & intravenous products
- Kit can be carried as a backpack using shoulder straps that stow away until needed

NAR 5 SAR KIT		NATO NSN*
OPTION	ITEM#	NSN #
Kit	80-0416	
Bags only	80-0399	8465-01-669-5796*

SPEC

antichoodod	OPTION	ITEM#	NSN#
SPECIFICATIONS H 24 in. x W 15 in. x D 7 in.	Kit	80-0416	
Weight: 14 lb 10 oz	Bags only	80-0399	8465-01-669-5796*
NAR-5 AID KIT:			





















4 Triangular Bandage















Naval USCG BO



For additional information about these products designed for emergency medical care on military & commercial afloat vessels, contact North American Rescue or visit www.NARescue.com

BOAT RESPONSE AID KITS

OPTION	ITEM#	NSN#
Naval Boat Response Aid Kit	80-0420	6545-01-648-3228
Naval Boat Response Aid Box Only	80-0354	8115-01-656-8316
USCG Boat Response Aid Kit	80-0353	6545-01-646-2623



Bag only option also

available at NARescue.com

NAVAL BOAT RESPONSE AID KIT:































USMC CLS KIT U.S. Marine Corps Combat Life Saver Kit



The Official Combat Lifesaver Kit of the U.S. Marine Corps

- Constructed of durable 1000 denier Cordura® solution-dyed nylon
- Sand and water resistant zippers with silent pulls & ITW Nexus Hardware

SPECIFICATIONS

Packaged: L12 in. xW10 in. xD8 in. Open Length: 37.5 in.

USMC CLS ReSupply (#80-0424)

Contains items for resupply in a clear resealable

plastic bag - excludes nylon bag

USMC CLS KIT

ITEM#	NSN #
80-0159	6545-01-571-4470

Weight: 6 lb 14 oz

USMC CLS KIT: \(\text{\text{\text{\text{\text{USMC}}}}\)



















Bag only #ZZ-0252 also available at NARescue.com



1 22F



1 Guedel



1 HyFin®

















JS Army CLS RESUPPLY KIT











US ARMY CLS RESUPPLY KIT

HEMOSTATIC	ITEM#
Combat Gauze	85-0197
ChitoGauze	85-0196

 Convenient kit containing the essential items to resupply the US Army CLS Combat Life Saver Kit delivered in a durable and easy-to-open resealable plastic bag

• Availble with 2 Hemostatic options: Combat Gauze or ChitoGauze

USMC CLS RESUPPLY:







































KIT CAPABILITIES INCLUDE: SQUAD RESPONDER K

• The NAR alternative to the US Army Combat Lifesaver (CLS) Kit

· Equipped with optimal medical gear to treat 2 to 4 casualties

• Constructed of Durable 1000 Denier Cordura® Solution-Dyed Nylon

 Large central compartment and multiple attachment points

• Sand and water resistant zippers with silent pulls

ITW Nexus Hardware

SQUAD RESPONDER KIT

ITEM #	NSN#
85-1764	

Bag only #80-0932 also available at NARescue.com















H 14 in. x W 14 in. x D 5 in Weight: 4 lb 14 oz



1 Surgical

Tape, 2 in.



4 Fox Eye

Shieldw/garter



4 Nitrile

Gloves, Lg





1 Trauma

Shears, Lg



3 Dry Sterile

Burn Dressing



1 Elastic

Bandage, 6 in.











SQUAD RESPONDER RESUPPLY

6 Alcohol

Prep Pads





US ARMY CLS RESUPPLY

ITEM#	NSN#
80-2495	

US ARMY CLS RESUPPLY:







NAR's alternative to the US Army Combat Lifesaver (CLS) Kit





















































Constructed of 500D Cordura® solution dyed nylon

• Sand- and water-resistant, self repairing zippers with silent 550 pull cords

• Separate lid compartment for the R.E.E.L. Non-Rigid Litter

• Six segregated removable pouches secured by hook and loop:

> 4 x Internal Bleed Pouches with Red pull handle

> 1 x Airway Pouch with Blue pull handle

>1 x Accessories Pouch with Coyote pull handle

• 2 internal lid pouches to secure additional medical products

• 2 external side pouches on either side of the kit for additional medical products

• 2 external straps with dual adjust MILSPEC IR reduced buckles to secure the kit in a vehicle or boat

SPECIFICATIONS

Case closed: H 17 in. x W 16.5 in. x D 10 in. Weight: 18 lb

ECRK

ITEM#	NSN#
80-0889	

Bag only **#80-0888** also available at NARescue.com



Lightweight platform to rapidly access and organize all your essential trauma management equipment

















































- Essential equipment for medical and non-medical personnel to address injuries common to shooting and training ranges
- Treatment capabilities include penetrating injuries, chest wounds, burns, eye injuries, minor lacerations, sprains and fractures

RANGE TRAUMA KIT

VΑ			

HEMOSTATIC	COLOR	ITEM#	NSN#
Combat Gauze	Orange	80-0299	
Chito Gauze	Orange	80-0213	6545-01-610-5411*

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in Weight: 5 lb 10 oz

















Bag only option also available at NARescue.com

KIT CAPABILITIES INCLUDE:























- Water-tight case maintains the integrity of your medical products when you need them most
- Kit contents arranged in a modular system featuring rugged lid storage organizer and removable nylon pouches, in the main compartment, with easy-to-read treatment labels

RANGE TRAUMA AID KIT

COLOR	ITEM #	
Orange	85-0889	

SPECIFICATIONS

H 9.5 in. x W 13.5 in. x D 6 in. Weight: 8 lb 12 oz











































K-9 HANDLER ROO M-



- One of NAR's most compact, versatile Mini First Aid Kits featuring rapid one handed access to your tourniquet from an easy-toopen front pouch
- Contains the First Responder's most requested critical point-of-wounding medical equipment for treating injuries typical to K-9s and their handlers operating in tactical or patrol environments
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Multiple elastic loops for secure gear organization
- Includes 1 in. Mini PVC Med Cross ID Patch

K-9 HANDLER ROO M-FAK

COLOR	LEVEL	ITEM#
Black	Basic	80-1067
Coyote	Basic	80-1068
OD Green	Basic	80-1069
Ranger Green	Basic	80-1070
Black	Combat Gauze LE	80-1071
Coyote	Combat Gauze LE	80-1072
OD Green	Combat Gauze LE	80-1073
Ranger Green	Combat Gauze LE	80-1074

Bag only option also available at NARescue.com

Rapid one handed access to your tourniquet without ever having to open the full kit!









1 C-A-T¹

Tourniquet



1 ETD[™] 4 in.

Mini Responde

SPECIFICATIONS L 7 in. x W 3.75 in. x D 4 in.

K-9 HANDLER ROO M-FAK:

Weight: 13 oz



Gauze

S-Rolled



1 Combat

Gauze® LE

KIT CAPABILITIES INCLUDE: .---

VAR K-9 LITTER & Carrying Strap Accessories





- K-9 Litter can be carried in cargo pocket or pouch
- Leg cut outs to safety carry the canine in the litter
- Six large and four small carry handles that allow for one to four person carries
- Optional K-9 Carry Strap & 2 Carabiners for over the shoulder carry of injured K-9 also available

K-9 LITTER

OPTIONS	ITEM#
K-9 Litter Only	60-0074W
Carry Strap & Carabiners	ZZ-1004
K-9 Litter with Carry Strap & Carabiners	60-0076

-9 HANDLER IFAK Designed especially for K-9 and Handler











8 Alcohol

Pad



















Bag only option also available at NARescue.com

• Designed for immediate lifesaving treatment to a wounded handler, operator or K-9

• Treatment capabilities range from major trauma typical to operating in tactical or patrol environments to K-9 emergencies

SPECIFICATIONS

H 7 in. x W 4 in. x D 2.75 in | Weight: 1 lb 3 oz

K-9 HANDLER IFAK

COLOR	HEMOSTATIC	ITEM#
Black		80-0439
Coyote		80-0440
Multi-Cam		80-0441
OD Green		80-0442
Black	Combat Gauze	80-0443
Coyote	Combat Gauze	80-0444
Multi-Cam	Combat Gauze	80-0445
OD Green	Combat Gauze	80-0446











• Designed for the handler or medic to deliver lifesaving treatment to a wounded K-9 at the point of injury

COMBAT GAUZE

- K-9 treatment capabilities from major trauma typical to operating in tactical or patrol environments to more common K-9 emergencies
- Provides the essential equipment to render care to the injured K9 when veterinary care is not immediately available

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in. Weight: 5 lb 2 oz

K-9 TACTICAL FIELD KIT

COLOR	HEMOSTATIC	ITEM #	
Black	Combat Gauze	80-0300	
Coyote	Combat Gauze	80-0301	
OD Green Combat Gauze		80-0303	
3 options online	ChitoGauze	See NARescue.com	

Bag only option also available at NARescue.com









































Eye



K-9 MINI MEDIC KIT Versatile, Multi-purpose Platform for K9 Trauma in the field

KIT CAPABILITIES INCLUDE: •

Inside Back

stowing extra equipment

Compartment for

- Compact, multi-purpose trauma aid bag with life-saving trauma equipment and K9 specific gear for point of wounding treatment of multiple traumatic injuries
- K9 specific gear enables handlers to render care to an injured canine when veterinary care is not immediately available
- Versatile platform features two external pockets with two external tourniquet holders, one external pouch with additional drop-down pouch for extra storage
- External MOLLE on sides and bottom of bag, internal MOLLE rows on inside back of bag, three clear tear away pouches | two internal mesh pouches on the inside main compartment, and one larger internal pouch on inner lid
- Internal tape roll strap also included

K-9 MINI MEDIC KIT:















4 in. Flat



2 28F Naso

Airway & Lube



2 HyFin®

Vent. Twin Pk























COLOR	ITEM#
Black	80-1075
Coyote	80-1076
OD Green	80-1077
Ranger Green	80-1078

SPECIFICATIONS H 14 in. x W 11 in. x D 5 in. (Depth expandable to 7.25 In.)

Weight: 12.5 lb

Bag only option also available at NARescue.com







KIT CAPABILITIES INCLUDE:

- Slim profile trauma aid pack designed for the K9 Handler or Tracking Team member
- Carries a minimum of two complete sets of IFAK contents along with K9 medical supplies
- Carriage options include single bandolier style sling strap as well as detachable waist strap
- Two tear away internal panels hold standard IFAK contents
- Two additional tear away clear accessory pouches to store additional equipment
- Internal shock cord area to hold K9 Litter
- External MOLLE on sides of bag for external tourniquet pouches
- Large exterior loop panel for ID patches

Bag only option also available at NARescue.com

K-9 RUF KIT:








































K-9 RUF KIT

OPTION	ITEM#
S-Rolled Gauze	80-1027
with Combat Gauze LE	80-1026







MULTI-MISSION TRAUMA PACK Customizable modular solution ideal for personalized gear organization spanning multiple Mission profiles

• A 3200-cubic inch backpack with (4) compartments

• Multi-mission belt attached for different configurations

• Internal Molle straps covered in hook and loop

• Detachable compartments

• Converts to trauma sheet for vehicles, boats, or aircrafts

• Internal channeling for electrical communications wiring

• Dedicated loops, providing rappelling or lowering of the bag

SPECIFICATIONS

H 25.5 in. x W 11 in. x D 9 in. Weight: 7 lb 8 oz

MULTI-MISSION TRAUMA PACK 'NATO NSN

COLOR ITEM#		NSN#	
Black 80-0950 Coyote 80-0949		6545-01-695-5545*	



Tiered Trauma Medical Bag **TTMB**







A. ReSupply Bag Item# 80-0271 (MTC only)



- 5-piece modular ensemble targeting the military medic's "Soldier-as-a-System" requisite
- Includes low-profile backpack-style aid bag, ample resupply bag, medical chest pouch, leg rig conversion panel & directional MOLLE converter
- Interchangeable internal pouches with removable labels compatible with the NAR-4 system or any hook and loop attachment
- Leg rig converter & directional change MOLLE panel allow the chest pouch to be mounted & carried using a drop-down, leg-attachment style platform

SPECIFICATIONS

- (A) ReSupply Bag: H 20 in. x W 15 in. x D 5 in.
- (B) NAR-4 Aid Bag: H 20 in. x W 15 in. x D 5 in.
- (C) NAR-4 Chest Pouch: H 7.75 in. x W 9.75 in. x D 4.5 in.
- (D) Leg Rig Conversion Panel: L 14 in. x W 6 in.
- (E) MOLLE Conversion Panel: L 7.5 in. x W 6 in.

TIERED TRAUMA MEDICAL BAG (TTMB)

COLOR	ITEM#
Multi-Cam	80-0274
Multi-Cam	80-0271
Multi-Cam	80-0241
Multi-Cam	80-0248
Multi-Cam	80-0272
Multi-Cam	80-0273
	Multi-Cam Multi-Cam Multi-Cam Multi-Cam Multi-Cam Multi-Cam

NOTE: Bags only all contents sold separately





♣ OCHO BAGS Low Visibility | Customizable Modular Solutions

 Low visibility aid bags purposely designed to resemble standard off-the-shelf bags to minimize the threat of theft and maximize low visibility discretion

• Featuring hidden compartments with interchangeable pouches backed with hook and loop fields

• Compatible with CrossBreed Holsters (not included)

OCHO MESH **AID POUCH**

SPECIFICATIONS

H 7 in. x W 4.5 in. x D 3 in. Weight: 1.87 oz

COLOR	ITEM#
Black/ODG	80-0259

OCHO ACCESSORY BAG

SPECIFICATIONS

H 7 in. x W 2.25 in. x D 3.25 in. Weight: 2.75 oz

COLOR	ITEM#	
Blue/Gray	80-0371	
Gray/Black	80-0370	
Red/Gray	80-0571	

OCHO Accessory Bag does NOT include Mesh Aid Pouches


COLOR	ITEM#
Black	80-0260

Gear sold separately









- Oversized pockets & large grommet attachment points facilitate rapid setup and equipment retrieval
- Quick attach and release straps allow for secure attachment to the inside wall of an aircraft or vehicle

SPECIFICATIONS

Open: H 41.5 in. x W 21.5 in. x D 2.5 in. Folded: H 20.25 in. x W 21.5 in. x D 5 in. Weight: 4 lb 12 oz

Bag only; contents are sold separately

MEDIC TRAUMA SHEET

COLOR	ITEM#	NSN#
Black	80-0067	6515-01-540-7612
Coyote	80-0069	6515-01-540-7603
Digital Camo	80-0094	6515-01-540-7579
OD Green	80-0068	6515-01-540-7585



Accessories NAR QUICK DETACH PANELS |



• Mil-Spec hook & loop and heavy duty snaps to ensure kit security

• Designed to be compatible with all of NAR IFAK's, including a **NEW horizontal mount option** in the medium version (W 4.5 x H 6.5 in.)

Small QUICK DETACH PANEL

(3 in. x 6 in.)

- Utilizes only 3 rows vest/bag, vertical mount configuration
- Designed specifically for NAR's most compact M-FAK kits

SPECIFICATIONS

W 3 in. x H 6 in. x D 0.5 in. Weight: 2.8 oz

SIZE	COLOR	ITEM #
3 x 6 in.	Black	80-0893
3 x 6 in.	Coyote	80-0894

by 2 columns on user

SMALL PANEL

Medium DUAL DIRECTION QUICK DETACH PANEL

• Utilizes 3 rows by 3 columns vertical mount or 2 rows by 4 columns horizontal mount

(4.5 in. x 6.5 in.)

• Designed for compatibility with all NAR IFAKs such as horizontal mount for M-FAK and vertical or horizontal mount of the TORK

SPECIFICATIONS

W 4.5 in. x H 6.5 in. x D 0.5 in .Weight: 4.5 oz

DUAL DIRECTION PANEL

SIZE	COLOR	ITEM#
4.5 x 6.5 in.	Black	80-0895
4.5 x 6.5 in.	Coyote	80-0896

Large QUICK DETACH PANEL

(9 in. x 6.5 in.)

- Utilizes 3 rows by 7 columns on user vest/bag, vertical mount configuration
- Designed specifically for use with NAR 4 Chest Pouch & Emergency Evacuation Litter (EEL)

SPECIFICATIONS

W 9 in. x H 6.5 in. x D 0.5 in.. Weight: 5.6 oz

LARGE PANEL

SIZE	COLOR	ITEM#
9 x 6.5 in.	Black	80-0897
9 x 6.5 in	Coyote	80-0898







D-BCRK Department of Defense Individual BLEEDING CONTROL RESPONSE KITS

- Individual Bleeding Control Response Kit designed to provide essential equipment that empowers DOD Forces and Federal Employee responders in public areas to take immediate action in stopping life threatening bleeding.
- Contained in a ruggedized, nylon carry bag or vacuumsealed pouch
- Provides intuitive and easy-to-use tools that are proven to help to save lives
- Available configurations: Basic, Intermediate, & Advanced

SPECIFICATIONS

Nylon Bag: H 7 in. x W 5 in. x D 2.75 in. Vacuum Sealed: H 7 in. x W 5 in. x D 2.75 in.

Nylon Bag Option:

- Easy to restock in cases of use
- Convenient to carry and easy to open in stressful situations
- Opens flat for ease of use and immediate identification of enclosed devices
- Easily recognizable red nylon bag with ample hook & loop field for attaching an ID patch

Vacuum Sealed Option:

- Cost-effective alternative to Nylon
- Compact and rugged Vacuum-Sealed pouch
- Large Identifying Card Vacuum-Sealed inside the pouch for ease of identification

D-BCRK INDIVIDUAL KITS

LEVEL	BAG	WEIGHT	ITEM#
Basic	Red Nylon	15 oz	85-2057
Intermediate	Red Nylon	16.5 oz	85-2058
Advanced	Red Nylon	17.4 oz	85-2059
Basic	Vacuum Sealed	11.7 oz	85-1749
Intermediate	Vacuum Sealed	13.25 oz	85-2066
Advanced	Vacuum Sealed	14 oz	85-2069

Also available at NARescue.com:

#ZZ-1000 Individual Clear Polycarbonate Wall Cabinet

Every D-BCRK INDIVIDUAL BLEEDING CONTROL RESPONSE KIT:











PLUS these Individual Components based upon configuration of choice:









D-BCRK 8-Pack Department of Defense Bleeding Control Response Kits (Nylon or Vacuum-Sealed)



• Eight (8) Nylon or Vacuum-sealed Individual Bleeding Control Response Kits inside an easy-toopen nylon carrying case

• Available in Basic, Intermediate, and Advanced configurations

SPECIFICATIONS

Each 8-Pack:: H 14 in. x W 11 in. x D 8 in. Individual Bleeding Control Kit: H 7 in. x W 5 in. x D 2.5 in.





D-BCRK 8-PACKS

LEVEL	BAG	WEIGHT	ITEM#
Basic	Red Nylon	8.5 lb	85-2072
Intermediate	Red Nylon	9.3 lb	85-2073
Advanced	Red Nylon	9.7 lb	85-2074
Basic	Vacuum Sealed	8.5 lb	85-1750
Intermediate	Vacuum Sealed	9.3 lb	85-2067
Advanced	Vacuum Sealed	9.7 lb	85-2070

D-BCRK 8-PACK - Vac





D-BCRK 8-PACK - Nylon



See page left for contents of the Individual Kits

D-BCRK Wall Station

Department of Defense Bleeding Control Response Kits (Nylon or Vacuum-Sealed)

- Wall Station contains One D-BCRK 8-Pack each containing Eight Nylon or Vacuum-sealed Individual Bleeding Control Response Kits
- Available in Basic, Intermediate, and Advanced configurations
- Includes a durable, clear polycarbonate cabinet with rounded corners provides both physical and environmental protection for medical equipment
- Includes a wire shelf to hold the 8-pack & an external eyelet that allows for a secondary safety seal to easily identify if the station has been opened or tampered with
- External mounting tabs and hardware enable easy installation on any wall
- Clear design allows visibility at any angle

SPECIFICATIONS

Clear Wall Case: H 18 in. x W 15.5 in. x D 10 in. Each 8-Pack: H 14 in. x W 11 in. x D 8 in. Individual Bleeding Control Kit: H 7 in. x W 5 in. x D 2.5 in.

D-BCRK CLEAR WALL STATION:

D-BCRK WALL STATIONS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	23.5 lb	85-2075
Intermediate	Red Nylon	24.3 lb	85-2076
Advanced	Red Nylon	24.7 lb	85-2077
Basic	Vacuum Sealed	23.5 lb	85-1751
Intermediate	Vacuum Sealed	24.3 lb	85-2068
Advanced	Vacuum Sealed	24.7 lb	85-2071









8 Individual

D-BCRK Vac-Pack







Also available at NARescue.com: #ZZ-0921 Semi-Recessed Metal Wall Cabinet







Weight: 15 lb 6 oz (6.97 kg) Static Payload: 1,800 lb (816 kg)

Ground Clearance: 1.5 in.

 ϵ

U.S. Patents: 6,842,923 & 5,598,592, Patent(s) Pending

RAVEN® 90C

Bi-fold Litter

- The official litter of the U.S. Navy and certified safe-to-fly by the U.S. Air Force
- Bifold design for reliable long-term storage
- Ergonomically designed handles
- Slip-resistant, chemical-resistant, rip-stop fabric tested & approved by Federal & private labs using live warfare agents
- Includes 6 IV-attachment points and 2 patient-securing straps

SPECIFICATIONS

Folded: H 45.5 in. x W 7.5 in. x D 8 in. Open: 90 in. x W 22.5 in. x D 6 in. Weight: 16 lb 8 oz (7.5 kg) Static Payload: 2,000 lb (909 kg) Ground Clearance: 1.5 in. U.S. Patents: 6,842,923 5,598,592, Patent(s) Pending

RAVEN 90C LITTER

ITEM#	NSN#
60-0001	6530-01-432-5114

T3™ LITTER

Lightweight, Assault

- Ultra light, ultra strong litter weighing less than 12 lb
- Specifically designed for Special Operations Forces
- Provides maximum maneuverability during TACEVAC in urban terrains

 ϵ

ITEM#

60-0002

NSN#

6530-01-504-9051

- Overall length of 81 in. enables floor loading on a UH-60 helicopter
- 2 straps for securing patients
- Ergonomically designed handles
- Quad-fold design for maximum cube space storage

SPECIFICATIONS

Folded: H 21 in. x W 8.75 in. x D 7.25 in Open: L 81 in. x W 22.5 in. x D 6 in. Weight: 11 lb 13 oz (5.34 kg) Static Payload: 1,200 lb (544.3 kg) Ground Clearance: 1.5 in.

T2 I ITTED

TO EITTER	
ITEM#	
60-0022	





• Lightest wheeled litter cart on the market for single person transport to and from MedEVAC vehicles

• Compatible with the Talon® II 90C, Raven® 90C, MedEvac & T3, in addition to all NATO litters

• Quick and easy to deploy or disassemble for storage

• Locking clamp system keeps the litter firmly in place and allows for different litter widths

 Performs as litter stand in the field due to its unique wheel locking system


 Also available with optional Carrying Case

SPECIFICATIONS

Deployed: L 56.29 in. x W 27.5 in. x H 36.2 in. Folded: L 42.1 in. x W 27.5 in. x H 9.84 in. Case: L 48 in. x W 29.5 in. x H 12.5 in.

Weiaht: Carrier only: 38 lb 8 oz (17.5 kg) Carrier with Case: 88 lb (39.9 kg)

Payload: 550 lb (250 kg)



OPTION	ITEM#	NSN#
Carrier with Case	60-0063	
Carrier only	60-0041	6530-01-663-2308*

LITTER STANDS

33 in. Collapsible

• Set of 2 lightweight, compact litter stanchions designed to collapse for convenient storage in a standard issue military carrying case

· Quickly locks open

• Allows the capacity to store a quad-folded Talon®II 90C (sold separately)

SPECIFICATIONS

H 17.95 in. x W 25.5 in. x D 6.92 in. Extended: 33 in. Base Width: 25.5 in. Weight: 10 lb 12 oz

COLLAPSIBLE 33 in. **LITTER STANDS**

Unit	ITEM#
Set of 2	60-0064

LITTER STANDS

33 in. Standard

- Set of 2 lightweight, foldable stanchions for standard sized litters
- Constructed of a lightweight aluminum frame that conveniently folds flat for ease in transport & storage

• Quickly locks open to safely secure the litter in place

SPECIFICATIONS

H 34 in. x W 25.5 in. x D 6.92 in. Extended: 33 in. Base Width: 23.5 in. Weight: 8 lb 8 oz

33 in. LITTER STANDS

UNIT	ITEM#	NSN#
Set of 2	60-0014	6530-01-523-3314

UNIT	ITEM#	NSN#
Set of 2	60-0014	6530-01-523-3314

BACKREST TALON® or RAVEN®

- Allows patients on a Talon® 90C or Raven® 90C litter to be placed at 90° or 45° angles
- Constructed of a lightweight aluminum frame with Nomex/Kevlar-blend fabric covering

SPECIFICATIONS

H 31.5 in. x W 29.25 in. x D 1.6 in. Weight: 2 lb 13 oz

LITTER BACKREST

ITEM#	NSN#
60-0027	6530-01-594-7444



IV POLE

Exclusively designed to fit the pre-molded holes of our framed litters

SPECIFICATIONS

H 19-33 in. x W 8.75 in. x D 0.6 in. Weight: 4.8 oz

IV POLE

	ITEM#	NSN#	
	60-0008	6530-01-504-9054	

TACEVAC | Litter Accessories

RIG SERIES CARRIER Assault Litter Carrier

- Enhanced carrier with additional strap over the top of the litter to hold it in place during drastic movement and FRIES (Fast Rope Insertion Extraction System) operations
- Upgraded shoulder straps provide comfort for long wear
- Includes pocket in the top flap for storage and three additional handles facilitate packing and unloading of the litter
- Accommodates all NAR quadfold litter platforms

 Equipment attaching webbing (MOLLE) along the outside of the carrier allows attachment of additional equipment

SPECIFICATIONS

H 22 in. x W 9 in. x D 9 in. Weight: 3 lb (1.36 kg)

RIG SERIES LITTER CARRIER

COLOR	ITEM#	NSN#
Multi-Cam	60-0057	6530-01-647-8571*
Coyote	60-0056	6530-01-643-9138*

*NATO NSN





TALON® CARRIER

Assault Litter Carrier

- Transport any NAR quad-fold litter (Talon, MedEvac & T3)
- MOLLE-style attachment points
- Includes thin back plates, drain hole, non-slip, adjustable shoulder straps & top-carrying handle

	TALON LITTI	ER CARRIE	R *NATO NSN
	COLOR	ITEM#	NSN#
	Coyote	60-0011	6530-01-531-6854
	OD Green	60-0010	6530-01-504-9055
	Black	60-0009	6530-01-531-6907
	Digital Camo	60-0025	6530-01-643-9134*

EQUIPMENT BAG

Store a casualty's personal items



LITTER BUCKLE Replacement Set

SPECIFICATIONS

L 3.75 in. x W 2.38 in. x D 0.5 in. Weight: 1 oz

ORDERING INFO

ITEM#	NSN#
ZZ-0032	6530-01-435-7337



ITEM#	NSN#	
80-0132	8465-01-575-6644	





MEDEVAC LITTER

- The ultimate, cost efficient casualty evacuation stretcher for First Responders operating in Active Shooter, Mass Casualty and Disaster incidents
- Compact, collapsible, quad-fold litter is easy to deploy and carry
- Two patient securing buckles, six IV-attachment points, ergonomically designed handles & patented joint configuration with locking leg crossbar

SPECIFICATIONS

 ϵ

Folded: H 21.5 in. x W 8.5 in. x D 7.5 in. Open: H 84.5 in. x W 22.5 in. x D 4 in. Weight: 13 lb 8 oz

Static Payload: 1,200 lb (544.3 kg) Ground Clearance: 1.5 in.

MEDEVAC LITTER

BED COLOR	FRAME COLOR	ITEM#	
Black	OD Green	60-0051	
Orange	Black	60-0052	

OSL™ LITTER OverSized Litter

- Safely evacuate larger patients with 30% larger surface area compared to standard-sized evacuation platforms
- Bi-fold configuration enables lower cube storage capability
- 6 litter stirrups, 6 padded carrying loops & 3 patient securing straps for firmer patient support as well as better weight distribution on the floor
- Flame-retardant cover & attached storage pouch made of a durable Nomex/Kevlar-blend
- Certified safe-to-fly by the U.S. Air Force

OSL OVERSIZED LITTER

ITEM#	NSN#	
60-0023	6530-01-552-7708	

SPECIFICATIONS Folded: H 45 in. x W 9 in. x D 8 in. Open: L 90 in. x W 31 in. x D 6 in. Weight: 22 lb 6 oz (10.2 kg) Static Payload:: 2,000 lb (909 kg) Ground Clearance: 1.5 in. U.S. Patents: 5,598,592 & Patent(s) Pending

OSL™ LITTER STANDS

for the OverSized Litter

- Set of 3 lightweight aluminum stanchions accommodate our oversized OSL™ Litter
- Quickly opens and securely locks the litter in place
- Easily folds for transport & storage

SPECIFICATIONS

CE

Folded: H 34 in. x W 25 in. x D 2.5 in. Extended Height: 32 in. Base Width: 23.5 in. Weight: 28 lb 8 oz (12.93 kg)



UNIT	ITEM#
Set of 3	60-0029

OSL™ BACKREST

for the OverSized Litter

- Position patient at 45° & 90° angles for maximized treatment options & increased comfort
- · Constructed of a lightweight aluminum frame with durable Nomex /Kevlar-blend fabric cover
- Compatible with our OSL™ litter

SPECIFICATIONS

H 31.5 in. x W 30.25 in. x D 1.6 in. Weight: 2 lb 15 oz (1.33 kg)

OSL BACKREST

CE

ITEM# 60-0024

TED Tactical Extrication Device

TACEVAC | Non-Rigid Litters

- Designed to allow two people to carry a casualty to safety
- Lighter, more compact non-rigid evacuation platform designed specially for use in high-threat environments & limited spaces
- Constructed of heavy duty vinyl with 4 sturdy carrying handles

 ϵ

Packaged in a sturdy, transport-friendly bag with adjustable shoulder strap

SPECIFICATIONS

Folded: H 5 in. x W 10.25 in. x D 3.5 in. Open: L 43.5 in. x W 19 in. Weight: 1 lb 8 oz (0.68 kg) Payload: 600 lb (272 kg)

TACTICAL EXTRICATION DEVICE

COLOR	ITEM#	NSN#
OD Green	60-0006	6530-01-492-4114
Black	60-0005	6530-01-531-6808



R.E.E.L. Rapid Expeditionary Evacuation Litter

- Cost-effective, non-rigid casualty evacuation platform designed for quick and effective casualty evacuation
- Rugged, lightweight, nylon construction
- 8 carry handles accommodate mission teams ranging from 2 to 8 rescuers
- Packaged in a compact, small cube, vacuum sealed pouch

SPECIFICATIONS

Packaged - L 13 in. x W 10 in. x D 0.75 in. Litter (Open) - L 78 in. x W 30 in. | Weight: 1 lb 3 oz Carrying Capacity: 700 lb

R.E.E.L. LITTER

ITEM#	
60-0075	

NAR K-9 LITTER & Carrying Strap Accessories

- Lightweight, compact litter for the movement of an injured canine
- Litter can be carried in cargo pocket or pouch
- Leg cut outs to safety carry the canine in the litter
- Six large and four small carry handles that allow for one to four person carries
- Optional K-9 Carry Strap & 2 Carabiners for over the shoulder carry of injured K-9 also available

SPECIFICATIONS

Litter: W 35 in. x H 48 in. Adjustable Carry Strap: L 9 to 54 in. x W 2 in. Carabiner open gate clearance: 27 mm Packaged: W 5.5 in. x H 10 in. x D 2.5 in. Weight: 1.38 lb

K-9 LITTER

OPTIONS	ITEM#
K-9 Litter Only	60-0074W
Carry Strap & Carabiners	ZZ-1004
K-9 Litter with Carry Strap & Carabiners	60-0076



• Ideal for decontamination roller systems & moving wet patients • Non-skid, flame-retardant

monofilament polypropylene fabric that is resistant to mildew, acids, alkalis and penetrating liquids • Tested by Federal & private labs

using live warfare agents Includes six heavy-duty carrying handles & 2

casualty securing straps

• Easy to store and clean

SPECIFICATIONS

Folded: H 20 in. x W 12 in. x D 2 in. Open: L 78 in. x W 24 in. Weight: 13 oz (0.37 kg) Payload: 1,000 lb (453.6 kg)

STINGRAY POLELESS LITTER

ITEM#	NSN#
60-0004	6530-01-492-4079



TACEVAC Spine Boards & Accessories

RESCUE SLEEVE II

- The industry's ultimate patient securing system refined for enhanced stability & ease of use
- Ideal for rescues requiring vertical and/or horizontal hoisting capabilities

SPECIFICATIONS

Folded: L 25.5 in. x W 13 in. x D 8 in. Open: L 73 in. x W 18 in. x D 2 in. Weight: 16 lb (7.26 kg) Payload: 1,000 lb (453.6 kg)

RESCUE SLEEVE II

COLOR	ITEM#	NSN#
Black	50-0011	6530-01-531-6649
Orange	50-0012	6530-01-477-8525

NOTE: Spine board sold sengrately, Hoisting straps & carabiners not

SPINE BOARD

- Thin & lightweight with tapered design
- Surface texture to prevent slippage & 14 raised, oversized hand-holds
- Black spine board includes ten rapid attachment pins

SPECIFICATIONS

L 72 in. x W 16 in. x D 1.75 in Weight: 14 lb (6.35 kg) Payload: 1,100 lb (498.96 kg) Float Capacity: 250 lb (113.4 kg)

SDINE POADD

SPINE BUAK	<u> </u>	
COLOR	ITEM#	NSN#
OD Green	50-0013	6530-01-490-2487
Orange	50-0014	6530-01-531-6771
Black w/Pins	50-0021	

NAR CASUALTY IMMOBILIZATION SYSTEM

- Constructed of black polypropylene webbing with hook & loop closures
- Featuring 6 attachment points: 2 "suspender" type, over-theshoulder straps & 4 adjustable torso/limb straps with adjustable center strap
- No choking "V" configuration by the patient's neck

CASUALTY IMMOBILIZATION SYSTEM

ITEM#	NSN#
50-0031	6530-01-586-2411



Pkg: L 16.5 in. x W 4 in. x D 2.75 in. Weight: 1 lb 1 oz

NAR HEAD **IMMOBILIZER**

- Quick, effective head block constructed of soft, closed cell foam with hook & loop closures and padded straps
- Ears and head are accessible without removing the device

SPECIFICATIONS

H 10 in. x W 16 in. x D 5.75 in. Weight: 2 lb 3 oz

HEAD IMMOBILIZER

ITEM #	NSN #
50-0020	6530-01-619-1777

NAR CASUALTY **RESTRAINT STRAPS**

Heavy-duty nylon webbing & convenient seat belt buckle utility

- A. 2 piece Plastic Buckle L5 ft x W 2 in. • Weight: 4 oz Black: 50-0023 · Orange: 50-0024
- B. 2 piece Plastic Buckle/Speed Clip L 5 ft x W 2 in. • Weight: 7.8 oz Black: 50-0026 · Orange: 50-0025
- C. 2 piece Metal Buckle L 5 ft x W 2 in. • Weight: 8 oz Black: 50-0027 · Orange: 50-0028
- D. 2 piece Metal Buckle/Speed Clip L 5 ft x W 2 in. • Weight: 11.7 oz Black: 50-0029 · Orange: 50-0030



CERVICAL COLLAR ACE®

CE

- Adjustable neck support with foam lining
- Extended back to ensure complete
- Adjusts to fit 16 adult neck sizes

SPECIFICATIONS

H 6.75 in. x W 20.5 in. x D 0.75 in. Weight: 5.2 oz

ACE CERVICAL COLLAR

COLOR	ITEM#	NSN#
OD Green	50-0010	6515-01-541-8147
White	50-1007	6515-01-650-0036









MIGHTYSAT™RX
Fingertip Pulse Oximeter from Masimo

• High Resolution Display

• Plethysmographic Waveform

features Masimo SET® Measurethrough Motion and Low Perfusion™ Pulse Oximetry which overcomes the challenges of low perfusion and motion artifact that limit conventional pulse oximetry.

This Signal Extraction Technology® (SET®) uses five parallel signal processing engines that measure through challenging conditions allowing users to accurately measure oxygen saturation and pulse rate even during movement or low blood flow.

> **SPECIFICATIONS** Device Dimensions:

H 2.9 in. x W 1.6 in. x D 1.2 in.

 $(7.4 \text{ cm} \times 4.1 \text{ cm} \times 3.0 \text{ cm})$

Weight (with batteries):

2.56 oz (73 g)

Display: 3.3 cm

patient's pleth waveform

 Touchpad Provides easy access of settings and device features

Rotating and adjustable color screen for

real-time display of all measurements

Provides high resolutions display of the

• Durable, Lightweight Design

Small and portable for use in clinical and home-care environments

• Flexible Pad

Easily conforms to finger shape for enhanced comfort & secure application to finger to minimize external interference

• Signal IQ (SIQ™)

An indicator of the SpO₂ measurement confidence

• Bluetooth Models (20-0094 & 20-0095) are compatible with free downloadable Masimo Professional Health App

KIT CONTENTS:

- 1 x MightySat™Rx Device
- 1 x Lanyard
- 2 x AAA Batteries
- 1x Carrying Case

MIGHTYSAT PULSE OXIMETER

OPTION	ITEM#
Standard	20-0093
Bluetooth LE	20-0094
Bluetooth LE/RRp/PVi	20-0095



FEATURES	20-0093	20-0094	20-0095
Masimo SET® Measure-through Motion and Low Profusion™ Technology	1	1	1
Measurements SpO ₂ (Oxygen Saturation)	✓	✓	✓
PR (Pulse Rate)	√	✓	/
PI (Perfusion Index)	✓	√	✓
PVI (Pleth Variability Index)			√
<i>Display</i> Pleth Waveform	✓	✓	✓
3.3 cm Color Display	√	✓	√
Rotational Screen	✓	✓	/
Bluetooth LE		√	✓
Compatible with Masimo Professional Health App*		√	1
Four Year Limited Warranty	1	√	✓

MED/SURG Diagnostics

FIELD CORPSMAN KITS

Portable ENT diagnostics

- Lightweight, durable ENT diagnostic sets for portable field care of the eyes, ears & nose
- A wide selection of high intensity lights & magnifiers for enhancing both day and night examinations

FIELD CORPSMAN KITS:































Additional items in Deluxe Kit only







SPECIFICATIONS

H 8 in. x W 4.25 in. x D 2.25 in. Weight: Basic Kit: 6.2 oz Deluxe Kit: 11.6 oz Ophthalmoscope only: 3.7 oz

FIELD CORPSMAN KITS

KIT OPTIONS	ITEM#	NSN#
Basic	20-0001	6515-01-458-6057
Deluxe (with Heine)	20-0002	6515-01-458-6178
Ophthalmoscope only	20-0007	6540-01-491-6651

BP/STETHOSCOPE COMBO KIT from NAR

- Ultra-compact, multi-pocketed BP kit that can be stored in most trauma kits
- Blood-pressure cuffs in 4 sizes (adult, large adult, child, infant).
- ADC stethoscope with extra-large bell & ultra-sensitive diaphragm
- Ambidextrous palm gauge


BP/STETHOSCOPE COMBO KIT

ITEM#	NSN#
20-0057	6515-01-542-7723

SPECIFICATIONS

H 6.75 in. x W 9 in. x D 5 in. Weight: 2 lb

HEINE MINI 3000 Ophthalmoscope





 Enhanced XHL Xenon Halogen Technology; concentrated white light for perfect illumination & crisp retinal detail

HEINE MINI 3000 SPECIFICATIONS

ITEM# NSN# H 6.75 in. x W 1.25 in. x D 0.9 in. 20-0007 6540-01-491-6651 Weight: 4 oz

TETHOSCOP

- Best full-range stethoscope for auscultation & comfort
- Superior acoustics for accurate auscultation of heart, lung & blood-pressure sounds

SPECIFICATIONS

Device Length: 30.5 in. Weight: 13.5 oz

STETHOSCOPE

ITEM#

20-0060





ARMADILLO

Medication & Storage Case

- Compact solution for keeping critical medications safe, secure and easily accessible
- Features Automatic Pressure Equalization Valve
- Hardware constructed of Marine grade 316 SS
- Case is water-resistant and crush-proof
- Hook & loop attachment system allows storage compartment to be tailored to hundreds of custom configurations

SPECIFICATIONS

Ext: L 7.95 in. x W 5.1 in. x D 2.14 in. Int: L 6.45 in. x W 3.87 in. x D 1.75 in. Weight: 10.5 oz

ARMADILLO	0	* NATO NSN
LID STYLE	ITEM#	NSN#
Black	ZZ-0081	8145-01-573-2533
Clear	ZZ-0259	6545-01-594-9399*

also available: Half Inch Hook Discs with pressure sensitive adhesive Set of 250 · Item# ZZ-0258 · NSN# 8145-01-617-9075*

I&A CASE Injury & Ailment Case

 Superbly engineered medical equipment organization system

Clear Lid

ZZ-0259 only includes Loop Linei

on the bottom

- Lid's clear pockets and hook & loop fasteners allow secure storage as well as superior inventory management
- Polymer O-ring provides dust and water-proof seal
- Automatic Gore-TEX® pressure equalization valve stops moisture from entering the case - also prevents vacuum-lock
- Double-throw latches engineered to tolerate up to 400 lb of pressure
- All-weather, sure-grip handle and reinforced padlock protectors



M.U.L.E. II

Multi-Use Luggage & Equipment Carrier II

- Extra large, heavy-duty, water-resistant 1000 denier nylon platform featuring 9,000 cubic in. of storage with large reinforced open compartment and multiple zippered mesh pockets for maximum storage
- Top & side web carry handles and 3 durable 3.75 in. diameter wheels for easy transport

Available in Black & Coyote

M.U.L.E. II COLOR ITEM# ZZ-1157T Black Coyote ZZ-1158T

SPECS Exterior: H 40 in x W 21 in. x D 15 in. Interior: H 40 in. x W 15 in. x D 15 in. Weight: 15.4 lb Capacity: 9,000 cubic in.

M.U.L.E. Multi-Use Luggage

& Equipment Carrier

- Heavy-duty, waterresistant nylon platform with large compartment and mesh pockets on the inside lid for maximum storage of critical items
- Outside end pockets feature open channels for expanded storage capability
- Suitcase-style carrying handles on top and side; Telescoping handle and wheels for easy transport
- Hidden rucksack shoulder pads



SPECIFICATIONS

H 30.5 in. x W 18.5 in. x D 11 in. Wgt: 11 lb 10 oz | Payload: 70 lb

M.U.L.E.	
ITEM#	
ZZ-0212	



• Compact, lightweight precision lighting solution with the rugged engineering required for use in field hospitals

- LED Technology with intense light fifteen times brighter than traditional halogen field lights
- Patented 3-pod head eliminates shadow interference common with single point light sources
- · Energy efficient, cool white color with an exceptionally long life
- Stays cool to the touch, projecting no heat on the patient
- Locking joints make it fully adaptable and secure with adjustable mounting clamp for attachment to hospital bed, litter or litter stand
- Includes a rugged storage case with wheels for easy transport

SPECIFICATIONS CE



Light in case: L 47.2 in. x W 16.5 in. x D 6.7 in. Light Head: 12.8 in. Diameter Light Arm/Post: L 78 in. x 2.8 in. Diameter Weight: Case & Light: 37 lb 8 oz Weight: Light with Arm/Post: 11 lb 8 oz

Light Output: 50,000 lux (4590 foot candles)

PAXLIGHT FIELD LIGH

60-0042 6530-01-549-108	9



Assembled: H 38 in. Weight: Base: 13.5 lb Assembled: 22 lb

WHEELED BASE

ITEM#	
60-0045	

SPECIFICATIONS

L 4.5 in. x W 1 in. x D 1 in. Shield Diameter: 4.25 in.

EZ LIGHT HANDLES

UNIT	ITEM#
Pack of 10	60-0046

MARK IV Field Hospital Bed

- Lightweight solution designed for patient comfort in the combat support or field hospital/treatment tent
- Durable & easy-to-open, supporting up to 850 lb of static weight
- Designed to be rapidly deployed with minimal training in less than 60 seconds by a single person

SPECIFICATIONS

Folded: L 42 in. x W 17 in. x D 32 in. Opened: L 79.25 in. x W 33 in. x H 35.25 in.

Weight:

Mattress and Frame: 54 lb Mattress, Frame & Accessories: 68 lb U.S. Patents Pending

- Constructed of a bi-fold aluminum frame that is adjustable to 6 positions
- 3 in. polyurethane foam pad with water-resistant decontaminable cover easily cleans with 10% bleach solution or other standard hospital cleaning products
- Includes IV Pole and a removable table platform that rests on its adjustable bed rails

MARK IV FIELD HOSPITAL BED

ITEM#	NSN#
66-0019	6530-01-624-7626



Wound Care | Sutures | MED/SURG





SPECIFICATIONS

43mm x W 12mm

Box Wt.: (24) 3.5 oz

X-SMALL: Ideal for small areas, such

as the fingers, toes, ears, and nose

Large

Device:

microMend 2 STERILE DEVICES XS 🔟

microMEND

Skin Closure Devices

- Easy-to-use alternative to sutures or surgical staples when skin closure is required
- Safe, secure, and rapid skin closure
- Equivalent holding strength to sutures
- Patent pending microstaples anchor the device to the skin

X-Small

Device:

microMend

L 28mm x W 7mm

Box Wt: (24) 3.5 oz

- Three times faster to apply than sutures
- More comfortable and easier to remove than standard invasive methods of skin closure
- Can be used on a variety of lacerations and incisions and with wounds of various lengths, as the product can be applied right next to each other

STEP 1:



Remove device from packaging by pulling TAB upwards

STEP 2:




Position device so the edge of narrow section on the TAB-side is over wound. Press MICROSTAPLE-side onto skin.

Wide

Device: L 38mm x W 28mm Box Wt: (24) 6.9 oz



STEP 3:



Pull TAB-side across the wound. Press MICROSTAPLE of tab-side down onto skin. Remove the tab with a twist.

microMEND[™]Wide

Small

Device:

microMend

L 32mm x W 12mm

Box Wt: (24) 3.5 oz

(48) 7.1 oz

- Three times wider for easier and faster closure of longer incisions with the same holding strength as other microMend sizes
- Ultra-flexible and elastic making it ideal for use in joint surgeries, Cesarean sections, and other procedures in areas of increased movement and varying skin architecture
- Available with apertures (Wide A) to allow drainage, or without apertures (Wide S) for a solid seal to keep the wound clean from contaminants

LARGE: For closing larger wounds on the extremities, back, chest, and abdomen. Great for closing hard to repair skin tears

microMend

microMEND DEVICES

DEVICE	UNIT (per box)	ITEM#
Large	24	30-0203
Small	24	30-0230
Small	48	30-0249
X-Small	24	30-0252
Wide A	24	30-0295
Wide S	24	30-0296

U.S. Patents Pending

Available with Aperature or Solid Seal

MED/SURG | Surgical Field Kits



R-SICK CONTENTS:

- 2 x Weitlaner Retractor, 3 x 4 Blunt (5.5 in.)
- 1 x Trousseau Trachael Dilator (5.25 in.)
- 2 x Adson Forceps, 1 x 2 (4.75 in.)
- 1 x Ewald Forceps 1 x 2
- 1x Plain Splinter Forceps (4.5 in.) Serr. Handle
- 1 x Kelly Forcep Heavy
- 1 x Converse Needle Holder (4.24 in.) Long
- 1 x Olsen-Hegar Needle Holder (7.25 in.) Serr
- 1 x Stevens Tenotomy Scissors Trbl/Bl
- 1 x Stevens Tenotomy Scissors Cvdbl/Bl
- 1 x Scissors, Mayo (5.5 in.) Str
- 1 x Mayo Dissecting Scissors, Crv (5.5 in.)
- 1 x Knife Handle (#3)
- 2 x Gigli Saw Blade (20 in.)
- 2 x Gigli Saw Handle (3.5 in.) Loop
- 1 x Hiupp Tracheal Hook, (6.5 in.)
- 2 x Forcep, Rochester-Pean (7.25 in.) Curv
- 6 x Scalpel #10 w/Handles

- 6 x Scalpel #11
- 6 x Scalpel #15
- 6 x Blade, #10 Pluss/S(Tefl)
- 10 x Blade, #15 Pluss/S(Tefl)
- 2 x Suture Ethilon 4-0 Blk Mono Fs-218
- 2 x Suture 3-0 Silk
- 2 x Black Braided Silk Suture, 2-0/2-60 in.
- 2 x 2-0 Suture Ct-2
- 4 x 4-0 Silk Suture Blk (30 in.)
- 1 x Sharps Shuttle
- 1x Farabeuf Retractor

P-SICK

IX-SICIX		
INSTRUMENTS	ITEM#	
Floor Grade	83-0016	
OR Grade	83-0017	
Bag Only	80-0871	



FBTK Field Blood Transfusion Kits

FBTK Basic:

Donor Module:

5 x Alcohol Prep Pad

2 x 2-pack Gauze Pad 4 in. x 4 in.

1 x IV Constricting Band - 1 in. x 18 in.

1 x Clear Adhesive Dressing - 4 in. x 4.75 in.

1 x Blood Bag Label

Recipient Module:

5 x Alcohol Prep Pad

1 x IV Catheter - 18 ga x 1.25 in.

2 x 2-pack Gauze Pad - 4 in. x 4 in.

1 x IV Constricting Band - 1 in.x 18 in.

1 x Clear Adhesive Dressing - 4 in. x 4.75 in.

1 x IV Catheter - 16 ga x 1.25 in.

Additional Supplies:

1 x roll 1 in. Surgical Tape

2 x pair Black Nitrile Gloves, Large

1x Black Permanent Marker

1 x Kelly Curved Hemostat - 5.5 in.

1 x Single Blood Bag (450ml)

1 x 550 Cord - 9.5 in.

1 x Beaded Security Tie - 11 in.

1 x IV Set Y-Type with Filter 170-260 UM

1x Questionnaire Card - DD572 Modified

FBTK Blood Surveillance Module:

Donor Module:

1 x IV Hypodermic Needle - 18 ga x 1.5 in.

1 x IV Syringe, 10cc

1 x Eldon Card

1 x Interlink Injection Site Needle

1 x Blood Transfer Device

2 x Vacutainer, Gold Top (5ml)

1 x Vacutainer, Purple Top (5ml)

2 x Vacutainer, Pearl Top (5ml)

Recipient Module:

1 x 3 Way Stopcock

1 x Eldon Card

1 x Interlink Injection Site Needle

*FBTK Complete:

Contains all components of

FBTK Basic & Blood Surveilleance Module

FIELD BLOOD TRANSFUSION

OPTIONS	ITEM #
FBTK Basic	83-0022
Blood Surveillance	83-0025
FBTK Complete*	83-0000

Surgical Kits | MED/SURG



MINOR SURGERYSets

MINOR SURGERY Contents:

1x Tray, Shallow

2 x Kelly Forceps, Straight 5.5 in.

2 x Kelly Forcep, Curved 5.5 in.

1x Mayo Hegar Needle Holder, 6 in.

1x Adson Forceps, Serr 4.75 in.

1x Forceps Tissue, 5 in.

1x Collier Needle Holder, Straight 5 in.

1x Handle Surgical Knife #3

1x Probe W. Eye 5 in.

1 x O.R. Scissors, Sb St 5.5 in.

1x Wrap, Kimguard 24 x 24 in.

SPECIFICATIONS

Packaged: L 9 in. x W 5 in. x D 1 in. Weight: 5 oz



MINOR SURGERY SETS

ITEM#	UNIT	NSN#
ZZ-0960	Single	
ZZ-0959	Pack of 20	

SURGICAL SUTURE Kits

SURGICAL SUTURE Contents:

1 x Tray, Shallow

1x Adson Forceps w/Teeth

1x Webster Needle Holder

1x Iris Scissors, Straight

1x Mosquito Hemostat, Curved

1x Poly Bag w/Twist Tie

1x Poly-backed Towel

1x Fenestrated Drape

1x Alcohol Prep Swab

1x PVP Prep Swab

1x Cotton Filled Sponge.

SPECIFICATIONS

Packaged: L 9 in. x W 5 in. x D 1 in. Weight: 6 oz



ITEM#	UNIT	NSN#
ZZ-0958	Single	
ZZ-0961	Pack of 20	6515-01-311-0360

TOENAIL REMOVALTrays

TOENAIL REMOVAL Contents:

1x Tray, Half Deep

1x CSR Wrap, Kimguard 24 in. x 24 in.

1x Kelly Forcep, Straight 5.5 in.

1x Kelly Forcep, Curved 5.5 in.

1 x Adson Forceps, 1 x 2 Teeth, 4.75 in.

1 x Iris Scissors, Straight 4.5 in.

1x Cuticle Nipper 4 in.

1x Probe w/Eye 5.5 in.

10 x 12-Ply Gauze, 4 in. x 4 in

1 x Kwiktear Header Bag, 10 in. x 15 in.

SPECIFICATIONS

Packaged: L 9.72 in. x W 5.4 in. x D 2 in. Weight: 11 oz



DENAU DEMOVAL TRAVE

IOENAIL REMOVAL IRATS			
ITEM#	UNIT	NSN#	
ZZ-0956	Single		
ZZ-0957	Pack of 10	6515-01-630-8050	




ED/SURG Suture Kits & Instruments



SUTURE REMOVAL KIT #1









ITEM#

ZZ-0664

SPECIFICATIONS

STERILE

Pkg: L 6 in. x W 3.5 in. x D 0.38 in. Weight: 2 oz



SUTURE REMOVAL KIT #2













ITEM#

ZZ-0665

SPECIFICATIONS

Pkg: L 6 in. x W 3.5 in. x D 0.38 in. Weight: 2 oz



SUTURE REMOVAL KIT #3











ITEM#

ZZ-0666

SPECIFICATIONS

STERILE

Pkg: L 6 in. x W 3.5 in. x D 0.38 in. Weight: 2 oz



SUTURE REMOVAL KIT #4





STERILE







Curved

6.25 in.

HALSTEAD

ITEM#

77-0667

SPECIFICATIONS

Pkg: L 6 in. x W 3.5 in. x D 0.38 in. Weight: 2 oz

SURGICAL

Precision standard surgical instruments available individually














TTEM #

ZZ-0702



STERILE

Laceration Trays | MED/SUR



FACIAL LACERATION TRAY

- Cost effective, sterile, single use kit for facial lacerations with a 3 year shelf life
- Conveniently packaged in a peel-to-open pouch for easy access



























LACERATION TRAY (MINOR



Scalpel

Towel

MINOR LACERATION TRAY

- Cost effective, sterile, single use kit for facial lacerations with a 3 year shelf life
- Conveniently packaged in a peel-to-open pouch for easy access











ITEM#	NSN#
ZZ-0684	6515-01-661-1786

SPECIFICATIONS

Pkg: L 9 in. x W 5 in. x D 2.75 in. Weight: 12 oz



1 Crvd Iris

TOENAIL REMOVAL KIT

- Cost effective, sterile, single use kit for facial lacerations with a 3 year shelf life
- Conveniently packaged in a peel-to-open pouch for easy access to components

ITEM# ZZ-0685 **SPECIFICATIONS** STERILE Pkg: L 9 in. x W 5 in. x D 2.75 in. Weight: 12 oz























PEGASUS ROPE®



8MM HAVOC Assault Rope

Innovative Synergy of 2 unique Fibers engineered for unmatched Confidence & Resilience from Ropes in Challenging Environments

North American Rescue & Edelrld. the industry leader in rope manufacturing with 150 years of experience, have joined forces with multi-disciplinary vertical specialists to innovate the *most tenacious* small diameter rope on the market

today - engineered to counter the most severe environmental pathologies found throughout full-spectrum operations... whether at work or at play.

The unique material engineering of the 8mm HAVOC core, specifically it's Innegra™ Braided Fiber Core **Integration**, provides **superior** durability, low fatigue properties, an extremely low specific gravity, and superior resistance to the harshest contaminates (acids, bases, and organic solvents). The novel tensioning and fiber refinement techniques utilized within the **12-braid Innegra™ core** amplfies its overall robustness and resilience providing an operational rope which acts as a counter-measure to anticipated environmental challenges.

This specialized core of the 8mm HAVOC is protected by a higher density, precision double braided **Technora sheath** to synergistically interface with the core for **optimized strength**, elongation, decreased creep, heat resistance, and operational safety. To reinforce the abrasion resistance, cut resistance, and overall durability, the 8mm HAVOC's sheath was bulked up to a 52% sheath ratio.





LENGTH	ITEM#
75 ft	90-0403
100 ft	90-1402
150 ft	90-0404
200 ft	90-2402
300 ft	90-3402
600 ft	90-6402

Specifications:

- Diameter: 8 mm
- Weight: 2.9 lb/100 ft
- Elongation: 300 lb: 3.58%
- 3-Sigma MBS: 5,615 lb (25 kN)
- Fibers (Core/Sheath): Innegra™/Technora®
- Sheath Ratio (by wt.): 52.5%

CORE

• *Innegra™* Lightest Weight Structural Fiber on the market with a very Low Specific Gravity

INNZGRAZ

CORE

52%

Technora®

Sheath

- Hydrophobic
- Resistant fo Harsh Contaminents
- Low Flex Fatigue
- Superior Durability
- Meets EN 564

SHEATH

• 52% Technora Sheath

Ratio (by weight) The most robust on the market!

- Resistant to Cuts. Abrasion & Heat
- Optimized Strength
- Decreased Creep
- Increased Durability & Safety



WARHORSE™ SERIES 100% Innegra™ fiber braided core

- High quality, multipurpose rope
- Static rope provides reliability, high performance and great handling
- Excellent knot-ability with high strength-to-weight ratio
- Low stretch/sheath slippage wet, dry or in any other climatic condition
- Ultra-lightweight, strong Innegra[™] core for reliable strength, durability & flotation
- Properties remain unchanged even when wet, dry or contaminated with chemicals
- Low dielectric and UV stable

SPECIFICATIONS

Diameter Weight: Tensile Strength: Fibers (core/sheath): Sheath Ratio: Ratings:

9.1 mm 4.88 lb/100 ft 5,670 lbf (22.7 kN) Innegra™/Polyester 63% Polyester NFPA 1983/2012 Throwline (Pending)

9.5 mm 3.34 lb/100 ft 5,100 lbf (22.7 kN) Innegra™/Polyester 28% Polyester NFPA 1983/2012 Throwline (Pending)



ORDERING INFORMATION

DIAMETER	COLOR	RUNNER	ITEM#
9.1 mm	Black	Gray	90-0207
9.5 mm	Black	Gray Reflective	90-0206
9.5 mm	Yellow	Blue Reflective	90-0205
9.5 mm	White	Blue	90-0208



OLYMPUS™ SERIES Engineered for rough terrain • Ideal for use in challenging and abusive environments • Durable, lightweight and great handling offers reliability in all weather conditions

• Technora sheath adds strength and abrasion resistance

• Buoyant even with a 62% Technora sheath (9.5mm)

SPECIFICATIONS

Diameter 9.5 mm 11 mm 3.82 lb/100 ft 4.99 lb/100 ft Weight: Tensile Strength: 5.495 lbf (24.4 kN) 6.834 lbf (30.4 kN) Innegra™/Technora Innegra™/Technora Fibers (core/sheath): Sheath Ratio: 62% Technora 60% Technora

ORDERING INFORMATION			
DIAMETER	COLOR	RUNNER	ITEM#
9.5 mm	Gold	Black/Red	90-0401

Black/Blue 90-0400

Gold

11 mm



STALLION[™] SERIES Low visibility - high performance

- 100 ft of this 9.5 mm rope offers 5,100 lb of strength under 3 lb - wet or dry
- Chemically resistant to acids, bases & common industrial chemicals
- Low sheath slippage, excellent strength-to-weight ratio and hydrophobic

SPECIFICATIONS

Diameter 9.5 mm Weiaht: 2.9 lb/100 ft 5,100 lbf (22.7. kN) Tensile Strength: Innegra™/Innegra™ Fibers (core/sheath): Sheath Ratio: 24% Innegra™

TRITON™ SERIES Specially designed for water rescue

- 10% lighter & 250% stronger than polypropylene alone
- 100% Innegra™ fiber braided core, the lightest structure fiber available, provides additional buoyancy, low elongation and consistent strength

SPECIFICATIONS

Diameter 75 mm 2 lb/100 ft Weight: Tensile Strength: 3,050 lbf (13.6. kN) Fibers (core/sheath): Innegra™/Polypropylene Sheath Ratio: 46% Polypropylene NFPA 1983/2012 Ratings: Throwline (Pending)

ORDERING INFORMATION

DIAMETER	COLOR	RUNNER	ITEM#
7.5 mm	Yellow	Blue Reflective	90-0303
7.5 mm	Red	Yellow	90-0305



R.A.C.E. SYSTEM Rapid Access Casualty Extraction

- Lightweight, multiple use casualty evacuation system enabling first responders/medical providers to ascend, descend, hoist & lower both rescuer & casualty to safety
- Combines 150 ft of lightweight durable 9.5 mm Olympus Pegasus Rope® with the latest in lightweight multipurpose hardware
- Capable of rigging 2:1, 3:1, 4:1, 5:1, 6:1, and 9:1 systems in different configurations for a wide variety of rescuer options with fully illustrated instruction booklet
- Front pocket allows for optional storage
- Carriage options include 4 Point, D-ring configurable modified side and/or backpack carry-style adjustable strap

R.A.C.E SYSTEM CONTENTS:















COLOR	ITEM#
Multi-Cam	90-0063
Coyote	90-0077















SPECIFICATIONS

H 15 in. x W 10 in. x D 7.5 in. Weight: 10 lb 14 oz

Water Rescue BELT THROW BAG

- Multipurpose water rescue solution designed for easy deployment and retrieval for rescuers
- Quick release reverse tapered throw bag for rapid packing of Triton rope and a slim profile for users



BELT THROW

ITEM#

90-0065

SPECIFICATIONS

Rope Bag: L 11 in. x W 4.5 in. x D 4 in. Gear Pouch: L 5.5 in. x W 5 in. x D 2 in. Total Weight: 2 lb 2 oz

50 Foot THROW BAG

- Must-have necessity for water rescue
- Durable, easy-to-use rescue device with highly visible orange nylon bag and mesh panels that drain quickly for use in multiple rescue operations
- Packed with 50 ft of 7.5 mm Triton Pegasus Rope® rated at 3.050 lbf
- Triton rope is strong enough for multiple uses: as a throw line, tension diagonal, tow line and for pinned kayak recovery





50 FT THROW BAG

ITEM#

90-0064

SPECIFICATIONS

Rope Bag: L 11 in. x W 4.5 in. x D 4 in. Gear Pouch: L 5.5 in. x W 5 in. x D 2 in. Total Weight: 2 lb 2 oz



HELO TAGLINE KIT

- Help control litter movement during hoist operations a requirement on all Army MEDEVAC helicopters
- Contains 250 ft of 7.5 mm Triton Pegasus Rope® featuring a highly visible polypropylene sheath with an Innegra[™] core











The Adaptable Vertical Response Rack was designed meet and exceed full-spectrum operational vertical requirements for urban through mountain / canyon access, rescue, and evacuation variables. The kit capabilities were driven utilizing NFPA 1006 (Operations & Technician Level) JPR's and the KPP's (Key Performance Parameters) within the USSOCOM SOF CASEVAC Program framework. The A-VRR Kit's organic versatility allows it to be scaled up to meet all relevant job performance requirements (JPR's) for complex team access and/or rescue or be scaled down for a

P-A-C-E methodology was engineered throughout the kit design, enabling multiple options for progressive capture and rope grab capabilities depending on the environmental pathology.

- Haul / Mechanical Advantage Capability Multiple simple, compound and complex configurations ranging from a 2:1 through 9:1
- Full spectrum Friction Configurations for Rappels, Bail-outs, And One or Two Person Lowering Systems
- Traditional and MA-assisted Ascending Techniques
- Multiple Basic and Anchor Configurations to Include Contingency and Releasable Systems
- Load Transfer Capabilities to include Knot Passing on Rappel and within a system and casualty Pick-Offs

1 Pulley with Rope Grab Carabiner_{wPulley}







1 8mm














Packaged: 12 in. x 15 in. Weight: 3.25 lb







A-VRR ITEM#

90-5000



R.A.C.E. SUB-T Rapid Access

Casualty Extraction Subterranean Kit

Specifically designed for full-spectrum vertical operations within harsh environments

- Ultra-Lightweight Components
- Hydrophobic Software
- Innate requisite variety (redundancy)
- Easy sustainment
- HRO (High Reliability Organizing) Engineered
- Scalable based on mission context, organic assets, & end-user
- Resilient to environmental pathology
- Capability focused up to Technician Level JPR's (Job Performance Requirements)
- Meets and exceeds operational requirements of SOF CASEVAC Program

Subterranean Deep Hydro | Cave Mountain | Rough Terrain Access | Evac Class B & C Canyon | Slot Mountain | Jungle

R.A.C.E. SUB-T CONTENTS:

















R.A.C.E. SUB-T

ITEM#

90-5002



















FLIGHT RISK KIT

- Ultra-lightweight, compact rescue, and emergency egress / bail-out modular system that easily fits within pocket, bag, or harness
- Auto-Locking F3 Descent device engineered with a double-action lever to give a brake assist in the event of a panic grab
- Gives the end-user a rapid, smooth, and safe egress capability
- Counter-balance casualty evacuation capability

FLIGHT RISK KIT

ITEM#

90-5003

FLIGHT RISK KIT CONTENTS:















SPECIFICATIONS H 6.5 in. x W 6 in. x D 3.5 in. Weight: 3.75 lb



For Multiple Operational Mission Sets:

- Austere | Remote Vehicle Recovery
- Heavy Lifting in Confined Space
- Subterranean
- Tunnel | Mine Environments

1 Squad

Nylon Bag-ODG

1 Pulley

ith Rope Grab

SMAS-H KIT Scalable Mechanical Advantage System-HEAVY

- Scalable Mechanical Advantage System designed for specialized heavy rigging situations constrained by limited amount of personnel
- Simple, rapid, and adaptable rigging configurations from a 3:1 up to a 23:1
- Developed and tested within Mobility, Vehicle Recovery, Subterranean, Confined Space, CBRN, Recovery, and Large Animal movement, stabilization, and rescue environments

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in. Weight: 5.7 lb

SMAS-H KIT

ITEM#

90-5004







3 HMS Auto-lock

Carabiner







1 H-Frame Auto-lock

Carabinerw/Pulley




2 Ascender

/Rope Clamp







- Simple to use yet expandable to the level of expertise of the rescuer and the complexity of the rescue
- Easily carried or packed/attached onto an aid bag, mountain pack, cargo pocket, harness gear loop, or fast-rope bag
- Traditional hauling operations (3:1, 4:1, 5:1, and 9:1)

IMAS KIT Individual Mechanical Advantage System

SPECIFICATIONS

H 7 in. x W 4.5 in. x D 2.75 in. Weight: 4.25 lb

IMAS KIT ITEM#

90-5005

IMAS CONTENTS:























High Threat Extraction | RESCUE



HIGH THREAT EXTRACTION KIT



that Care Under Fire casualty treatment may be performed as soon as possible

• Facilitates immediate casualty extraction so

- Self-contained kit with all of the essential tools for tactical operators and rescue personnel
- Includes a C-A-T® tourniquet to immediately address extremity hemorrhage at the point of wounding
- Rugged nylon bag attaches to the duty belt as a drop down leg rig or to the assault vest/body armor with PALS/MOLLE-style connectors

















HIGH THREAT EXTRACTION KIT

COLOR	ITEM#
Black	80-0139
OD Green	80-0140
Coyote	80-0141
Digital Camo	80-0142

CUE Rescue Gear

IRUPU HARNESS



UIAA: Yes Certification: EN 12277 Type C Waist: 65-110 cm | Leg Loops Max: 70 cm Weight: 21.16 oz (600 gm)

- Multi-disciplinary rugged harness for use in operational environments (Mountain evac, Sub-T, Canyons, & Caves)
- 2 robust, floating gear loops
- Enhanced stitching and replaceable seat provide uniform protection in environments with abrasion hazards
- Stitched clip-in buckles with tension straps for easy on & off
- Left/ right leg loop markings make it easier to put the harness on correctly
- Large, textile tie-in point

IRUPU HARNESS

ITEM# 90-4020

The first & only foldable climbing helmet on the market

EDELRID

SPECS

CE marking: CE 0123 Weight: 13.76 oz (390 gm)

MADILLO HELMET

ITEM#

90-4019

HELMET

 Robust and impactresistant injectionmoulded ABS shell

 Innovative combination of 3 different foam types (EPS, EPP and EVA) deliver optimum safety and comfort

- Fully adjustable chinstrap for optimal fit
- 4 robust headlamp clips
- Detachable, washable padding
- Side/front/back impact tested to EN 12492

RAT STRAP

Rescue Assault Tether

• For head/feet first linear drags during shield-assisted, openfield/indoor-corridor rescues or single-rescuer extractions

• Innovative Daisy Chain configuration offers multiple reinforced attachment points

SPECIFICATIONS

Pkg: H 1 in. x W 4 in. x D 5.5 in. Extended: L 62 in. Minimum Break Strength: 2,702 lb Weight: 12.7 oz

RAT STRAP

ITEM#	NSN#
90-0002	6530-01-632-7872

Under the Vest Harness

high threat extraction

• Built-in handle system for immediate

• Operator harness to aid in the rescue

of downed personnel without the

need to holster or sling weapons

casualty's ballistic vest, reducing

exposed, unprotected body regions

• Prevents the "ride up" of the

Height: 20 in. from waist to non-deployed handle

Length: 14 in. from harness

attachment point to top of handle

SPECIFICATIONS

Weight: 1 lb 1.6 oz. Waist Size: 30 in. to 45 in.



Point of Wounding Extraction

- Point of wounding extraction
- Back pack carry
- Hasty harness for traditional and confined space lowering and hauling

SPECIFICATIONS

Pkg: H 18 in. x W 3.5 in. x D 1 in. Extended: L 58 in. Weight: 9.3 oz

EX-TRACT STRAP

ITEM#

90-0082

DRAGON DRAGON HANDLE

ESCUE High Threat Extraction

Tool offering head/feet first drags

- Maintain an anatomically correct posture for rapidly extracting greater weight
- Decreases manpower needs in 'zeroed in" locations
- Black Omega, non-locking "D" carabiner rated for 31kN (6.968.8 lb).

SPECIFICATIONS

H 6.3 in. x W 4.1 in. x D 2.25 in. Weight: 11.3 oz Device Length: L 42 in. US Patent(s): Pending

DRAGON HARNESS

ITEM#

90-0003

DRAGON HANDLE

ITEM#	NSN#
90-0004	6530-01-598-9595

HASTY HARNESS Multi-purpose Harness

- Indispensable gear for every webbing equipment toolbox
- Resourceful, multi-purpose harness for unlimited applications & maximum versatility
- Constructed of 22 ft loop of 1 in. Tech Tape (MBS: 4,400 lb) with opposing color thread to aid in visual safety inspection

SPECIFICATIONS

Loop: L 10 ft 11 in. x W 1 in. Circumference - 21 ft 10 in. Min. Break Str: Tech Tape (4,400 lb) Sewn Sling (6,000 lb)





TIBLOC

Ascender/Rope Grab by Petzl®

Ultra-light Compact Ascender/Progress Capture Device for rope ascents, hauling systems or friction knot replacement in selfrescue situations



Rope Size: 8-11 mm Weight: 1.23 oz (35 g)

ITEM#

ZZ-0547



ROPE GRAB PULLEY-SPOC ROLLCLIP PULLEYCarabiner

Ultra-light pulley designed for lifting loads, crevasse rescue operations, or as an emergency ascender

SPECS

Rope Size: 7-11 mm Weight: 2.12 oz (60 g) Diameter Pulley: 20 mm

ITEM#

90-4005

Pulley-carabiner with a non-locking gate opening

SPECS

Rope.: 7 to 13 mm Major axis: 20 kN Minor axis: 8 kN Open gate: 7 kN 5.25 in. x 3 in. x 1.25 in. Weight: 3.76 oz

ITEM#

90-4006



PIRATE CARABINER

by Rock Exotica

Ideal as a belay biner at the harness, or as a connector in a system

SPECS

Gate Opening: 25 mm MBS Major: 26 kN MBS Minor: 12 kN MBS Open: 7 kN H 107 mm x W 71 mm Weight: 3.09 oz (87 g)

ITEM#

90-4003

ROCKO AUTO LOCK

Carabiner by Rock Exotica

3-stage Autolocking carabiner with intuitive pull-down and twist action

SPECS

Gate Opening: 24 mm MBS Major: 24 kN MBS Minor: 11 kN MBS Open: 6 kN H 107 mm x W 71 mm Weight: 2.88 oz (81.65 g)

ITEM#


90-4001



MBS Major: 24 kN MBS Minor: 9 kN MBS Open: 6 kN H 107 mm x W 71 mm Weight: 2.96 oz (84 g)

ITEM#

90-4002



ROCKO NON LOCKING

Carabiner by Rock Exotica® Non locking enhanced Oval biner with Snagresistant keylock nose design

SPECS

Gate Opening: 25 mm MBS Major: 22 kN MBS Minor: 10 kN MBS Open: 6 kN H 107 mm x W 71 mm Weight: 2.36 oz (67 g)

ITEM#

90-4022



OMEGALITE 4 WIREGATE

ISO Cold Forged design

SPECS

Gate Opening: 1.10 in. (28 mm) MBS Major: 26 kN MBS Minor: 9 kN MBS Open: 8 kN L 4.25 in. x W 2.75 in. Weight: 0.12 lb

ITEM#

ZZ-0314



QUICK LINK SCREW LOCK

5mm Oval Stainless Steel Quick Link Carabiner

SPECS

Rod stock Diameter: 5 mm Gate opening: 6 mm Major Axis: 25 kN H 1.97 in. x W 0.98 in. Weight: 0.7 oz

ITEM#

90-4015



QUICK LINK SCREW LOCK

7mm Oval Stainless Steel Quick Link Carabiner

SPECS

Rod stock Diameter: 7 mm Gate opening: 16 mm Major Axis: 25 kN Minor Axis: 10 kN H 3.25 in. x W 1.25 in. Weight: 2.4 oz

ITEM#

ZZ-0544



ALPHA FIX Y SLING

30/60 cm by Edelrid® Y-Shaped self belay sling for edge restraint, DCD extender, knot passing on rappel, edge transition safety, & hasty pick-offs

SPECS

Load capacity: 23 Kn CE 0123 Weight: 4.16 oz (118 gm)

ITEM#

90-4021



ARAMID CORE BLUE SLING

120 cm by Edelrid

High-strength sewn Aramid sling ideal as a full-spectrum anchor, harness extension, alternate prusik, & canyon quickdraw for Rig to Raise conversions

SPECS

Load capacity: 22 Kn CF 0123 Width: 0.24 in. (6 mm) Weight: 0.78 oz (22 gm)

ITEM#

90-4011



BASIC ASCENDER

bv Petzl^o

Compact ascender with excellent grip for rope ascent on a fixed line or as progress capture in a hauling system

SPECS

Rope Size: 8-11 mm Weight: 3 oz (85 g)

ITEM#

ZZ-0554



SPECS

bv Petzl®

Rope Size: 8-11 mm Weight: 6 oz (170 g)

REVERSO 4

amount of braking

Belay/Rappel device by Petzl

braking capability

Compact, ultra-light

belay device allowing excellent control during

GRIGRI 2 BELAY DEVICE

ITEM#

ZZ-0545

Regulates the

friction on the

rope according

to diameter and

condition

SPECS

ITEM#



MINI MACHINED PULLEY

by Rock Exotica

Mini pulley with recessed axle screw, machined sideplates & integrated axle for a stronger, compact, snag-free profile

SPECS

Rope Size: up to 11 mm Weight: 2.9 oz (83 g)

ITEM#

ZZ-0548

TRIRIG RIGGING PLATE

by Rock Exotical

Lightweight. high performance 3 hole rigplate with numbered holes to assist with detailed rigging

SPECS


Break Strg: 36 kN Weight: 1.8 oz (51 g)

ITEM#

ZZ-0553

FOOTAPE

Adjustable Foot Loop by Petzl®



SPECS

Weight: 2.3 oz (65 g)

ITEM#

ZZ-0551



77-0546

ROCKD CARABINER by Rock Exotical

Weight: 2.08 oz (59 g)

Narrow asymmetric "D" shape machined frame with snag resistant keylock nose and large, ergonomic angled gate opening

SPECS

Gate Opening: .24 mm Weight: 2.6 oz (73 g)

ITEM#

ZZ-0543



120 cm

ANNEAU SLINGS

by Petzl®

Polyester sling to set up a belay or an anchor

SPECS

W 19 mm L 60 or 120 cm Wgt: 60 cm - 2 oz 120 cm - 3.5 oz





Rescue Knives



- Perfectly suitable to shorten static lines or ropes
- Large swivel hole for attaching to harness or backpack with carabiner
- Sharp, locking blade (partly serrated)
- High-quality Stainless Steel

SPECIFICATIONS

L 3 in. x W 1 in. (closed blade) Blade Length: 3 in. Weight: 1.69 oz (48 gm)

ROPE TOOTH KNIFE

ITEM#







HOOK SAFETY CUTTERS by Benchmade®

- Hook cutters designed to quickly and reliably cut through a wide array of materials
- Extensively tested and widely used by the US Military for medical and egress purposes

ITEM#	STYLE	NSN#
ZZ-0652	7 Hook	
ZZ-0250	8 Hook	

SPECS

7 Hook: L 4.3 in. 8 Hook: 6.5 in. Weight: 1.5 oz

SOCP RESCUE TOOL by Benchmade®

- Fixed 440C high-chromium stainless steel blade hook for rescue performance and less-than-lethal self-defense
- Includes carbide glass breaker, oxygen wrench, single-piece handle with finger loop & molded sheath
- Designed by Special Operations Combatives Program (SOCP)

_	_	_	
		6	CS
	г		CO

Blade: 2.5 in. Overall: 6.75 in. Weight: 1.88 oz

ITEM #	NSN#
ZZ-0760	



NAVY SAR SOML 112 SBK by Benchmade®

Approved for Navy SAR -" So Others May Live"

- N680 ultra stainless steel blade with opposing fixed bevel blade style with blunt tip, integrated rescue hook and carbide glass breaker
- Boltaron® sheath designed to fit current Navy SAR kits

Blade: 3.5 in. Overall: 8.02 in. Weight: 3.2 oz

ITEM#	NSN#
90-0078	1095-01-677-1013

TRIAGE 917 BK bv Benchmade®

- Black drop point style CPM-S30V steel blade with a smooth edge
- 440C folding hook & carbide glass breaker

ITEM#	NSN#
90-0092	

SPECS

Blade: 3.5 in Overall: 8.3 in. Weight: 5.24 oz

TRIAGE 916 SBK bv Benchmade®

- Dynamic triple utility tool featuring three essential rescue tools in a single compact design
- N680 high corrosion resistant, partially serrated, opposing bevel blade, 440C folding

hook & carbide glass break

SPECS

Blade: 3.5 in Overall: 8.2 in. Weight: 5.1 oz

ITEM#	HANDLE	NSN#
90-0096	Orange	
90-0060	Black	



Leatherman® RAPTOR" MULTI-TOOL

High-performance multi-tool equipped with a variety of mission critical tools essential to emergency medical practitioners in the field

- 1. 420HC Stainless Steel Folding Medical Shears
- 2. Strap Cutter
- 3. Oxygen Tank Wrench
- 4. Rina Cutter
- 5. Carbide Glass Breaker
- 6.5 cm Ruler
- 7. Holster with Clip

HANDLE	ITEM#	NSN#
Black	70-0014	5110-01-626-8766
Blue	70-0264	
Coyote	70-0265	
Red	70-0266	

See our Treatment Support Section on page 87 for the compatible Scissor Leash

SUPERTOOL 300

by Leatherman®

Heavy-duty, all purpose multi-tool featuring 19 devices in one easy to carry configuration

SPECIFICATIONS

Closed: L 4.5 in. Primary Blade: L 3.2 in.. Weight: 9.6 oz

SUPERTOOL 300

ITEM#

ZZ-0542

- 1. Needlenose Pliers
- 2. Regular Pliers
- 3. Replaceable Wire Cutters
- 4. Hard-Wire Replaceable
- 5. Stranded-Wire Cutters
- 6. Electrical Crimper
- 7. Wire Stripper
- 8.420HC Knife
- 9. 420HC Serrated Knife



- 11. Awl with Thread loop
- 12.9 in. Ruler
- 13. Can Opener
- 14. Bottle Opener
- 15. Wood/Metal File
- 16. Phillips Screwdriver
- 17. Large Screwdriver
- 18. Medium Screwdriver
- 19. Small Screwdriver



MUT® by Leatherman®

Tactical and practical multi-tool for military, law enforcement and civilian shooters - features 16 devices

- 1. Needlenose Pliers
- 2. Regular Pliers
- 3. Replaceable Wire Cutters
- 4. Hard-Wire Replac. Cutters
- 5. Stranded-Wire Cutters
- 6. Electrical Crimper
- 7.420HC Knife
- 8. Saw
- 9. Replaceable Cutting Hook

- 10. Hammer
- 11. Bolt Override Tool
- 12. Bronze Carbon Scraper
- 13. #8-32 Cleaning Rod/ Brush Adapter
- 14. Firearm Disassembly
- 15. Carabiner/Bottle Opener Punch
- 16. Large Bit Drive

SPECIFICATIONS

Closed: L 5 in. Primary Blade: L 3 in. Weight: 11.2 oz

MUT MULTITOOL

ITEM#

ZZ-0539



TACTIKKA HEADLAMP

ITEM#

ZZ-0552

SPECIFICATIONS

Light Output: Min. 5 Lumens Max: 120 Lumens (160 Lumens on Boost mode) Weight: (with batteries) 3 oz

 ϵ

TACTIKKA + RGB

Headlamp by Petzl®

- Ambient Mode: Wide beam with reduced brightness, avoids blinding others in a group situation and provides a long burn time
- Close Range Vision Mode Mixed beam with wide, uniform flood or 3 colored LED lighting modes (red, green & blue) to remain discreet
- Boost Mode for temporary access to maximum power, 160 lumens (60 m)
- Strobe mode for signaling one's presence, accessible at any time
- Reserve Mode- Burn time management automatically switches on by signaling with blinking light when batteries are almost drained. Red lighting activates to preserve a minimal amount of usage

NAR HEADLAMP

- 100 Lumen LED Light
- Dual color (white and red)
- High and low lighting modes (White)
- Sturdy and simple to rotate single arm bracket
- Large push button for operation
- Splash and water resistant (IPX4)
- Comes with 2 AAA batteries
- 58 hour burn time

NAR HEADLAMP

ITEM#	NSN #	
ZZ-0106	6515-01-542-7937	

SPECIFICATIONS

Size: W 2.25 in. x D 1.5 in. Weight: 2.25 oz





Streamlight HEADLAMP

- Ring Sealed construction -Waterproof to 1M for 1 Hr
- Up to 6 hrs, runtime on high up to 24 hrs. on low
- Impact resistant ABS housing
- Includes (2) 1.5V "AAA" batteries
- Solid State Power Regulation provides maximum light output throughout battery life
- Shatter resistant polycarbonate lens

SPECIFICATIONS

Faceshield Diameter: 1.3 in. Maximum headlight: W 2.25 in. x D 1.6 in. Weight: (with batteries & headband) 2.75 oz



ITEM#

ZZ-1090

HEADLAMP

GECKO GRIP®

Multi-purpose Tape

Laminated polyethylene cloth tape featuring an aggressive natural rubber based adhesive formulation and perfect for heavy-duty applications

SPECIFICATIONS

Pack of 6: L 12.25 in. x W 0.9 in. x D 0.7 in. Individual Roll: W 2 in. x L 5 yd Weight: Roll: 2.7 oz Pack: 1 lb 8.7 oz

GECKO GRIP

UNIT	ITEM#
6 Pack	80-0131



RESCUE | Tactical Lights

TOMAHAWK NIGHT VISION

by First Light

- Right angle tactical light with rotating head and 3 button control panel for one-handed thumb actuation
- Provides White, IR, and Red+Green LED mix for highlighting blood during casualty treatment
- Featuring Integrated high intensity IR signaling beacon
- Includes momentary and three intensity levels of constant-on (all output colors)
- Can be used handheld or hands-free with MOLLE TRS Blade Mount

TOMAHAWK NIGHT VISION

CASE COLOR	TASK LIGHT	ITEM#
Coyote	UV Task	70-0015
OD Green	UV Task	70-0112
Coyote	Blue Task	70-0111
OD Green	Blue Task	70-0110

SPECIFICATIONS

Output Wavelengths: White, IR, & Red+Green (Blood Highlighting), IR Signaling Beacon

Maximum White Light Output: 425 lumens Max Output Runtime: 2 hours (ANSI Tested)

Power Source: CR123 Batteries (2) | Material: Aerospace-Grade Aluminum With mount: H 3.4 in. x L 2.4 in. x W 1.5 in. | Weight: (with mount & batteries) 8 oz



TOMAHAWK MULTI-COLOR

by First Light®

- Right angle tactical light with rotating head and 3 button control panel
- Provides White, Red. Green, and Red+Green LED mix for highlighting blood during casualty treatment
- Includes momentary and three intensity levels of constant-on

SPECIFICATIONS

Output Wavelengths: White, Red, Green, & Red+Green (Blood Highlighting) Maximum White Light: 425 lumens Max Output Runtime: 2 hours (ANSI Tested) Power Source: CR123 Batteries (2) Material: Aerospace-Grade Aluminum With mount: H 3.4 in. x L 2.4 in. x W 1.5 in. Weight (with mount & batteries): 7.2 oz

TOMAHAWK MULTI-COLOR

CASE COLOR	TASK LIGHT	ITEM#
Coyote	None	70-0016
OD Green	None	70-0117
Coyote	UV Task	70-0116
OD Green	UV Task	70-0115
Coyote	Blue Task	70-0114
OD Green	Blue Task	70-0113

TORQ by First Light®

- Right angle tactical light, rotating head and 3 button control
- White, Red. Green, and Red+Green LED mix for highlighting blood
- Includes momentary and three intensity levels of constant-on (all output colors)
- Features integrated high intensity white tactical strobe and red/white/green flashing safety beacon
- 10 hour Runtime on high output

TORQ

ITEM#	NSN#	
70-0118	6515-01-646-2636	

SPECIFICATIONS

Wavelengths: White, Red, Green, & Red+Green (Blood Highlighting) Max Runtime: 10 hours (ANSI Tested) With mount: H 3.9 in. x L 2.6 in. x W 1.3 in. Weight (with mount & batteries): 7.1 oz









TREATMENT SUPPORT | Hyperthermia

POLAR SKIN™

- The latest cooling technology for effective Human Thermo-regulation
- Polar Skin Mini-Mesh,™ a unique smart textile with a hydrofunctional polymer that allows for liquid absorptive capacity 350% greater than base weight
- ArcticAqua,™ a proprietary solution that prevents rigidity at temperatures well below 0°F while providing the requisite moisture to support evaporative, convective and conductive cooling
- Adheres to current U.S. Army and National Athletic Trainers' Association protocols for heat illness treatment
- Protracted direct cooling of high blood flow regions

COMPREHENSIVE COOLING SYSTEM

- Highly-effective portable cooling system for the treatment of heat-related illnesses
- Scientifically proven to reduce core temperature by .09°F/ minute in 100°F heat when used according to protocol

COMPREHENSIVE COOLING

ITEM# ZZ-0674

SPECIFICATIONS

L 11 in. x W 10 in. x H 4 in. Weight: 8 lb 5.3 oz



COMPLETE CARE PACK

- Portable, highly insulated option
- Contains one Polar Skin Comprehensive Cooling System and one heavy-duty, leak & sweat proof cooler that keeps ice from melting for 24 hours in +100°F heat

COMPLETE CARE PACK

ITEM# 80-0451

SPECIFICATIONS

W 14 in. x L 12 in. x D 7 in. Weight: 10 lb 3.2 oz

COMPLETE CARE BACKPACK

Available in either Single Patient (one comprehensive cooling system) or Double Patient (two cooling systems) configurations

COMPLETE CARE BACKPACK

SIZE	ITEM #
Single	80-0543
Double	80-0544

SPECIFICATIONS

W 19.5 in. x L 15 in. x D 7.5 in. Single Weight: 12 lb 13 oz Double Weight: 21 lb 13 oz



Available in Standard, Large and Backpack COOLER

Standard size accommodates one Polar Skin™ Comprehensive Cooling System and Large or Backpack Size accommodates two cooling systems

POLAR SKIN COOLERS

COLOR	SIZE	ITEM#
Red	Standard	ZZ-0682
Black	Standard	ZZ-0759
Black	Large	ZZ-0756
Black	Backpack	ZZ-0794



Standard: L 14 in. x W 12 in. x H 7 in. Weight: 1 lb 3.2 oz Large: L 18.5 in. x W 13 in. x H 9.5 in Weight: 2 lb 13 oz Backpack: L 19.5 in. x W 15 in. x D 7.5 in. Weight: 3 lb 13 oz

ICE SHEETS Standard

- Principle cooling component of Polar Skin™ cooling systems
- Adheres to current National Athletic Trainers' Assoc and U.S. Army protocols for heat illness treatment





STANDARD ICE SHEETS

UNIT	ITEM#	NSN#
Single Pack	ZZ-0668	6532-01-653-6054
Double Pack	ZZ-0669	
Triple Pack	ZZ-0755	

SPECIFICATIONS

Single: L 4.5 in. x W 7.75 in. x H 0.5 in. Wgt: 9.6 oz Double: L 9.5 in. x W 8 in. x H 0.5 in. Wgt: 1 lb 6.4 oz Triple: L 12.5 in. x W 7.5 in. x H 0.5 in. Wgt: 2 lb 3.2 oz

ICE SHEETS Large

Same characteristics as the standard size Ice Sheets, but in a larger 36 in. x 60 in. size

LARGE ICE SHEETS

ITEM#

ZZ-0758

SPECIFICATIONS

L 8.5 in. x W 5.5 in. x H 1.5 in. Weight: 1 lb 8 oz

C2E ICE PACK BUNDLE

- Four Comfort Gel Packs
- Facilitates Conductive, Convective, and Evaporative (C2E) cooling of high blood flow (large artery) regions of the body

C2E ICE PACK

ITEM#

ZZ-0675

SPECIFICATIONS

Indiv. Gel Pk: W 6 in. x L 10 in. Pkg: L 8 in. x W 9.5 in. x H 1.25 in. Weight: 4 lb 4 oz

COOLING COLLAR



- Two gel packs encased in integrated pouches
- Designed to align with the carotid artery

COOLING COLLAR

ITEM#

ZZ-0670

SPECIFICATIONS

Indiv. Gel Pk: W 6 in. x L 10 in. Pkg: L 8 in. x W 9.5 in. x H 1.3 in. Weight: 4 lb 4 oz



- Multi-capable device for splinting a wide variety of major and minor fractures
- Lightweight and compact for easy field carry with 6 easily adjustable interconnected telescoping poles and 4 wide elastic bands with hook & loop securing
- Includes clearly marked directions on the device as well as an easy-to-use, step-by-step instructional insert



EATMENT SUPPORT | Splinting

TACTICAL TRACTION SPLINT

(TTS) Compact Mechanical Traction Device

SPECIFICATIONS Packaged: L 11 in. x W 7 in. x D 3 in.

Weight: 1 lb 5 oz

TACTICAL TRACTION SPLINT

COLOR	ITEM #	NSN#
Black	50-0015	6515-01-346-9186
Coyote	50-0033	
Digital Camo	50-0034	
OD Green	50-0032	

• Relieve pain, increase patient comfort and reduce the risk of vascular injury or nerve damage in femur fractures

 Lightweight and portable kit with clearly marked elastic straps for the proper fit to the thigh, knee & shin

TACTICAL TRACTION SPLINT (TTS):













SAM® SPLINT II

- Versatile, malleable splint easily molds to fit most fractures
- Lightweight, comfortable & strong

• Waterproof, Reusable

SPECIFICATIONS

CE

H 9.25 in. x W 4.25 in. x D 1 in. Unfolded: L 36 in. x W 4.25 in. Weight: 4.5 oz

SAM SPLINT

ITEM#	NSN#
50-1005	6515-01-494-1951

SAM® PELVIC SLING for Pelvic Ring Fractures

- Force-controlled, ideal compression every time without the guesswork
- Rapid, simple, three-step application, usually in less than one minute



SAM PELVIC SLING

ITEM#	NSN#
50-0004	6515-01-509-6866

SPECIFICATIONS

H 7.25 in. x W 7.75 in. x D 4 in. Hip Circumference: 32 in. to 50 in. Weight: 9.8 oz

SAM[®] SOFT SHELL SPLINT

SPECIFICATIONS

L 15 in. x W 4.25 in. x D 0.5 in. Weight: 3.2 oz

SAM SOFT SHELL SPLINT

UNIT	ITEM#	NSN#
Each	50-0008	6515-01-523-5440
Box of 20	50-0009	6515-01-507-8074

TRAUMA SHEARS

SPECIFICATIONS

SMALL: L 6.25 in. x W 2.75 in. x D 0.3 in.

Weight: 1.5 oz

LARGE: L 7.25 in. x W 3.7 in. x D 0.4 in.

Weight: 2 oz

TRAUMA SHEARS

* NATO NSN

SIZE	ITEM#	NSN#
S 6.25 in.	ZZ-0064	6510-01-622-1624*
L 7.25 in.	ZZ-0063	6515-01-538-9276



SCISSOR LEASH

SPECIFICATIONS

H 6 in. x W 1.6 in. x D 0.6 in. Cable Length - 36 in. Weight: 1.6 oz

SCISSOR LEASH

ITEM#	NSN#
70-0005	6515-01-540-7226

6 oz MOLLE Retractable

GEAR LEASH

- MOLLE Mount compatible / 360° rotating
- High break strength, salt water resistant self-flushing
- 36 in. Extension / Spectra Nylon Line

CEADIEACH

GEAR LEASH	
ITEM#	
70-0155	



SPECIFICATIONS

Retracted: H 3.5 in. x W 1.6 in. x D 1 in. Cable Length - 36 in. Weight: 2.0 oz

BEAR CLAW



- Form-fitting, ambidextrous design with textured fingertips
- Thick palm ensures durability
- Textured for improved wet and dry gripping
- · Longer and thicker than standard nitrile gloves



Nitrile gloves available in Black or Tan colors and come

either as pre-packaged pairs (kits) or in bulk (box)



ITEM#

ZZ-0218

70-0001

ZZ-0217

70-0002

ZZ-0216

70-0003

77-0215

Weight: 1 lb 1 oz Box: L 11 in. x W 5 in. x D 3.75 in. Weight: 2 lb 4 oz

NSN#

6515-01-515-0195

6515-01-672-7977*

6515-01-515-0200

6515-01-632-0059

6515-01-515-0197

6515-01-672-7985*

6515-01-515-0190

BEAR	BEAR CLAW GLOVES CE *NATO NSN		
SIZE	PKG QTY	ITEM#	NSN#
S	Kit 25 pair	ZZ-0246	6515-01-619-5350
S	Box 50 pair	70-0009	6515-01-600-6059*
М	Kit 25 pair	ZZ-0247	6515-01-619-5689
М	Box 50 pair	70-0010	6515-01-600-6057
L	Kit 25 pair	ZZ-0248	6515-01-600-0262
L	Box 50 pair	70-0011	6515-01-600-6039
XL	Kit 25 pair	ZZ-0249	6515-01-619-5367

NAR RESPIRATOR

- NIOSH N95 certified particulate barrier
- Adjustable nose bridge for a secure, comfortable fit

70-0012

SPECIFICATIONS

Packaged: L 5 in. x W 4.25 in. x H 0.75 in. Mask size: Adult Weight: 0.9 oz

Box | 45 pair

NAR RESPIRATOR

ITEM #	NSN #
70-0006	6515-01-537-1531



BLACK TALON GLOVES

Kit | 25 pair

Box | 50 pair

Kit | 25 pair

Box | 50 pair

Kit | 25 pair

Box | 50 pair

Kit | 25 pair

SIZE PKG | QTY

S

М

L

L

ΧI

1 N95 certified Face Mask & 2 pair Bear Claw Nitrile Gloves

SPECIFICATIONS

L 3.5 in. x W 4 in. x H 1.5 in. Weight: 2 oz

ITEM#	NSN#
70-0007	6515-01-537-1534

NAR PPE KIT

ITEM #	NSN #
70-0007	6515-01-537-1534



6515-01-600-6052*

Designed to maximize personal protection when handling and managing detainees, refugees, etc.

H 6 in. x W 8 in. x D 1 in. Weight: 4.6 oz

SPECIFICATIONS

BIOLOGICAL PERSONAL PROTECTION KIT:





















BIOLOGICAL PPE KI	
ITEM#	NSN#
80-0488	6515-01-388-1351

C€ *NATO NSN



REATMENT SUPPORT | Eye Care & Protection

PES Polycarbonate Eye Shield

- Meets CoTCCC treatment protocols for penetrating eye trauma
- Convex eye shield accommodating either eye
- Translucent orange shell features perforations to allow air circulation



POLYCARBONATE EYE SHIELD

UNIT	ITEM #	NSN#
Single	30-0142	6515-01-647-5456
Box of 100	30-0143	6515-01-630-5551

NAR **IPRO**® Tactical Goggle System

High performance eye protection with Hexo-vents for maximum ventilation, fine particle filtration and fog-free viewing

SPECIFICATIONS

Packaged: L 2.6 in. x W 7 in. x D 3 in. Weight: 8 oz

NAR IPRO

ITEM#

ZZ-0109





Glasses: W 5.5 in. x H 2 in. x D 4.75 in. Case: L 6.5 in. x W 4.3 in. x H 3 in. Weight: 9 oz


 Superior ballistic impact protection per US MIL-PRF-31013 standard clause 3.5.1.1 and 100% protection from UVA/UVB rays

Includes clear & gray interchangeable lenses

AEGIS ECHO II

ITEM#

90-0067



SPECIFICATIONS

3 in. x 2.38 in. Weight: 0.2 oz

FOX EYE SHIELD

OK ETE STILLED			
STYLE	UNIT	ITEM#	NSN #
No Garter	Box of 12	30-0144	6515-01-449-1016
Garter	Single	30-0055	
Garter	Box of 50	30-5550	6515-01-598-1877

OTW GOGGLES Outside the Wire - with Case

• 100% UVA/B/C protection with maximum ballistic protection

• Tapered, neutral tint gray lenses with anti-fog and scratch resistant coating, reduced glare and correct distortion

OTW GOGGLES

ITEM#

90-0066

SPECIFICATIONS

Goggles: W 7.5 in. x H 2.75 in. Case: L 8.75 in. x W 3.5 in. x H 3 in. Weight: 10.1 oz

GRAY MAN

EyePro

Low visibility EyePro option with classic, subtle styling

GRAY MAN EYEPRO

ITEM#

90-0068



SPECIFICATIONS

Glasses: W 5. in. x H 2 in. x D 5 in. Case: L 6.625 in. x W 3 in. x H 2.5 in. Weight: 9 oz



BURNTEC® Hydrogel Dressings

- Sterile non-stick hydrogel dressings for treating 1st, 2nd & 3rd-degree burns (thermal, chemical, electric)
- Also treats other skin injuries, such as sunburns, abrasions, skin irritations, bruises and swelling (insect bites)
- Applied dressings provide a cool, soothing effect for up to 24 hours
- Easy to handle without spilling or damaging the hydrogel
- Permeable to water vapor and oxygen but impermeable to bacteria
- Biocompatible and hypo-allergenic; Reduces risk of hypertrophic scars

BURNTEC DRESSI	C€	
SIZE	ITEM#	NSN#
2 in. x 2 in.	30-0114	6510-01-631-1913
2.5 in. x 5 in.	30-0115	6510-01-631-2080
4 in. x 4 in.	30-0110	6510-01-631-1883
5 in. x 5 in.	30-0116	6510-01-631-2138
5 in. x 10 in.	30-0117	6510-01-631-1907
8 in. x 8 in.	30-0118	6510-01-631-1897
8.5 in. x 11 in.	30-0119	6510-01-631-1890
4 in. x 16 in.	30-0111	6510-01-631-2124
8 in. x 16 in.	30-0112	6510-01-626-2990
16 in. x 24 in.	30-0120	6510-01-631-1886
10 in. x 10 in. Face	30-0121	6510-01-627-1834
12 in. x 16 in. Face	30-0113	6510-01-627-1855
16 in. x 24 in. Face	30-0122	

DRY STERILE BURN DRESSING

- Dry sterile, non adherent burn dressing & bandage
- Includes two No. 3 closed safety pins and rubber band individually packaged

SPECIFICATIONS

Bandage: 45 in. x 45 in. x 63 in. Pkg: W 4.5 in. x H 4 in. x D 0.5 in. Weight: 1 oz



ITEM#	NSN#
30-0162	6510-01-587-6579



• Convenient array of BurnTec® dressings in a handy carrying case

• Kit contains (1) each of the following dressings:

2 in. x 2 in. 5 in. x 10 in. 8 in. x 16 in. 8 in. x 8 in. 10 in. x 10 in. (face) 2.5 in. x 5 in. 4 in. x 4 in. 4 in. x 16 in. 12 in. x 16 in. (face) 5 in. x 5 in. 8.5 in. x 11 in.

BURNTEC BURN KIT

OPTION	ITEM#	NSN #
Kit	80-0361	6515-01-630-9177
Bag Only	30-0123	6510-01-626-6079



online video available ∧ scan to view

SPECIFICATIONS


L 12 in. x W 15.5 in. x D 3.5 in. Weight: 6 lb 4 oz

BurnTec® MINOR BURN DRESSING KIT

 ϵ

• Variety of BurnTec® dressings in a handy carrying case

• Kit contains the following: Two 4 in. x 4 in. dressings Two 4 in. x 16 in. dressings One 8.5 in. x 11 in dressing



Bag Only 30-0150

SPECIFICATIONS

L 14 in. x W 8.5 in. x D 2.5 in. Weight: 2 lb 13 oz



BurnTe



REATMENT SUPPORT | Triage System

T2® TAG NAR

• Triage & treatment tag utilizing the S.A.L.T. Mass Casualty Triage Standard

• To aid in documenting critical patient information throughout the entire mass casualty response process

• Liquid-resistant in all settings; made of rugged materials contained in a plastic sheath with latex-free band

T2 TAG

UNIT	ITEM#	NSN#
20 Pack	20-0031	6515-01-624-7148
1000 Pack	20-0032	6545-01-592-3261









TACTICAL TRIAGE KIT 12

- · Compact incident command system
- Streamlines medical triage
- Contained in a lightweight, rugged vacuum-sealed pouch

SPECIFICATIONS

L 5.25 x W 7 in. x D 1.25 in. Weight: 4.9 oz

TACTICAL TRIAGE KIT

ITEM #	NSN#
20-0040	6545-01-591-2134

TACTICAL TRIAGE KIT:







COMBAT CASUALTY CARD

- Streamlined triage card that meets all TCCC standards for casualty management
- State-of-the-art ruggedized paper is tear resistant, waterproof & designed to fit into the IFAK
- Includes grommet reinforced punch with latex-free band for attachment
- Currently in use by Special **Operations Forces**

SPECIFICATIONS

H 6 in. x W 4.5 in. each Weight: 0.22 oz ea | 5.5 oz 25-Pack

COMPAT CASHALTY CADD

COMBAI CASCALI I CARD		
UNIT	ITEM#	NSN#
25 Pack	ZZ-0262	6515-01-587-1700

COMBAT CASUALTY REFERENCE CARD

- Rugged reference card for emergency management featuring the current TCCC Guidelines
- Clear, concise instructions for medics and warfighters
- Information divided into four sections:
 - Care Under Fire
 - Tactical Field Care
 - TACEVAC Protocols
- 9-line MEDEVAC request/notes
- Pocket sized for convenience
- Eighteen pages



COMBAT CASUALTY REFERENCE CARD

ITEM#	NSN#
ZZ-0057	6515-01-537-4161

SPECIFICATIONS H 6 in. x W 5 in. x D 0.1 in. Weight: 2.3 oz



SPECIFICATIONS

Pkg: H 10 in. x W 8 in. x D 5 in. Weight: 1 lb 14 oz

- equipped with T2/SALT triage cards
- Includes a laminated S.A.L.T. Triage Flowchart &

TRIAGE OFFICER KIT

COLOR	ITEM#	NSN#
OD Green	20-0023	6545-01-537-4189
Digital Camo	20-0075	
Red	20-0076	

COMMAND KIT 72°

- Organized, rapid response system for Incident Commanders during a mass casualty/tactical triage incident
- 3 laminated, detachable Command Boards
- Pad of 50 Patient-Tracking Boards
- 3-ring, clear storage sheets for bar codes & mesh pockets to store extra accessories

SPECIFICATIONS

Case Closed: L 21 in. x W 17 in. x D 6 in. Weight: 4 lb 14 oz





COLOR	ITEM#	NSN#
OD Green	20-0025	6545-01-537-4197
Digital Camo	20-0074	
Red	20-0026	



- Comprehensive triage system ideal in disaster response, field hospitals & military troop concentrations.
- Lightweight, modular, man-portable triage kit.
- includes user-friendly tools for every aspect of the uniformed S.A.L.T. model for disaster triage

SPECIFICATIONS

H 21 in. x W 17 in. x D 9 in. Weight: 15 lb 2 oz

MCIM KIT

COLOR	ITEM #	NSN#		
OD Green	20-0029	6545-01-537-3138		
Digital Camo	20-0073			
Red	20-0077			



REATMENT SUPPORT | Apparel & Accessories

OPERATORS CAPNAR

- 100% cotton
- Adjustable strap
- 3.5 in. by 2.5 in. cut loop on the front designed for attaching patches (sold separately)
- 1 in. square of cut loop on the top for attaching an IR square

SPECIFICATIONS

One size fits most Weight: 3.2 oz

OPERATORS CAP

ITEM#

ZZ-0092

MEDIC CROSS PATCHES PVC

- Soft flexible PVC rubber patch
- Backed with hook material for easy attachment and removal
- Ideal for use on medical kits such as the M-FAK and T.O.R.K.
- 1 in. x 1 in. available at NARescue.com

MEDIC CROSS PATCHES

COLORS	ITEM#
Coyote & Black	ZZ-0797
Black & Gray	ZZ-0798
OD Green & Black	ZZ-0800
Black & Orange	ZZ-0799



SPECIFICATIONS

L 2 in. x W 2 in. x D 0.125 in. Weiaht: 0.6 oz

EMBROIDERED LOGO CAPNAR • All American Cap

- Made in the USA
- Comfortable, attractive Fit-4-You™ design fits most wearers
- American flag on the front, NAR logo on back

SPECIFICATIONS

Fit-4-You™ S/M Weight: 3.2 oz

EMBROIDERED LOGO CAP

ITEM#

ZZ-0176

IR FLAGSNAR

- Clearly visible under IR light
- Hook-and-loop attachment
- Sold in Packs of 2

SPECIFICATIONS

L 2 in. x W 2 in. x D 0.125 in. Weight: 0.6 oz

IR FLAGS

COLOR	STARS	ITEM#
Full Color	L	ZZ-0093
Coyote	L	ZZ-0094
Coyote	R	ZZ-0235
OD Green	L	ZZ-0107



MED ID PATCH LUMINOUS

- Soft flexible custom rubber PVC
- Ultimate glow-in-the-dark identification technology
- 1 in. block letters
- 2 medical patches per package

SPECIFICATIONS

Each: L 3 in. x W 2 in. x D 0.125 in. Weight: 0.4 oz

LUMINOUS MED ID PATCH

ITEM#

ZZ-0121



SQUARES NAR

- Infrared patch printed on honeycomb reflective material
- Black hook fasteners on the back for easy attachment
- Clearly visible under IR light illumination
- Six squares per pack

SPECIFICATIONS

Each: 1 in. x 1 in. Weight: 0.2 oz.

IR SQUARES

ITEM#

ZZ-0108



NAR **Training**



NORTH AMERICAN RESCUE

Education and Training

North American Rescue Education and Training (NAR Training)

has roots that run deep and wide in the world of pre-hospital medicine and technical rescue. Our team was built from the ground up on the vision and experience of seasoned military medics seeking to fill the training gaps commonly encountered in the operational environment.

Since then, we have continued to broaden our vision to include other invaluable first responder perspectives such as pre-hospital healthcare providers in civilian, law enforcement, emergency medical services and community public safety initiatives.

CORE Competencies

- Continuing Medical **Education (CMEs)**
- Prolonged Casualty Care
- Austere Medicine
- Technical Rescue
- Emergency Medical Services
- Critical Care Flight Medicine


ONLINE



designed to develop both competence and confidence. Using **Professional Role** Players, Moulage Artists, Subject Matter Experts, High Fidelity Manikins, Re-animated Cadavers and State of the Art Facilities, our team delivers unparalleled medical instruction to today's heroes.







ACCREDITATIONS:





INSTRUCTOR











➤ Phone: 888.689.6277 ➤ Email: NARTraining@NARescue.com

COURSES

PUBLIC ACCESS

NAR BCON

6-hour course oriented towards civilians; students are exposed to realistic trauma scenarios and taught how to treat major hemorrhaging injuries, receiving critical knowledge on how respond to preventable bleeding deaths prior to arrival of EMS. Stop the Bleed certification issued upon completion.

First Responder - Work Place Safety

8-hour course designed for civilians. Students will learn CPR, bleeding control, and basic first aid to the most common workplace injuries. AHA BLS and NAR First Responder certification upon successful test completion.

Civilian Medical

3-day medical course designed for civilians. Students will learn basic emergency medicine, med bag packing, and med supply procurement. AHA BLS and NAEMT TECC certification upon successful test completion.

First Responder - Trauma Care

8-hour course for law enforcement and first responders that will provide students the ability and confidence to apply trauma care in the field. Upon completion of testing, students receive NAEMT TECC certification.

NATIONAL REGISTRY

EMT Hvbrid

Self-paced (over period of one year), online learning designed to prepare entry-level student for written and practical application testing for National EMT testing and certification.

AEMT Hybrid

Self-paced (over period of one year), online learning designed to prepare EMT student for written and practical application testing for National AEMT testing and certification.



42-hour inhouse course for IBSC critical care course certification or recertification (FP-C or CCP-C) and national registry recertification. Students will learn or refresh all skills and knowledge needed to certify/re-certify, and will earn 42 hours of IBSC and National Registry Continuing Education.

Vent Hero

Supporting literature that provides an introduction into the theory of successful ventilation management geared towards the critical care provider and clinician responsible for managing complex patient ventilation support.

Essentials In Critical Care 4.0

Online refresher course which exceeds the minimum requirements (above 16 hrs) for approved IBSC certifications.

LAW ENFORCEMENT

> TECC-LE

8-hour course specifically designed for law enforcement. Covers hemorrhage control, airway control, splinting, extraction techniques, and stop the bleed scenarios. Upon passing test, students receive NAEMT TECC-LE certification.

FEDMED

5-day course designed specifically for federal law enforcement officers. Students learn how to diagnose and treat various illnesses and injuries utilizing role players and simulations. Upon successful completion, AHA BLS and NAEMT TECC certifications are achieved; for EMTs, this course also serves as a full EMT re-certification.



NORTH AMERICAN RESCUE®

Education and Training

➤ Phone: 888.689.6277 ➤ Email: NARTraining@NARescue.com

COURSES



MILITARY

TCCC-AC

This 40-hour course is designed for non-medical military personnel that provides a 4 year NAEMT TCCC certification upon completion. Students will learn emergency battlefield medicine focusing on preventable causes of death.

> TCCC-MP

This 16-hour course is designed for non-medical military personnel that provides a 4-year NAEMT TCCC certification upon completion. Includes standard training with addition of scenarios utilizing role players and moulage.

► EGT-AC

This 5-day course is designed for small military units that provides a 4-year NAEMT TCCC-AC certification. Provides students with TCCC medicine and drills, field maneuvers, and trauma lanes with role players and CASEVAC trainers.

➤ EGT-MP

This 5-day course is designed for military medical provider units that awards a 4-year NAEMT TCCC-AC certification upon completion. Provides students with TCCC medicine and drills, field maneuvers, and trauma lanes with role players and CASEVAC trainers.

► SOFMED

This 8-day course is specifically designed for Special Operations Medics that provides the following certifications: BLS, ACLS, PALS, TCCC-MP, and a full NCCP 60-hour CEU/ NRP refresher. Students focus on medicine while enhancing tactics via use of simulations, role players, and trainers.



► SOTM

This 8-day course is specifically designed for Special Operations Medics that provides the following certifications: BLS, ACLS, PALS, TCCC-MP, and a full NCCP 60-hour CEU/NRP refresher. Students focus on medicine while enhancing tactics via use of perfused cadavers, simulations, role players, and trainers.

CP MED

Toxic MED 101

2-day course for anyone who may be involved in a CBRN incident. Students will receive critical knowledge and skills on how to respond to incidents involving chemical/ biological/radiological/ nuclear contaminants. **NAEMT TECC certification** awarded upon completion.



Toxic MED 201

5-day course designed for medical providers/first responders. Provides background, knowledge and methodologies associated with responding to or supporting CBRN operations with casualties. Students will have the opportunity to receive the AHLS certification if desired. All students will earn NAEMT TECC certification upon successful completion.

ADVANCED


Austere Rescue Course

5-day course aimed at military medical providers and paramedics. Students will receive periods of instruction covering austere medicine, remote rescue, and patient transport in rugged terrain with limited resources. Course will involve treatment and movement in remote environments coupled with challenging problems involving emergency medicine, transport, and elements of PCC.



Austere Surgical Course

5-day course exclusively designed for military surgical teams. Students will utilize perfused cadavers to hone surgical skills in a variety of environments; battlefield conditions and equipment will be incorporated into training to increase proficiency and team skills. SOPs will be tested in preparation for full mission profile event.



5-day course that gives military medical providers the knowledge and skills necessary to treat patients in an austere environment for longer periods of time with limited resources when patient evacuation is not tenable.

Students will be instructed on how to keep patients stable with whatever resources are available until transport arrives.

Explore training opportunities online at NARTraining.com

Tactical Operation Medical Manikin SERIES

TOMManikin[™]

BASIC | BLAST | GSW | WATER | CONFINED SPACE | BURN | CBRNE | FEMALE

Rugged, state-of-the-art training simulators, NARS TOMManikin Series provides the comprehensive solution for enhanced training in Tactical Combat Casualty Care (TCCC) from Point of Injury to transfer of higher care!

- All-in-one, hyper-realistic trauma simulation solution enabling instructors to build multiple realistic traumatic casualty scenarios with reactive patients Ideal for use in classroom settings, trauma lanes and field training scenarios
- Multiple, interchangeable extremities present a variety of wounding patterns to provide full mission profiles in combat scenarios and disaster/mass casualty re-creations
- Solid metal bones and joints constructed of welded steel for unparalleled durability and quick exchange leg joints for flexibility in the most austere environments
- Proprietary silicone skins available in light, medium or dark skin tones that will not break down during prolonged UV exposure or extreme temperatures
- 180 lb manikin provides extremities capable of withstanding up to 6,000 lb of pressure -Ideal for building collapse/entrapment training scenarios beyond conventional exercises
- Ruggedized remote controls provide simultaneous long-range communications with multiple NARS devices all from a single device!



- MILSPEC Wireless Remote Instructor Tablet & Software
- Stand alone point-to-point communications with .5 mile line of site range

Flexible Wound Configurations:

Basic | Blast | GSW | Water | Confined Space | Burn | CBRNE Custom configurations including Female Anatomy on request

Clinical Interventions available on ALL TOMManikins: ◀

MASSIVE HEMORRHAGE:

SPECIAL FEATURES

Interactive Animatronics:

• Variable Radial Pulse Rates & Pressure

High Fidelity Bleeding from Multiple Wounds

• Customizable, Multi-language Voice Tracks

Rugged Softech® Skin in a variety of skin tones

Lifelike Breathing Mechanisms

Functional Anatomical Airway

- Tourniquet Application (proximal, junctional & abdominal)
- Active Bleeding Control by Direct **Pressure** on Pressure Points

Wound Packing

ADVANCED AIRWAY:

- Nasopharyngeal/Oral/ Supraglottic Airway Insertion
- Endotracheal Intubation
- Surgical Airway and Tissue **Dissection** with Realistic Bleeding

NORTH AMERICAN RESCUE®

Simulation







RESPIRATIONS:


- Chest Seal Application
- Bilateral Needle Decompression
- Chest Tube Capable

CIRCULATION:

- **IO** (sternal & humeral)
- Antecubital IV Access or uninjured upper limbs

RFV072021





TOMManikin Interventions

- Tourniquet Application (proximal, junctional, and abdominal tourniquets)
- Active bleeding controllable by direct pressure on pressure points
- Wound packing with realistic bleeding control
- Nasopharyngeal and oral airways, supraglottic airway insertion and endotracheal intubation
- Surgical airway with realistic bleeding and tissue dissection 35 operative sites included
- Remote controlled breathing rates: Slow, Normal, and Fast
- **Bilateral Needle Decompression** at either 2nd/3rd intercostal (midclavicular), or 5th/6th intercostal (mid-axillary) locations
- Realistic chest seal application & adhesion
- Chest Tube capable (bilateral)
- Radial pulses in uninjured upper extremities & bilateral carotid pulses
- Variable pulse rate & pressure (Normal, bounding, thready)
- IO (sternal & humeral) and antecubital IV access on uninjured upper limbs



All Manikins include a MILSPEC-rated instructor tablet

BASIC TOMManikin™

BASIC Wounding Patterns:

- 1. Partial Arm Amputation below the elbow with Massive bleeding and dramatic soft tissue damage, exposed hard bone and discernible tissue layers Responds to tourniquet application or brachial artery direct pressure application
 - 2. Full Amputation Leg with massive bleeding, extensive soft tissue damage, exposed hard bone, visible tissue layers. Responds to direct junctional pressure and tourniquet application

2



BASIC TOMMANIKIN*

ITEM#

93-0035

*Other configurations are available upon request

BLAST TOMManikin™

BLAST Wounding Patterns:

- 1. Open Jaw Fracture Head Skin with massive maxillofacial bone and soft tissue injuries, actively bleeds. Supports OP/NP/ET/Surgical airways
 - 2. Packable Abdominal Evisceration with intricate handcrafted attention to anatomical detail with visible tissue layers, and 10 in. of eviscerated bowel with peripheral ecchymosis surrounding large open wound. Massive bleeding with "clampable" vessel
- 3. Shrapnel Wound Arm with active bleeding, multiple soft tissue wounds, responds to tourniquet application
- 4. Full Amputation Leg with massive bleeding, extensive soft tissue damage, exposed hard bone, visible tissue layers. Responds to direct junctional pressure and tourniquet application

BLAST TOMMANIKIN ITEM#

93-0036

*Other configurations are available upon request



Supports chest seal application and pleural n decompression

2. *GSW to Groin* with GSW anterior entry/ posterior exit with soft tissue trauma to penis and avulsed scrotum with prolapsed testicle. Wounds present with peripheral wound site ecchymosis, exposed sub-q fat layer and provide attriction.

muscle striation. Entry wound is packable and responds to both peripheral and junctional tourniquet application and femoral pressure point application

3. *GSW Leg* with through & through GSW to leg above the knee. Entry wound presents with small entrance hole with circular powder burn. Exit wound has exposed fat and wound site ecchymosis. Variable bleeding rate. Responds to both peripheral and junctional tourniquet application and femoral pressure point application. Can be used as an impaled object using steel rebar, or wooden stake

4. Full Amputation Leg (Option)** with massive bleeding, extensive soft tissue damage, exposed hard bone, visible tissue layers. Responds to direct junctional pressure and tourniquet application

Other configurations are available upon request

**Leg Amputation Option 93-0369/93-0947

GSW TOMMANIKIN*

OPTION	ITEM#
GSW TOMM (Basic)	93-0040
GSW TOMM (w/Leg Amputation)	93-0369
GSW TAMikin*** (Basic)	93-0039
GSW TAMikin***(w/Leg Amputation)	93-0947

***Female version

Addressing injury locations based on US Army study suggesting females who suffer penetrating injuries to chest, abdomen, pelvis have a higher mortality rates than male colleagues

1

TAMikin gsw

Female version of TOMManikin GSW, with anatomically correct female features and voice

- Injuries located in sensitive areas such as under breast & pubic area as well as the leg
- Parallels GSW TOMManikin's capabilities & wounding patterns



*Other configurations are available upon request

WATER TOMManikin™

Wounding Patterns:

Partial Leg Amputation below the knee with massive soft tissue trauma and exposed hard bone, muscle and tissue layers.

Options:

Sucking Chest Wound with anterior GSW entry & posterior exit. Actively bleeds from both locations. Anterior wound mixes with air from onboard compressor to create bubbling bleeding effect. Supports chest seal application and pleural needle decompression

-or-

Packable Abdominal Evisceration with intricate handcrafted attention to anatomical detail with visible tissue layers, and 10 in. of eviscerated bowel with peripheral ecchymosis surrounding large open wound. Massive bleeding with "clampable" vessel

Closed Fracture Arm with ecchymosis and palpable crepitus.

Designed to be submerged up to 33 feet in either fresh or salt water

WATER TOMMANIKIN*

OPTION	ITEM#
Basic (Sucking Chest Wound)	93-0041
ABD Evisceration Closed Arm Fracture	93-0108







*Other configurations are available upon request

Wounding Patterns:

Sucking Chest Wound with anterior GSW entry & posterior exit. Actively bleeds from both locations. Anterior wound mixes with air from onboard compressor to create bubbling bleeding effect. Supports chest seal application and pleural needle decompression

Partial Leg Amputation below the knee with massive soft tissue trauma and exposed hard bone, muscle and tissue layers.

Shrapnel Wound Arm with active bleeding, multiple soft tissue wounds, responds to tourniquet application

CONFINED SPACE TOMMANIKIN*

ITEM#

93-0110







CONFINED SPACE TOMManikin™



Wounding Patterns:

Burn Face* with 1st, 2nd, 3rd-degree burns to face that supports all airway adjuncts. Head skin is interchangeable with other TOMManikins

BURN TOMManikin™

Burn Arms with 1st, 2nd, 3rd-degree circumferential burns to both arms

BURN TOMMANIKIN*

ITEM#

93-0037

*Other configurations are available upon request

**Head skin is interchangeable with other TOMManikins

CBRNE TOMManikin™

Wounding Patterns:

Burn Face* with 1st, 2nd, 3rd-degree burns to face that supports all airway adjuncts. Head skin is interchangeable with other TOMManikins

Sucking Chest Wound with anterior GSW entry & posterior exit. Actively bleeds from both locations. Anterior wound mixes with air from onboard compressor to create bubbling bleeding effect. Supports chest seal application and pleural needle decompression

Abdominal Evisceration (Wearable Moulage) with intricate handcrafted attention to anatomical detail with visible tissue layers. and 10 in. of eviscerated bowel with peripheral ecchymosis surrounding large open wound. Massive bleeding with 'clampable' vessel.

GSW to Groin with GSW anterior entry/posterior exit with soft tissue trauma to penis and avulsed scrotum with prolapsed testicle. Wounds present with peripheral wound site ecchymosis, exposed sub-q fat layer and muscle striation. Entry wound is packable and responds to both peripheral and junctional tourniquet application and femoral pressure point application



CBRNE TOMMANIKIN'

ITEM#

93-0038

*Other configurations are available upon request




K9 Features:

- Redesigned all-metal frame
- Expanded joint range of motion
- CPR chest
- 5 interchangeable limbs with various injuries: Complete Amputation, Open Fracture, Soft Tissue Injuries, Road Rash/Burns, & Gunshot Wound
- Improved airway
- Barks, Bleeds & Breathes
- 1L onboard blood capacity
- Wireless operation via NARS tablet









Features include:

- Pelican 1650 case with custom-cut foam to protect all system components
- 2 x Cameras
- 2 x 72 in. Telescopic Tripods with camera mounts
- 3 x Terabyte stand-alone recorder
- 1 x NARS mil-spec rated Instructor Tablet for wireless communications
- Proprietary software with instructor bookmark/replay feature, simulator status capture, trainee biometric capture, trainee evaluation software (exportable reports)
- 2 x 25 ft. CAT-6 retractable cord reels
- Trainee-worn stress monitor
- Features streaming participant heart and respiratory rate



MOBILE FEEDBACK SYSTEM

ITEM#

93-0945



Wearable Trainee

Stress Monitor

NEW H-60 ▲

Cockpit with Realistic Pilot/Co-Pilot Seats and Graphics Panel

CASEVAC HELOS Simulators HELOS

Features:

- Available in H60, CV22 & MH6 configurations
- Reusable, all-weather systems
- Includes loudspeakers for inflight noise, smoke generator to replicate inflight emergencies, wind machine to simulate rotor wash
- DOT-rated for over-the-road towing
- Platform-specific options include HD cameras (with playback system), FRIES bar, litter stanchions, blackout window drapes, ICS Kit (with headsets and cords)



RAINING | TCCC All Combatant Training Modules



TCCC SKILLS SYSTEM

All Combatant

- All the necessary tactical medical and casualty extraction/ evacuation equipment needed when teaching the TCCC - All Combatant Skills required tasks for completing the Tactical Combat Casualty Care (TCCC) course
- A turnkey solution for training centers and educational institutions designed to equip up to twenty-five (25) students
- Complete system consists of three self-contained Hard Storage Case modules: Skill Stations | IFAKS | Evacuation

SPECIFICATIONS

Hard Storage Case: External – \bar{L} 32.5 in. x W 21 in. x D 13 in. Internal - L 29 in. x W 18 in. x D 10 in. Weights: TCCC Skills Set: 55 lb TCCC IFAK Skills Set: 70 lb TCCC Evacuation Skills Set: 50 lb

TCCC SKILLS SYSTEM		
	UNIT	ITEM#
	3 Modules	80-0970

TCCC SKILLS SET All Combatant

• Train up to twenty-five (25) students at individual skill stations in the treatment of: extremity hemorrhage, wound packing, chest seal application, needle chest decompression, basic airway adjuncts, splinting, hypothermia prevention & management, and eye injuries

• Self-contained, wheeled, rigid case with internal dividers to keep the equipment organized and easy to inventory

SPECIFICATIONS

Hard Storage Case: External - L 32.5 in. x W 21 in. x D 13 in. Internal - L 29 in. x W 18 in. x D 10 in. Weight: 55 lb

TCCC SKILLS SET

MODULE	ITEM#
TCCC Skills Set	80-0968

































All Combatant TCCC IFAK SKILLS SET

SET MODULE:



Kit contents below ▼

- **TCCC IFAK SKILL** Twenty-five (25) pre-packed Training IFAKs for use during trauma lanes or training scenarios
 - IFAKs come packaged in a self contained wheeled hard, case that allows for easy transport and storage

SPECIFICATIONS

External Case -L 32.5 in. x W 21 in. x D 13 in. Internal Case -L 29 in. x W 18 in. x D 10 in. Weight: 70 lb

TCCC IFAK **SKILLS SET MODULE**

MODULE	ITEM#
TCCC IFAK Skills Set	80-0969



TCCC IFAK SKILLS KIT All Combatant

SPECIFICATIONS

H 7 in. x W 4.5 in. x D 2.75 in. Weight: 2 lb 8 oz

TCCC IFAK SKILLS KIT

OPTION	ITEM#
TCCC IFAK Skills Kit	85-1958
TCCC IFAK ReSupply*	80-1028









Tourniquet,Traine



1 NARETD"

6 in. Flat



Gauze

S-Rolled



Combat

Gauze® Traine

* TCCC IFAK ReSupply Kit contains all components of the IFAK SKILLS Kit, excluding the Nylon Bag

TCCC EVACUATION SKILLS SET

• Various high threat extraction devices and casualty evacuation systems to train students in different options to move casualties during the all three TCCC Phases of Care

• Packaged in a self-contained, wheeled, rigid case that allows for easy transport and storage

TCCC EVACUATION SKILL SET MODULE:







SPECIFICATIONS External Case -L 32.5 in. x W 21 in. x D 13 in. Internal Case -

L 29 in. x W 18 in. x D 10 in. Weight: 50 lb















TRAINING | TCCC Naval Training Kits

NAVAL FIRST AID BOX RESPONSE KIT All Service Member Trainer

- Designed to meet the TCCC requirements of the U.S. Naval Ships, Military Sealift Command ships and afloat commercial vessels
- Provides the All Service Members and first responders with rapid access and organization of critical lifesaving equipment for rendering on scene emergency care at sea in a training environment
- Lightweight, sand/water resistant, slimline pouch constructed of Berry Amendment 500D Cordura solution dyed nylon, Ghille Tex hardware, self repairing zippers & silent 550 cord pulls
- Pouch features securable anti-pilferage seals and a clear window to readily view contents
- Easily fits into the standard Naval First Aid Box

NAVAL FIRST AID BOX RESPONSE KIT TRAINER

ITEM#	NSN#
85-3204	6910-01-695-9363

SPECIFICATIONS

L 15 in. x H 6 in. x D 4 in. Weight: 2 lb

















1 28F

Naso Airway



1 CPRMask,

Res-Cue Key



1 Survival

Blanket





















NAVAL GUN RESPONSE AID KIT Trainer

- Designed to meet the TCCC requirements of the U.S. Navy for the emergency care equipment onboard U.S. Naval Ships, Military Sealift Command Ships and afloat commercial vessels
- Ergo dynamic, low profile, casualty response system featuring 2 identical removable equipment pouches that attach via hook and loop to a lightweight, durable outer bag
- Meets the requirements for the U.S. Navy's Naval Gun Response Aid Kit - Trainer
- Lightweight, sand/water resistant, slimline pouches constructed of Berry Amendment 500D Cordura solution dyed nylon, Ghille Tex hardware, self repairing zippers & silent 550 cord pulls
- Both outer bag & 2 inner pouches are securable using anti-pilferage seals
- Slim-line, lightweight kit stores easily on a wall, can be carried over the shoulder or bandolier-style

NAVAL GUN RESPONSE AID KIT TRAINER

ITEM#	NSN#
85-3205	6910-01-695-9287

SPECIFICATIONS

Outer Bag: H 16 in. x W 11 in. x D 5.5 in. Inner Pouch: H 15 in. x W 6 in. x D 4 in. Weight: 1 lb 15 oz

















































RAINING Individual Components

C·A·T®TRAINER

Combat Application Tourniquet

Fully functioning C-A-T® tourniquet colored blue to distinguish it as a training device as opposed to a personal/operational life-saving device



Tactical & Rescue versions of the C·A·T®

CAT TRAINER

ITEM#	NSN #
30-0033	6910-01-560-2972

SPECIFICATIONS

H 6.5 in. x W 2.4 in. x D 1.5 in. Weight: 2.7 oz

HYFIN® TRAINER Twin Pack

Unique adhesive gel formulation does not damage mannequins or leave behind excessive residue

SPECIFICATIONS

Packaged: L 7.5 in. x W 4.5 in. x D 0.25 in. Chest Seal: H 6 in. x W 6 in. x D 0.1 in.

Weight: 2.5 oz

See page 12 for other HyFin® Chest Seal products

HYFIN TRAINER

ITEM#	NSN #
10-0032	6910-01-908-8679

HUF MICHEST SEAL TRAINER



COMBAT GAUZE®

Moulage Trainer

Completely inert (no hemostatic properties) - training purposes only in a re-sealable package for multiple uses

SPECIFICATIONS

Packaged: H 6 in. x W 4.5 in. x D 1.25 in. Gauze: L 4 yd x W 3 in. Weight: 0.75 oz

COMBAT GAUZE TRAINER

ITEM#

30-0063

Responder

TRAINING ETD

Non-sterile, elastic pressure bandage for an economic, low cost training solution




ITEM#

30-0145

SPECIFICATIONS

Bandage: L 60 in. x W 4 in. Pad: H 6 in. x W 4 in. Packaged: L 4.5 in. x W 2.5 in. x D 1.75 in. Weight: 2.0 oz



SAM IO TRAINING KIT Single Bone

SAM IO SINGLE BONE

ITEM#

30-0256

• All of the equipment needed to familiarize trainees with the SAM IO Access hardware & usage

 Featuring single training humerus bone with tissue block & rugged nylon carrying case

SPECIFICATIONS

H 7 in. x W 11.25 in. x D 5.5 in. | Weight: 2.32 lb

SAM IO TRAINING KIT - SINGLE BONE:



















SAM IO FIELD STORAGE CASE

ITEM#

30-0254

SPECIFICATIONS H 7 in. x W 11.25 in. x D 5.5 in.

SAMIO

Weight: 1 lb





G3 Molded Hard Side Carrying Case **#30-0077** sold separately

G3 POWER DRIVER TRAINER

by EZ-IO®

- Battery-powered revolutionary technology
- Provides access within 10 seconds
- Multiple insertion sites with three needle sets
- Trigger guard prevents accidental activation of the device when stored in kits
- Standard Luer-Lock catheter for fast, easy removal
- · Lightweight, compact design
- Red coloring to immediately indicate it as a training only device

G3 POWER DRIVER TRAINER

ITEM#

30-0076

SPECIFICATIONS

H 4.5 in. x W 6.5 in. x D 2.5 in Weight: 15 oz

G3 EDUCATOR KIT by EZ-IO®

• Includes all of the Arrow® EZ-IO® System training equipment

• Allows clinicians to train in the correct use of the EZ-IO® System on adult and pediatric patients

KIT CONTENTS:

- G3 Training Power Driver
- EZ-IO® Training Needle Sets: 15mm Pediatric, 25mm Adult 45mm Large Adult
- 4 Wristbands
- 4 EZ-Connects
- G3 Vascular Access Pack
- 4 EZ-Stabilizers
- Training Bones: Infant, Tibia. Adult Proximal. Humerus Adult Tibia
- Softside Nylon Case
- Training CD

EZ-IO Products are restricted to Federal Purchases only.

Please contact **Customer Service** for details

G3 EDUCATOR KIT

ITEM#

30-0074



EZ-IO® G3 Power Driver

TRAINING NEEDLES

Designed for use with any EZ-IO® **Power Driver**



G3 TRAINING NEEDLES

SIZE	ITEM#
15 mm (Pediatric)	30-0082
25 mm (Adult)	30-0084
45 mm (Adult, large)	30-0086

T.A.L.O.N. MANUAL TRAINING NEEDLE

Train in manual I/O insertion at all seven optional sites, including the sternum

T.A.L.O.N. **MANUAL** TRAING NEEDLE

ITEM# 30-0133

MANUAL TRAINING NEEDLE by EZ-10

For training in six different FDA cleared sites for easily accessible insertion

EZ-IO TRAINING NEEDLE

30-0079

ITEM#



TRAINING | Part Task Trainers

WOUND MANAGEMENT **SIMULATOR**

- Rugged part-task trainer for learning essential lifesaving skills of tourniquet application, wound packing and applying pressure bandages
- Multiple wounding patterns to enhance training objectives

SPECIFICATIONS

L 15 in. x Diameter 6 in. Weight: 9.3 lb

- Simulates severe laceration. penetrating trauma (entrance and exit wounds), along with the ability to apply tourniquets
- Proprietary silicon that is rugged and durable and will not break down due to environmental exposure
- Can be used in the classroom or in field training scenarios

TRUECLOT Large Caliber GSW TRAINER:



TRUECLOT®



deep wound packing techniques



• Large caliber gunshot wound with hidden bleed points and a

simulated fracture/defect hard plastic bone insert to aid in teaching

• Sturdy silicone construction with no moving parts, reinforced tubing

mount, and a more realistic wound pattern and outer anatomy



Large Caliber GSW with Bone Insert

Hemorrhage Training Kit[™]









SPECIFICATIONS

ITEM#

93-0020

TRUECLOT LARGE CALIBER GSW **TRAINER**

Ideal for the repeated use of a training environment





5 TrueClot®

Flow Valves



















TRUECLOT® Hemorrhage Training Kit™

with 1 Gun Shot Wound & 1 Knife/Laceration

- Task Trainer for teaching realistic hemorrhage control in wounds including a large caliber gunshot & knife/laceration
- Training capabilities include general wound packing techniques & advanced hemostatic gauze training

TRUECLOT GSW & LACERATION TRAINER

ITEM#

93-0932



SPECIFICATIONS H 17 in. x W 20 in. x D 7.5 in. Weight: 14 lb



HEMORRHAGE CONTROL TRAINING KIT

with Bio Feedback

- Training tool designed to teach effective wound packing skills with Immediate results for both trainer and student
- Verifies correct or incorrect technique

SPECIFICATIONS

Case: 23 in. x 8 in. x 15 in. | Weight: 17 lb 3.2 oz

HEMORRHAGE CONTROL TRAINING KIT

Combat Gauze®

- Comes with 10 Combat Gauze Trainers & 5 Combat Gauze LE Demo Samples
- Includes a "training leg" that

COMBAT GAUZE

ITEM#

93-0017

93-0019

features an
entry/exit bullet
wound and a deep
laceration

SPECIFICATIONS

Case: H 15 in. x W 23 in. x D 8 in. Weight: 15 lb 8 oz



CUSTOMER SUPPORT

- COMMITMENT: We pride ourselves in providing an exceptional customer experience beginning with your initial inquiry up to order placement and beyond. Our customer service staff members are trained by certified instructors on all NAR products and are standing by to help with any questions you may have
- SUPPORT: Our in-house technical support professionals consist of an elite group of Nationally Registered Paramedics and Special Operations Medics from military, law enforcement, tactical law enforcement, fire-rescue and EMS backgrounds

ORDERING

ONLINE: www.NARescue.com

SPECIAL PRICING: automatically available ONLINE for all military customers with a .mil e-mail address! Register your account in just minutes by creating an online account using your .mil e-mail address

ESTABLISHED CUSTOMERS:

Want to order online with special pricing? Contact Customer Service at 888.689.6277 to have your online account set up

NEW CUSTOMERS: Simply create an online account and you will have access to customer specific pricing – anytime, anywhere you are logged in to NARescue.com

QUOTE REQUESTS: Online accounts have no waiting to receive your custom pricing to place an order or, we will gladly accept your submittal for a pricing quote request

• CONTRACT:

Tax ID: 27-1024029 Cage Code: 06ST7 DUNS: 83-24-26782 SAM Registered: Yes

ECAT: dmmonline.dscp.dla.mil

Prime Vendor: dmmonline.dscp.dla.mil/primevendor/primehome

• EMAIL: info@narescue.com

• **PHONE**: 888.689.6277

ADDRESS:

35 Tedwall Ct | Greer, SC 29650

• FAX: 864.675.9880

SHIPPING

- In-stock merchandise ships within two business days
- Need it ASAP? Ask Customer Service about overnight or second day delivery

