

PRODUCTS WITH A MISSION®



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

*Over twenty years of helping
to reduce preventable causes
of death on the battlefield
and in our home towns...*



Items marked with this icon have been specifically branded with NAR's signature Red Tip Technology® packaging. Use these items with confidence, knowing that they have been engineered for less-than-optimal, often hostile conditions that exist when our equipment is needed the most. *Accept no imitations!*

table of contents

volume 6 • issue 1

4 New Product Preview
Our newest product offerings from throughout the catalog including new kits, public access products and ballistic protection accessories

6 Armor and Ballistic Protection
Designed to help safeguard rescue personnel in "warm zone" situations

11 M.A.R.C.H. Products
Massive Hemorrhage, **A**irway Management, **R**espiratory Support, **C**irculation/Vascular Access and **H**ypothermia Prevention

22 Kits
Kits designed to meet the needs/challenges of our customers and to address the leading causes of potentially preventable death

22 Public Access Bleeding Control
Individual kits, packs and stations designed to equip bystanders and first responders with the essential medical gear needed to help save lives

48 Training
Products designed to help train users on TCCC/TECC protocols and decreasing preventable death at the point of wounding

54 TAC-EVAC
The most effective and highest quality casualty evacuation platforms and accessories for use in the tactical medical environment

58 Treatment Support
From personal protection products to splinting and immobilization products to help from the point of injury throughout the entire continuum of care

66 Fire Suppression
Patented CAFS (Compressed Air Foam System) technology that brings increased capabilities and performance in a cost effective solution

70 Rescue Equipment
Includes a wide array of products ranging from lightweight and versatile Pegasus® Ropes to tactical lights and carabiners

81 Index
Product offerings indexed by name, item number and NSN number for quick reference and information

OUR NEWEST product

SOFT LITTER CARRIER



See page 9 for more information

GEAR BAGS



See page 9 for more information

ETD™ 8x10

Emergency Trauma Dressing



See page 14 for more information

AMBU® RDIC

For use in chemical or biological conditions



See page 16 for more information

PUBLIC ACCESS

5-Pack Kits and Stations



See page 23 for more information

PUBLIC ACCESS

Semi-Recessed Metal Stations



See page 25 for more information

preview

RIGID C·A·T® HOLDERS



See page 12 for
more information

EMERGENCY SURVIVAL WRAP



See page 21 for
more information

PATROL VEHICLE TRAUMA KIT



See page 46 for
more information

ANKLE KITS



See page 27 for
more information

Public Access **BLEEDING CONTROL SKILLS TRAINING KITS**



See page 48 for
more information

BALLISTIC PROTECTION

AGENT RTF KIT

Lightweight and compact system

- An innovative, lightweight, compact AGENT Armor Carrier System incorporating Level IIIA soft armor protection and medical pouch all contained in a small, discreet soft-sided briefcase
- Medical pouch includes essential trauma equipment to provide immediate point-of-wounding care to injured casualties in accordance with TECC Guidelines
- Agent RTF kit includes two Rescue Task Force ID Placards and one Luminous Med ID Patch







ORDERING INFO

COLOR	ITEM #
Black	80-0410
Red	80-0411

SPECIFICATIONS

Stored: H 12 in. x W 15 in. x D 5 in.
Maximum Waist Measurement: 60 in.
Total Weight: 7 lb 4 oz

AGENT RTF MEDICAL KIT COMPONENTS:

 QTY 2 C-A-T® Tourniquet	 QTY 2 ETD 6 inch	 QTY 2 Gauze S-Rolled	 QTY 2 Naso Airway & Lube	 QTY 1 HyFin® Vent, Twin Pk	 QTY 1 Trauma Shears, Lg	 QTY 1 Scissor Leash™	 PAIR 3 Nitrile Gloves, Lg	 QTY 4 Combat Casualty Card	 QTY 2 Marker Permanent
---	--	--	--	--	---	---	---	--	--

Agent RTF kit includes two Rescue Task Force ID Placards and one Luminous Med ID Patch

AGENT ARMOR CARRIER SYSTEM

Lightweight and compact system

- Available as an individual carrier and Level IIIA soft armor system
- Lightweight, compact armor carrier system is easily deployed yet conveniently stores in vehicle compartments or passenger spaces
- Fully adjustable hook and loop internal soft armor and hard plate stabilizers
- Two integrated armor pouches that accommodate 14 in. x 11 in. Level IIIA soft armor and/or additional DKX Rifle Plate (10 in. x 12 in., Level III, sold separately)

SPECIFICATIONS

Stored: H 12 in. x W 15 in. x D 4.5 in.
Vest Front: H 12 in. x W 15 in. x 0.5 in.
Maximum Waist Measurement: 60 in.
Weight: Carrier only: 2 lb 6.4 oz Carrier with Soft Armor: 5 lb 3.2 oz

ORDERING INFORMATION

COLOR	STYLE	ITEM #
Black	Carrier & Soft Armor	70-0070
Red	Carrier & Soft Armor	70-0071
Black	Carrier Only	70-0072
Red	Carrier Only	70-0073



BALLISTIC PROTECTION

RESCUE TASK FORCE VEST KIT

Available with or without side armor

- Level IIIA Soft Armor Ballistic Protection
- Front and back hard armor plate pockets able to receive multiple size plates
- Optional cummerbund-style Side Armor Protection with additional plate pouch
- Mil-Spec industrial reinforced webbing
- Includes Medical Kit that addresses the leading causes of preventable death
- Adjustable hook and loop side straps – expandable up to 62 in.
- Rescue Handle – 1,200 lb tensile strength reinforced strap
- One of the largest soft armor coverage areas in the industry

SPECIFICATIONS

Vest with Medical Kit (stored): H 19 in. x W 13.5 in. x D 5 in.

Vest, Medical Kit & Side Armor (stored): H 19 in. x W 21 in. x D 6 in.


ID Placards: H 2 in. x W 11 in.

Maximum Waist 62 in.

Vest & Kit Weight: 9 lb

Vest, Kit & Side Armor Weight: 13 lb

RESCUE TASK FORCE VEST MEDICAL KIT CONTENTS:



Kit also includes two Rescue Task Force ID Placards & Luminous Med ID Patch

ORDERING INFORMATION

COLOR	STYLE	ITEM #
Black	Vest, Side Armor & Medical Kit	80-0365
Blue	Vest, Side Armor & Medical Kit	80-0448
Coyote	Vest, Side Armor & Medical Kit	80-0366
Red	Vest, Side Armor & Medical Kit	80-0385
Yellow	Vest, Side Armor & Medical Kit	80-0386
Black	Vest and Medical Kit	80-0363
Blue	Vest and Medical Kit	80-0447
Coyote	Vest and Medical Kit	80-0364
Red	Vest and Medical Kit	80-0383
Yellow	Vest and Medical Kit	80-0384

VEST,
SIDE ARMOR
& MEDICAL KIT

VEST &
MEDICAL KIT

Black, Blue and Yellow kits come standard with Black C-A-T® Holder
Coyote kits come standard with Coyote C-A-T® Holder
Red kits come standard with Red C-A-T® Holder

ARMOR

TRAINING

PARC

KITS

TRAINING

TAC-EVAC

SUPPORT

RESCUE

LIFESTYLE

INDEX

BALLISTIC PROTECTION

ARMOR

M-A-R-C-H

PABC

KITS

TRAINING

TAC-EVAC

SUPPORT

RESCUE

LIFESTYLE

INDEX

TACTICAL RESPONDER VEST

MKII - Available with and without side armor

VEST
& SIDE
ARMOR

VEST
ONLY

- Level IIIA Soft Armor Ballistic Protection
- Front and back hard armor plate pockets able to receive multiple size plates
- Optional cummerbund-style Side Armor Protection with additional plate pouch

SPECIFICATIONS

Vest (stored): H 19 in. x W 13.5 in. x D 1 in.

Weight Vest Only: 6 lb 3 oz

Vest & Side Armor (stored): H 19 in. x W 13.5 in. x D 3 in.

Maximum Waist Measurement: 62 in.

Weight with Side Armor: 10 lb 9 oz

ORDERING INFORMATION

COLOR	STYLE	ITEM #
Black	Vest & Side Armor	70-0019
Blue	Vest & Side Armor	70-0097
Coyote	Vest & Side Armor	70-0020
Red	Vest & Side Armor	70-0045
Yellow	Vest & Side Armor	70-0046
Black	Vest Only	70-0017
Blue	Vest Only	70-0094
Coyote	Vest Only	70-0018
Red	Vest Only	70-0041
Yellow	Vest Only	70-0042

See page 9 for accessories

RESCUE RESPONDER SIDE ARMOR

- Level IIIA Soft Armor ballistic protection
- Cummerbund-style side armor set compatible with the Tactical Responder Soft Armor Vest/Rescue Task Force Vest Kit
- Mil-Spec industrial reinforced webbing

SPECIFICATIONS

Each Panel: H 8 in. x W 28 in. x D 0.5 in.

Maximum Waist: 62 in. Weight: 4 lb 7 oz

ORDERING INFORMATION

COLOR	ITEM #
Black	70-0021
Blue	70-0095
Coyote	70-0022
Red	70-0043
Yellow	70-0044

BALLISTIC PROTECTION

SOFT LITTER HOLDER

Lightweight and low profile



- Designed to carry two (2) NAR QuikLitters or one (1) NAR QuikLitter XL (each sold separately)
- Litters can be quickly removed from either the bottom or top of carry device for self or buddy deployment
- Featuring a MOLLE-style attachment system that connects to the rear of protective vest leaving front of vest open for other mission critical equipment

SPECIFICATIONS

Folded: H 14 in. x W 14 in. x D 0.25 in.
Weight: 9 oz

See page 56 for
QuikLitter & QuikLitter XL

ORDERING INFO

COLOR	ITEM #
Black	80-0556
Blue	80-0557
Red	80-0558



DKX LEVEL III RIFLE PLATE

- Level III Hard Plate Level III Compliant
- Multi-Hit Capable and one of the lightest hard armor plates on the market
- Rugged, drop resistant plate that does not require x-ray examination

ITEM

70-0024

SPECIFICATIONS

H 12 in x W 10 in x D 11 in
Weight: 2 lb 14.4 oz

RESCUE TASK FORCE PLACARDS



- Integrates seamlessly with the Rescue Task Force Vest Kit and the Tactical Responder Soft Armor Vest

ORDERING INFO

COLOR	ITEM #
Black	ZZ-0509
Coyote	ZZ-0510
Red	ZZ-0537

SPECIFICATIONS

H 2 in. x W 11 in.
Weight: 3.2 oz

GEAR BAGS

- Single large compartment for maximum capacity and simplicity
- Large enough to hold a Tactical Responder Vest & Ballistic Helmet

SPECIFICATIONS

L 30 in. x W 14 in. x D 14 in.
Weight: 1 lb 4 oz



ORDERING INFO

COLOR	ITEM #
Black	80-0547
Blue	80-0477

BALLISTIC PROTECTION

BATLSKIN® VIPER HELMETS

Viper A1, Viper A3 and Viper P2

- Full cut, low profile series of ACH combat helmets designed to stop NIJ IIIA level threats (9mm and 44 mag) (NIJ 0106.01)
- Unsurpassed impact protection and ballistic rated hardware
- Tested for ballistic performance from -40°F to 158° F
- Compatible with NVGs, Chem/Bio Masks and communication devices
- Easily adjusts for a comfortable fit featuring movable comfort pads for impact protection, a 4-point attachment system with supportive neck pad

ORDERING INFORMATION

Size	Style	ITEM #
S	Viper A1	70-0025
M	Viper A1	70-0026
L	Viper A1	70-0027
XL	Viper A1	70-0028
S	Viper A3	70-0048
M	Viper A3	70-0049
L	Viper A3	70-0050
XL	Viper A3	70-0051
S	Viper A3-MSS	70-0081
M	Viper A3-MSS	70-0082
L	Viper A3-MSS	70-0083
XL	Viper A3-MSS	70-0084
S	Viper P2	70-0033
M	Viper P2	70-0034
L	Viper P2	70-0035
XL	Viper P2	70-0036

Batlskin® Viper A1

- Ultimate performance combat helmet built in the ACH shape to exact military specifications

Batlskin® Viper A3

- Standard A3 easily adjusts for fit with movable comfort pads for impact protection & a 4-point attachment suspension system
- A3 with MSS features the Modular Suspension System with 4-point integrated harness & unique buckle system

Batlskin® Viper P2

- Premium, lightweight ACH built to exact military specs
- Infinitely adjustable nape pad to ensure helmet stability



Coyote Helmets available see
NARescue.com for more information

OTW GOGGLES

Outside the Wire - with Case

- 100% UVA/B/C protection with maximum ballistic protection
- Tapered, neutral tint gray lenses with anti-fog and scratch resistant coating, reduced glare and correct distortion

ORDERING INFO

ITEM #
90-0066

SPECIFICATIONS

Goggles: W 7.5 in. x H 2.75 in.
Case: L 8.75 in. x W 3.5 in. x H 3 in.
Weight: 10.1 oz



IPRO® by NAR

Tactical Goggle System

- High performance eye protection with Hexo-vents for maximum ventilation, fine particle filtration and fog-free viewing

SPECIFICATIONS

Packaged: L 2.6 in. x W 7 in. x D 3 in.
Weight: 8 oz

ORDERING INFO

ITEM #
ZZ-0109



MASSIVE HEMORRHAGE

C·A·T® COMBAT APPLICATION TOURNIQUET

- **Single route buckle** decreases application time to help decrease blood loss, allows for slack removal for proper TQ application
- **Applied same way** to both upper and lower extremities
- True **one-handed tourniquet**
- **Scientifically proven** 100% effective in occluding blood flow in both upper and lower extremities by the USAISR and recommended by CoTCCC

SPECIFICATIONS

Packaged: H 6.5 in. x W 2.4 in. x D 1.5 in.
Device: L 37.5 in.
Weight: 2.7 oz

ORDERING INFORMATION

COLOR	STYLE	ITEM #	NSN #
Black	Tactical	30-0001	6515-01-521-7976
Orange	Rescue	30-0023	6515-01-612-7129
Blue	Trainer	30-0033	6910-01-560-2972



online video available

< scan to view



See page 49 for more about the Blue Trainer version of the C·A·T®

JUNCTIONAL EMERGENCY TREATMENT TOOL **JETT®**



- Immediately treat high extremity wounds in the thigh/groin area that is not amenable to TQ application.
- Unilateral or bilateral treatment solution
- Simultaneous occlusion of blood flow to both lower limbs that will not impede respiration

SPECIFICATIONS

Packaged: H 5.5 in. x W 7 in. x D 3.75 in.
Overall Device: L 56.75 in.
Buttocks Pad: L 19 in. x W 6.5 in.
Compression Pad: L 4.85 in. x W 2.5 in. x D 11 in.
T-Handle: H 4.62 in. x W 4.04 in. x D 0.63 in.
Weight: 1 lb 9.6 oz

ORDERING INFORMATION

ITEM #	NSN #
30-0088	6515-01-616-5841

JETT Pack (a nylon pouch designed specifically for carrying the JETT) is also available (Item# 80-0355 & 80-0356)
See NARescue.com for more information

online video available

scan to view >



MASSIVE HEMORRHAGE

RIGID C·A·T® TQ CASE GEN 7

- Allows you to carry the Gen 7 C·A·T® Tourniquet on your duty/patrol/ pistol belt so you have life saving equipment on your person at all times
- Blade-Tech TekLok Belt Attachment fits up to a 2.25 in. belt
- Can be mounted vertically or horizontally
- Injection Molded Nylon Case stands up to rugged use

ORDERING INFORMATION

Color	ITEM #
Black	80-0479
Black, Basket Weave	80-0560
Black, Polished	80-0562
Blue	80-0480
Coyote	80-0481
Gray	80-0484
Multi-Cam	80-0561
Orange	80-0483
Ranger Green	80-0482



SPECIFICATIONS

Packaged: H 4.75 in. x W 2.25 in. x D 2 in.
Weight: 3.8 oz

TQ 911 C·A·T® HOLDER

- Solid one-piece injection molded design
- Allows fast and easy deployment when needed

SPECIFICATIONS

Case: H 4 in. x W 2.25 in. x D 2.25 in.
Weight: 1.6 oz

ORDERING INFO

ITEM #
30-0160



C3 - CONCEALED C·A·T® CARRIER

- Soft armor, under-the clothing concealed carrying option
- Easily mounts to body armor in multiple number of locations

SPECIFICATIONS

Carrier: L 6.5 in. x W 2.25 in. x D 0.5 in.
Weight: 0.13 oz

ORDERING INFO

ITEM #
30-0137



RIGID C·A·T® TQ CASE & COVER

- Designed to stand up to rugged use
- Cover is made of Hypalon One-Piece Construction

SPECIFICATIONS

L 6.75 in. x W 2.25 in. x D 2.5 in.
Weight: 4.5 oz

ORDERING INFO

Color	ITEM #
Black	80-0485
Coyote	80-0486
Ranger Green	80-0487



C·A·T® HOLDER

- Tailored specifically for your C·A·T® to protect it from the elements

SPECIFICATIONS

L 6.25 in. x W 1.75 in. x D 1.5 in.
Weight: 1.26 oz

ORDERING INFO

NATO NSN*

Color	ITEM #	NSN #
Black	30-0057	----
Coyote	30-0044	----
Multi-Cam	30-0059	----
OD Green	30-0053	6530-01-605-9432*
Red	30-0128	----



See page 27 for our new Ankle Trauma Holster and Kit

MASSIVE HEMORRHAGE

COMBAT GAUZE® & COMBAT GAUZE® LE

- Sterile, non-woven Z-fold gauze impregnated with kaolin for temporary external control of bleeding
- Effective, simple and safe to use on wounds
- Does not generate heat
- Features an embedded x-ray detectable strip

SPECIFICATIONS

Standard and LE: L 4 yd x W 3 in.

Pkg: H 5.75 in. x W 4.5 in. x D 0.25 in. Weight: 0.8 oz

XL: L 4 yd x W 4 in. (2 ply)

Pkg: H 6 in. x W 5 in. x D 0.75 in. Weight: 1.75 oz

ORDERING INFORMATION

Product	ITEM #	NSN #
Standard	30-0039	6510-01-562-3325
XL	30-0070	6510-01-600-4705
LE	30-0109	----
Trainer	30-0063	----

See page 50 for more information on the Trainer Version



ARMOR

M.A.R.C.H

PABC

KITS

TRAINING

TAC-EVAC

SUPPORT

RESCUE

LIFESTYLE

INDEX

QUIKCLOT® BLEEDING CONTROL DRESSING



- Advanced hemostatic agent combined with the simplicity of gauze available without medical device authorization (MDA)
- Utilizes the same inorganic component as QuikClot Combat Gauze® used by the US Military to control severe bleeding

SPECIFICATIONS

H 5.75 in. x W 4.5 in. x D 0.65 in.

Gauze: L 4 yd x W 3 in. Weight: 0.8 oz

ORDERING INFO

ITEM

30-0161

CHITOGAUZE®

- Proven hemostatic properties of chitosan & the familiar convenience of non-woven, polyester/rayon blend medical gauze

SPECIFICATIONS

Gauze: L 4 yd x W 3 in.

Packaged: H 6 in. x W 5 in. x D 0.65 in.

Weight: 0.8 oz

ORDERING INFORMATION

ITEM #	NSN #
30-0049	6510-01-591-7740



CELOX™

- Stops life-threatening bleeding resulting from shock as found in severe traumatic injuries

SPECIFICATIONS

5ft Gauze: L 5 ft x W 3 in.

Packaged: H 5.25 in. x W 4.25 in. x D 0.5 in. Weight: 1.1 oz

10ft Gauze: L 10 ft x W 3 in.

Packaged: H 6.75 in. x W 4.25 in. x D 0.75 in. Weight: 2.1 oz

ORDERING INFORMATION

SIZE	ITEM #	NSN #
5 Foot	30-0064	6510-01-623-9910
10 Foot	30-0065	6510-01-619-8636



MASSIVE HEMORRHAGE

COMPRESSED GAUZE

- Low-cube, ultra-compact sterile gauze
- Premium rolled 100% cotton gauze for maximum absorption

SPECIFICATIONS

Gauze:
W 4.5 in. x L 4.1 yd
Packaged:
L 3 in. x W 2.5 in. x D 1 in.
Weight: 12 oz



ORDERING INFORMATION

ITEM #	NSN #
30-0052	6510-01-503-2117

S-ROLLED GAUZE

- Unique dispenser, allows controlled linear feed application
- Premium 100% cotton

SPECIFICATIONS

Gauze:
W 4.5 in. x L 4.1 yd
Packaged:
L 3 in. x W 2.5 in. x D 1 in.
Weight: 12 oz



ORDERING INFORMATION

ITEM #	NSN #
30-0003	6510-01-605-2234

ETD™ EMERGENCY TRAUMA DRESSING

- Premium cotton construction for max quality & absorption
- Durable C-clasp securing device



ORDERING INFORMATION

SIZE	STYLE	ITEM #	NSN #
4 inch	Standard	30-0031	6510-01-558-4108
4 inch	Flat	30-0093	----
6 inch	Standard	30-0032	6510-01-558-4114
6 inch	Flat	30-0094	----
Abd/Stump	----	30-0012	6510-01-541-8121
8x10 inch	----	30-0165	6510-01-654-6801

See NARescue.com for ETD™ specifications

RESPONDER ETD™ EMERGENCY TRAUMA DRESSING

- Sterile pressure bandages that deliver extremely effective pressure and excellent absorption
- Suitable for applying immediate direct pressure



ORDERING INFORMATION

SIZE	STYLE	ITEM #
4 inch	Flat	30-0171
6 inch	Rolled	30-0168

See NARescue.com for Responder ETD™ specifications

CRAVAT (TRIANGULAR BANDAGE)

Bandage: 37 in. x 37 in. x 52 in.
Pkg: 3.5 in. x 2 in. x 0.63 in. Wgt: 1.6 oz

ORDERING INFORMATION

ITEM #	NSN #
30-0089	6510-00-201-1755



TACTICAL COMPRESSION WRAP

Wrap: L 38 in. x W 5 in. Wgt: 6oz
Pkg: 5.4 in. x 2.7 in. x 2 in.

ORDERING INFORMATION

ITEM #	NSN #
30-0016	6510-01-614-0938



AIRWAY MANAGEMENT

TACTICAL CRICKIT® Tactical Surgical Cricothyroidotomy Kit

- Contains all features required by CoTCCC for a surgical airway kit in the tactical setting
- Ideal length & size tracheal tube with securing device
- Configured in the sequence of use, ideal for high-stress, low-light, austere environments
- Lightweight, ruggedized and minimal in cube size
- Included rolled pouch provides a clean working surface

INCLUDED IN THE TACTICAL CRICKIT:



SPECIFICATIONS

Packaged: H 6.5 in. x W 2.5 in. x D 1.25 in.
Unrolled: H 6.5 in. x W 7.5 in. x D 1.25 in.
Weight: 3.7 oz

ORDERING INFORMATION

ITEM #	NSN #
10-0017	6515-01-540-7568

NASOPHARYNGEAL AIRWAY & LUBRICANT

- Packaged with surgical lubricant attached

L 6 in. x W 16 in. x D 0.5 in.
Weight: 0.5 oz

ORDERING INFORMATION

SIZE	ITEM #	NSN #
28F	ZZ-0034	6515-01-529-1187

22F-36F available at NARescue.com

NAR TACTICAL SUCTION DEVICE

Suction Capacity:
1000 mL
Vacuum Force:
100 mmHg
Weight: 7.3 oz
L 9.75 in. x W 3.75 in. x D 3 in.

ORDERING INFORMATION

ITEM #	NSN #
10-0018	6515-01-540-7206

CYCLONE® BVM - IMPROVED POCKET



- Unique collapsible design packaged 25% smaller than conventional BVMs
- Requires minimum pressure for compression of the bag to ventilate

Diameter: 5.3 in. x H 2.9 in.
Weight: 17.4 oz

ORDERING INFORMATION

ITEM #	NSN #
10-0019	6515-01-568-0193

AMBU® SPUR II BVM



L 11.63 in. x W 5 in. x H 5 in.
Weight: 11.2 oz

ORDERING INFORMATION

ITEM #	NSN #
10-0043	6515-01-519-3612

NAR TRACHEAL HOOK



L 4.2 in. x W 0.5 in. x D 0.8 in.
Weight: 0.5 oz

ORDERING INFORMATION

ITEM #	NSN #
10-0020	6515-01-562-9734

AIRWAY MANAGEMENT

KING LTS-D™ Supraglottic Airway Kit



- Unique dual channel design, anatomical shaped distal tip/cuff & open second lumen
- Tapered ramp design provides additional ventilation outlet

SPECIFICATIONS

Folded: H 8 in. x W 11.43 in. x D 1.92 in.
Unfolded: H 8 in. x W 13.36 in. x D 1.22 in.
Weight: 5.7 oz

INCLUDED WITH KING LTS-D™



QTY
1 Sterile
Lubricant



QTY
1 Gastric
Tube, 18F



QTY
1 Syringe
Luer-Lock 60cc

ORDERING INFORMATION

CLR	SIZE	ITEM #	NSN #
YEL	3	10-0025	----
RED	4	10-0026	6515-01-591-1578
PUR	5	10-0027	6515-01-591-0397

AMBU® RDIC Mark III Resuscitator



- Genuine Ambu Mark III reusable resuscitator equipped with impermeable elastomeric components on the compression bag and patient valve (butyl rubber)
- Specialized inlet filtration system, including C2 canister, for filtering contaminated air

ORDERING INFORMATION

ITEM #	NSN #
10-0047	6515-01-338-6602

SPECIFICATIONS

L 14 in. x W 7.5 in. x H 7 in.
Weight: 5.4 lb

Special Operations LARYNGOSCOPE SET



- Perform direct visual laryngoscopy even in the most severe environments
- Lightweight, compact kit to equip medical care providers in the field with tools to secure an airway
- Contains 2 single-use fiber optic blades, textured handle & removable light source all rolled up in a rugged nylon pouch

SPECIFICATIONS

Open: L 8 in. x W 5.5 in. x D 2 in.
Rolled: L 8 in. x W 2 in. x D 2 in.
Weight: 10 oz

INCLUDED IN THE LARYNGOSCOPE SET:



QTY
1 Pouch,
Nylon



QTY
1 Light,
Pocket Exam



QTY
1 #2 Blade,
Miller



QTY
1 #3 Blade,
Mac



QTY
1 Handle
Laryngoscope



QTY
1 Spare
Bulb



QTY
2 Size AA
Batteries

ORDERING INFORMATION

ITEM #	NSN #
10-0001	6515-01-450-9790

NAR FIBER OPTIC LARYNGOSCOPE BLADES

- Disposable, safe and convenient

ORDERING INFORMATION

AVAILABLE SIZES	ITEM #	NSN# available on all sizes offered
Miller #1, #2, #3 & #4	See NARescue.com	See NARescue.com for details
Mac #2, #3 & #4	See NARescue.com	See NARescue.com for details




RESPIRATORY SUPPORT

HYFIN® CHEST SEALS



HyFin® Vent Twin Pack

- Two individually packaged, vented/valve-type chest seals in one unique package designed for treatment of both entry/exit wounds or multiple penetrating injuries
- Patented 3 channel pressure relief vents
- Ruggedized dual perforated pouch designed to allow the provider to open only one dressing at a time, as needed
- Advanced adhesive technology designed to provide a superior seal



ORDERING INFORMATION

ITEM #	NSN #
10-0037	6510-01-642-6212

SPECIFICATIONS

Device: 6 in. x 6 in. (excludes tab) Wgt: 2.5 oz
 Packaged: H 7.5 in. x W 4.5 in. x D 0.25 in.
 Unfolded: H 7.5 in. x W 9 in. x D 0.13 in.



online video
available

< scan to view

HyFin® Vent Individual

ORDERING INFORMATION

ITEM #	NSN #
10-0029	6510-01-624-0840

SPECIFICATIONS

Device: 6 in. x 6 in. (excludes tab) Wgt: 1.31 oz
 Packaged: H 4.63 in. x W 7.5 in. x D 0.1 in.



HyFin® Occlusive Individual

ORDERING INFORMATION

ITEM #	NSN #
10-0015	6510-01-532-8019

SPECIFICATIONS

Device: 6 in. x 6 in. (excludes tab) Wgt: 1.31 oz
 Packaged: H 4.63 in. x W 7.5 in. x D 0.1 in.



HyFin® Vent Compact, Twin Pack

- Two Chest Seals packaged approximately 25% smaller than our standard HyFin®
- 3-channel pressure relief vents
- Two Chest Seals designed for the treatment of entry/exit or multiple penetrating wounds
- Large red tab designed for single step, peel and apply application and for burping the wound if necessary

ORDERING INFO

ITEM #
10-0042

SPECIFICATIONS

Device: 4.75 in. x 4.75 in. (excludes tab)
 Wgt: 1.55 oz
 Packaged: H 3.88 in x W 6.6 in x D 0.13 in.



Also available at
www.NARescue.com

HyFin® Twin Pack Trainer
 Item# 10-0032
 See page 49 for more information

RESPIRATORY SUPPORT

ARS[®] NEEDLE DECOMPRESSION KIT

- Penetrate the pleural space in 99% of casualties presenting the signs and symptoms of a tension pneumothorax
- Quality, reliable devices with the appropriate length needle/catheter
- Capless flash chamber provides immediate confirmation of correct placement
- Patented design packaged in an easy-to-open, crush-proof protective tube featuring easy-to-identify textured twist top with clip

SPECIFICATIONS

Pkg: H 5.5 in. x W 0.7 in. x D 0.8 in.
Weight: 0.5 oz

ORDERING INFORMATION

GAUGE	ITEM #	NSN #
14 ga	ZZ-0056	6515-01-541-0635
10 ga	ZZ-0298	----

U.S. Design Patents: D584,409 S; D595,847S; D584,410S; U.S. Patent 7,874,426B2; © North American Rescue,® LLC., 127820, 127822, 127823, 001002372-0001, 001013940-0001, 001013940-0002



online video available

< scan to view

NAR NEEDLELESS SALINE LOCK KIT

- IV safety catheter & needle-free valve port reduce or eliminate accidental needle-stick exposures

ORDERING INFORMATION

ITEM #	NSN #
30-0046	6515-01-607-1943

SPECIFICATIONS

Pkg: H 5.5 in. x W 3.25 in. x D 1 in.
Weight: 1.3 oz



NEEDLELESS SALINE LOCK KIT CONTENTS:



Meets all requirements for US Army Task: 081-831-1011

NAR SALINE LOCK KIT

- Essential components for establishing a saline lock for safer, quicker and more effective IV cannulation in the most demanding environments

ORDERING INFORMATION

ITEM #	NSN #
30-0013	6515-01-537-4094

SPECIFICATIONS

Pkg: H 5.5 in. x W 3.25 in. x D 1 in.
Weight: 1.3 oz



SALINE LOCK KIT CONTENTS:



BOA® IV CONSTRICTING BAND

- Direct peripheral IV access utilizing a revolutionary rolling action for transfusing venous blood to the distal veins while still allowing arterial flow
- Quick Click connector & two-finger Quick Release improves "time on task"
- XL version stretches approximately 50% longer than the standard BOA®

SPECIFICATIONS

Standard: L 7.5 in. to 18 in. XL: L 10 in. to 27 in.
Packaged: L 4 in. x W 3 in. x D 0.9 in. Weight: 1.5 oz

US Patents: 5,304,202 & 5,607,448



online video
available

< scan to view



ORDERING INFORMATION

SIZE	ITEM #	NSN #
Standard	30-0009	6515-01-537-2611
XL	30-0071	----

CIRCULATION ACCESS

SHARPS SHUTTLE

SPECIFICATIONS

L 6.38 in. x Diameter 1.11 in.
Weight: 15 oz

ORDERING INFORMATION

UNIT	ITEM #	NSN #
Single	30-0014	6515-01-518-8551
6 pack	30-00146	6515-01-619-5671



RAPTOR® IV Securing Device

SPECIFICATIONS

Extended: L 19 in. x W 1 in.
Weight: 19 oz

ORDERING INFORMATION

UNIT	ITEM #	NSN #
6 pack	30-0006	6515-01-523-3317



G3 POWER DRIVER by EZ-IO®

- Battery-powered solution for providing rapid vascular access to humerus & tibia in less than 10 seconds (Not for sternal use)
- Trigger guard prevents accidental activation
- Molded Hard Side Carrying Case (sold separately)



NOTE: EZ-IO Trainers also available - see page 52

ORDERING INFORMATION

UNIT	ITEM #
Driver	30-0075
Case	30-0077

SPECIFICATIONS

H 4.5 in. x W 6.5 in. x D 2.5 in.
Weight: 15 oz

MANUAL EXTREMITY NEEDLE SET

Manual I/O insertion at 3 extremity sites by EZ-IO®

MANUAL NEEDLE SET CONTENTS:



ORDERING INFO

SIZE	ITEM #
25mm	30-0078



SPECIFICATIONS

Device: L 6 in. x W 3.25 in.
Weight 1.4 oz

POWER DRIVER NEEDLE SET

for use with the EZ-IO® G3 Power Driver

ORDERING INFORMATION

CLR	SIZE	ITEM #
Pink	15mm/PD	30-0081
Blue	25mm/AD	30-0083
Yellow	45mm/LD	30-0085

SPECS

Packaged:
H 4.75 in.
W 4.25 in.
D 0.7 in.
Weight: 15 oz

POWER DRIVER NEEDLE SET:



T.A.L.O.N.™ MANUAL NEEDLE SET

Manual I/O insertion at 7 sites – including the sternum

ORDERING INFORMATION

SIZE	ITEM #
38.5mm	30-0132

SPECIFICATIONS

Device: L 6 in. x W 3.25 in.
Weight 1.5 oz



HYPOTHERMIA PREVENTION

HPMK® HYPOTHERMIA PREVENTION & MANAGEMENT KIT

- Includes reinforced Heat Reflective Shell (HRS™) of 4-ply, composite fabric with a protected non-conductive reflective layer
- Built-in hood, integrated fluid absorption pad & top-to-bottom taper providing excellent thermal insulation
- Continuous 1.5 in. Velcro® closures provide rapid 360° access to the casualty with minimal exposure to the elements
- Self-heating, 4-cell liner (Ready-Heat™ Blanket) provides 10 hours of continuous dry heat up to 125°F (52°C)

SPECIFICATIONS

Shell Packaged: H 6.75 in. x W 10.5 in. x D 5.5 in.

Shell Open: L 78 in. x W 43 in. US Patent: 7,766,950 B2

Weight: 3 lb 8 oz

HPMK® CONTENTS:



ORDERING INFORMATION

ITEM #	NSN #
80-0027	6515-01-532-8056



HRS® HEAT REFLECTIVE SHELL



- As an individual unit, reinforced 4-ply, composite shell with built-in hood, weight tested to support 250+ lb.

ORDERING INFORMATION

ITEM #	NSN #
ZZ-0051	6532-01-627-7888


SPECIFICATIONS

Shell Packaged: L 10 in. x W 6.75 in. x D 2.5 in.

Open: L 78 in. x W 43 in. US Patent: 7,766,950 B2

Weight: 1 lb 13 oz

READY HEAT™ BLANKET



- As an individual unit, the self-heating, 4-cell liner provides 10 hours of continuous dry heat after reaching peak temperature approximately 15-20 minutes after opening

ORDERING INFORMATION

ITEM #
ZZ-0052

SPECIFICATIONS

Open: L 48 in. x W 36 in.

Weight: 1 lb 11 oz

EMERGENCY SURVIVAL WRAP

- Designed to reflect heat back to the body

ORDERING INFORMATION

ITEM #
ZZ-0834

SPECIFICATIONS

Unfolded: 60 in. x 90 in.

Weight: 3.7 oz





PUBLIC ACCESS KITS

INDIVIDUAL BLEEDING CONTROL KITS • NYLON BAGS or VACUUM-SEALED

Each Individual Bleeding Control Kit Contains These Items:



PLUS, choose from 4 different levels, each adding additional capabilities to your kit (see level choices below)

Plus, these Individual Components Based Upon Your Desired Level:

 QTY 1 Gauze Compressed BASIC LEVEL also includes an additional Compressed Gauze (2 total)	or	 QTY 1 Gauze Compressed QTY 1 HyFin® Compact, Twin INTERMEDIATE LEVEL also includes additional Compressed Gauze (2 total) & Hyfin® Vent Compact, Twin Pack	or	 QTY 1 QuikClot Dressing QTY 1 HyFin® Compact, Twin ADVANCED BCD LEVEL also includes QuikClot Dressing & Hyfin® Vent Compact, Twin Pack	or	 QTY 1 Combat Gauze® LE QTY 1 HyFin® Compact, Twin ADVANCED LEVEL also includes Combat Gauze LE & Hyfin® Vent Compact, Twin Pack
---	----	---	----	--	----	---

NYLON BAG PUBLIC ACCESS INDIVIDUAL KITS

- Essential equipment needed for persons to take action as immediate responders in helping to stop and/or control life threatening bleeding
- Easy to restock in cases of use
- Convenient to carry and easy to open in stressful situations
- Opens flat for ease of use and immediate identification of enclosed devices
- Easily recognizable red nylon bag with attached PVC Bleeding Control Patch



ORDERING INFORMATION

Kit Level	ITEM #	Weight
Basic	80-0452	15.0 oz
Intermediate	80-0453	16.5 oz
Advanced BCD	80-0523	17.4 oz
Advanced	80-0454	17.4 oz

SPECIFICATIONS

H 7 in. x W 5 in. x D 2.75 in
Weight: see table to the left

VACUUM-SEALED PUBLIC ACCESS INDIVIDUAL KITS



- Essential equipment needed for persons to take action as immediate responders in helping to stop and/or control life threatening bleeding
- Cost-effective alternative to nylon, compact and rugged
- Large identifying card vacuum-sealed inside the pouch for easy kit identification

ORDERING INFORMATION

Kit Level	ITEM #	Weight
Basic	80-0465	11.7 oz
Intermediate	80-0466	13.3 oz
Advanced BCD	80-0522	14.0 oz
Advanced	80-0467	14.0 oz

SPECIFICATIONS

H 7 in. x W 5 in. x D 2.75 in
Weight: see table to the right



online video
available

< scan to view

PUBLIC ACCESS PACKS

BLEEDING CONTROL PACKS



5-PACK VAC

BLEEDING CONTROL VACUUM-SEALED

- Five Individual Bleeding Control Vacuum-Sealed Pouch Kits contained inside a zippered nylon tote bag
- Also contains one NAR QuikLitter for moving casualties
- Tote Dimensions: H 7.5 in. x W 13 in. x D 8 in.



5-PACK VAC CONTENTS:

- 5x Individual Bleeding Control Kits (Vacuum-Sealed Pouches)
- 1x NAR QuikLitter
- 1x Nylon Zippered Tote Bag

ORDERING INFORMATION

Kit Level	ITEM #	Weight
Basic	80-0573	5 lb 11 oz
Intermediate	80-0574	6 lb 8 oz
Advanced BCD	80-0576	6 lb 8 oz
Advanced	80-0575	6 lb 8 oz

8-PACK VAC

BLEEDING CONTROL VACUUM-SEALED

- Eight Individual Bleeding Control Vacuum-Sealed Pouch Kits contained inside an easy to open nylon carrier
- Also contains one NAR QuikLitter for movement of casualties
- Folded Dim: H 16 in. x W 16 in. x D 8 in.



8-PACK VAC CONTENTS:

- 8x Individual Bleeding Control Kits (Vacuum-Sealed Pouches)
- 1x NAR QuikLitter
- 1x Nylon Folding Carrier

ORDERING INFORMATION

Kit Level	ITEM #	Weight
Basic	80-0468	6 lb 13 oz
Intermediate	80-0469	7 lb 10 oz
Advanced BCD	80-0524	8 lb 1 oz
Advanced	80-0470	8 lb 1 oz

8-PACK NYLON

BLEEDING CONTROL NYLON BAG

- Eight Individual Bleeding Control Nylon Bag Kits contained inside an easy to open nylon carrier
- Also contains one NAR QuikLitter for moving casualties
- Kits can be easily removed via hook and loop for individual deployment
- Folded Dimensions: H 16 in. x W 16 in. x D 8 in.



8-PACK NYLON CONTENTS:

- 8x Individual Bleeding Control Kits (Nylon Kit Bags)
- 1x NAR QuikLitter
- 1x Nylon Folding Carrier

ORDERING INFORMATION

Kit Level	ITEM #	Weight
Basic	80-0460	8 lb 8 oz
Intermediate	80-0461	9 lb 5 oz
Advanced BCD	80-0525	9 lb 11 oz
Advanced	80-0462	9 lb 11 oz



PUBLIC ACCESS STATIONS

BLEEDING CONTROL CLEAR POLYCARBONATE WALL STATIONS



SPECIFICATIONS

H 16 in. x W 16 in. x D 8 in.

Empty Case Weight: 10 lb

Complete Wall Station Weight: see individual tables

Empty Case with Placards: Item# ZZ-0822

CLEAR POLYCARBONATE • 8-PACK WALL STATION

- Includes Individual Bleeding Control Kits packaged in either red nylon bags or vacuum-sealed pouches stored within an easy-to-open nylon pack carrier and one (1) NAR QuikLitter for moving injured casualties

8-PACK WALL STATION CONTENTS:

- 8x Individual Bleeding Control Kits (Nylon Bags or Vacuum-Sealed)
- 1x NAR QuikLitter
- 1x Nylon Folding Carrier
- 1x Polycarbonate Wall Station Case with identifying placards

8-PACK STATION ORDERING INFORMATION

Kit Level	#Kits	Kit Packaging	ITEM #	Weight
Basic	8	Vac-Sealed	80-0471	16 lb 13 oz
Intermediate	8	Vac-Sealed	80-0472	17 lb 10 oz
Advanced BCD	8	Vac-Sealed	80-0526	18 lb 1 oz
Advanced	8	Vac-Sealed	80-0473	18 lb 1 oz
Basic	8	Nylon Bags	80-0474	18 lb 8 oz
Intermediate	8	Nylon Bags	80-0475	19 lb 5 oz
Advanced BCD	8	Nylon Bags	80-0527	18 lb 1 oz
Advanced	8	Nylon Bags	80-0476	19 lb 11 oz

CLEAR POLYCARBONATE • 5-PACK WALL STATION

5-PACK WALL STATION CONTENTS:

- 5x Individual Bleeding Control Kits (Vacuum-Sealed Pouches)
- 1x NAR QuikLitter
- 1x Nylon Tote Bag
- 1x Polycarbonate Wall Station Case with identifying placards

5-PACK STATION ORDERING INFORMATION

Kit Level	#Kits	Kit Packaging	ITEM #	Weight
Basic	5	Vac-Sealed	80-0581	15 lb 10 oz
Intermediate	5	Vac-Sealed	80-0582	16 lb 8 oz
Advanced BCD	5	Vac-Sealed	80-0584	16 lb 8 oz
Advanced	5	Vac-Sealed	80-0583	16 lb 8 oz



Polycarbonate Clear Wall Stations also available with Alarmed Cases visit www.NARescue.com for ordering information



online video
available
< scan to view

PUBLIC ACCESS STATIONS

BLEEDING CONTROL SEMI-RECESSED WALL STATIONS



ARMOR

M-A-R-C-H

PABC

KITS

TRAINING

TAC-EVAC

SUPPORT

RESCUE

LIFESTYLE

INDEX



SEMI-RECESSED METAL • 8-PACK WALL STATION

- Includes Individual Bleeding Control Kits packaged in either red nylon bags or vacuum-sealed pouches stored within an easy-to-open nylon carrier and one (1) NAR QuikLitter™ for moving injured casualties

8-PACK WALL STATION CONTENTS:

- 8x Individual Bleeding Control Kits (Nylon Bags or Vacuum-Sealed)
- 1x NAR QuikLitter
- 1x Nylon Folding Carrier
- 1x Semi-Recessed Metal Wall Case with identifying placards

8-PACK STATION ORDERING INFORMATION

Kit Level	#Kits	Kit Packaging	ITEM #	Weight
Basic	8	Vac-Sealed	80-0861	16 lb 13 oz
Intermediate	8	Vac-Sealed	80-0862	17 lb 10 oz
Advanced BCD	8	Vac-Sealed	80-0863	18 lb 1 oz
Advanced	8	Vac-Sealed	80-0864	18 lb 1 oz
Basic	8	Nylon Bags	80-0865	18 lb 8 oz
Intermediate	8	Nylon Bags	80-0866	19 lb 5 oz
Advanced BCD	8	Nylon Bags	80-0867	18 lb 1 oz
Advanced	8	Nylon Bags	80-0868	19 lb 11 oz

SPECIFICATIONS

Outside: H 17.5 in. x W 17.5 in. x D 1 in.
Inside: H 14.13 in. x W 14.13 in. x D 6.88 in.
Empty Case Weight: 10 lb
Complete Wall Station Weight: see individual tables
Empty Case with Placards: Item# ZZ-0833

SEMI-RECESSED METAL • 5-PACK WALL STATION

5-PACK WALL STATION CONTENTS:

- 5x Individual Bleeding Control Kits (Vacuum-Sealed Pouches)
- 1x NAR QuikLitter
- 1x Nylon Tote Bag
- 1x Semi-Recessed Metal Wall Case with identifying placards

5-PACK STATION ORDERING INFORMATION

Kit Level	#Kits	Kit Packaging	ITEM #	Weight
Basic	5	Vac-Sealed	80-0577	15 lb 10 oz
Intermediate	5	Vac-Sealed	80-0578	16 lb 8 oz
Advanced BCD	5	Vac-Sealed	80-0580	16 lb 8 oz
Advanced	5	Vac-Sealed	80-0579	16 lb 8 oz



Semi-Recessed Metal Wall Stations also available with Alarmed Cases visit www.NARescue.com for ordering information



M-FAK

- Contains the essential point-of-wounding equipment needed for treating penetrating, blast and other traumatic injuries
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access

SPECIFICATIONS

H 6 in. x W 3 in. x D 3 in.


Weight: 13 oz

ORDERING INFORMATION


COLOR	LEVEL	ITEM #
Black	Basic	80-0494
Coyote	Basic	80-0495
Gray	Basic	80-0872
Multi-Cam	Basic	80-0496
OD Green	Basic	80-0497
Black	with Combat Gauze	80-0498
Coyote	with Combat Gauze	80-0499
Gray	with Combat Gauze	80-0873
Multi-Cam	with Combat Gauze	80-0500
OD Green	with Combat Gauze	80-0501
Black	Advanced	80-0502
Coyote	Advanced	80-0503
Gray	Advanced	80-0874
Multi-Cam	Advanced	80-0504
OD Green	Advanced	80-0505




BASIC M-FAK CONTENTS:



BASIC WITH COMBAT GAUZE M-FAK CONTENTS:



ADVANCED M-FAK CONTENTS:



See our Lifestyle Section on page 78 for the compatible Medic Cross Patch

Bag only option also available at NARescue.com



ANKLE TRAUMA KIT

- Very lightweight, low profile, versatile carrier that allows for carry of essential trauma equipment
- Designed to accommodate a C-A-T® Tourniquet, Emergency Trauma Dressing and Gauze/Hemostatic
- Advanced breathable materials that provide user comfort and long-term wearability
- Protective, water resistant pouches for each carried component with strong hook and loop closures
- Quick, single-handed, ambidextrous access to components
- One size fits most with convenient hook and loop adjustability

SPECIFICATIONS

L 17.65 in. x H 6.5 in. x D 1.5 in.

Weight: 8.25 oz

ORDERING INFORMATION

Kit/Holster Option	ITEM #
Kit with Combat Gauze LE	80-0870
Kit with S-Rolled Gauze	80-0869
Ankle Holster Only	30-0177T

ANKLE TRAUMA KIT • S-ROLLED



ANKLE TRAUMA KIT • COMBAT GAUZE LE





EAGLE IFAK

- Can connect to the assault vest/body armor or to a duty/pistol belt
- Contains the essential equipment needed for self-aid/buddy-aid to address the leading causes of potentially preventable death in a tactical environment
- Quick Pull Tab easily opens the kit with just one hand
- Unique, detachable shock-cord inner panel that can be adjusted to your components
- Inner panel is removable for easy access to your kit

SPECIFICATIONS


H 8 in. x W 4 in. x D 3 in.
Weight: 1 lb 2 oz

ORDERING INFORMATION


COLOR	LEVEL	Hemostatic	ITEM #
Black	Basic	None	85-0177
Coyote	Basic	None	85-0178
Multi-Cam	Basic	None	85-0179
OD Green	Basic	None	85-0180
Black	Advanced	Combat Gauze	85-0189
Coyote	Advanced	Combat Gauze	85-0190
Multi-Cam	Advanced	Combat Gauze	85-0191
OD Green	Advanced	Combat Gauze	85-0192
4 options online	Advanced	Chito Gauze	See NARescue.com
4 options online	Advanced	Celox Gauze	See NARescue.com

Bag only option also available at NARescue.com

BASIC EAGLE IFAK CONTENTS:



ADVANCED EAGLE IFAK CONTENTS:



RIG SERIES EAGLE IFAK

- Tourniquet attached for quick deployment by banding it to inside of front opening flap
- Protective hood with quick pull latch release helps protect IFAK components and includes sewn in embroidered medical cross identifier
- Detachable (hook & loop backed) shock-cord inner panel adjusts to secure internal components and is easily removable for rapid access

SPECIFICATIONS

H 8 in. x W 4 in. x D 3 in.
Kit Packed Weight: 1.18 lb
Empty Bag Weight: 7 oz

ORDERING INFORMATION


COLOR	LEVEL	ITEM #
Black	Basic	80-0509
Coyote	Basic	80-0510
Multi-Cam	Basic	80-0511
OD Green	Basic	80-0512
Black	with Combat Gauze	80-0517
Coyote	with Combat Gauze	80-0518
Multi-Cam	with Combat Gauze	80-0519
OD Green	with Combat Gauze	80-0520
Black	Advanced	80-0513
Coyote	Advanced	80-0514
Multi-Cam	Advanced	80-0515
OD Green	Advanced	80-0516

Bag only option also available at NARescue.com




BASIC RIG SERIES EAGLE CONTENTS:



RIG SERIES EAGLE WITH COMBAT GAUZE CONTENTS:



ADVANCED RIG SERIES EAGLE CONTENTS:



NAR-4™ CHEST POUCH

- Streamlined chest pouch that provides quick access to essential equipment
- Designed to be secured to the front of the body armor using a MOLLE-style attachment system
- Ambidextrous trauma shear accessibility



NAR-4 CHEST POUCH CONTENTS:

 QTY 2 C-A-T® Tourniquet	 QTY 2 ETD™ 6 inch	 QTY 2 Gauze S-Rolled	 QTY 2 Choice of Hemostatic
 QTY 2 ARS® Needle, 14g	 QTY 2 Naso Airway & Lube	 QTY 1 CricKit Tactical	 QTY 1 HyFin® Vent, Twin Pk
 QTY 1 Trauma Shears, Lg	 QTY 1 Scissor Leash™	 PAIR 3 Nitrile Gloves, Lg	 QTY 2 PES Eye Shield
 QTY 2 Combat Casualty Card	 QTY 1 Marker Permanent	SPECIFICATIONS H 7.75 in. x W 9.75 in. x D 4.5 in Weight: 2 lb 10.3 oz	

ORDERING INFORMATION

COLOR	Hemostatic	ITEM #
Black	Combat Gauze	80-0173
Coyote	Combat Gauze	80-0175
Digital Camo	Combat Gauze	80-0176
Multi-Cam	Combat Gauze	80-0247
OD Green	Combat Gauze	80-0289
5 options online	Chito Gauze	See NARescue.com
5 options online	Celox Gauze	See NARescue.com

Bag only option also available at NARescue.com

BELT TRAUMA KIT (BTK)

- Rugged black pouch fits a belt up to 2.25 in. and can be worn rotated 90° to make room for other gear
- Mesh pockets allow quick and easy contents identification

SPECIFICATIONS

BASIC: H 4.5 in. x W 6 in. x D 2.5 in

Weight: 6.6 oz





COMPLETE: H 4.5 in. x W 6.75 in. x D 3.75 in.

Weight: 10.5 oz

ORDERING INFORMATION

Level	ITEM#
Basic	80-0426
Complete	80-0427

BASIC BELT TRAUMA KIT CONTENTS:

 QTY 1 ETD™ 4 in., Flat	 QTY 1 Combat Gauze LE	 QTY 1 Shield, CPR Face	 PAIR 1 Nitrile Gloves, Lg
--	---	---	---

COMPLETE BELT TRAUMA KIT CONTENTS:

 QTY 1 C-A-T® Tourniquet	 QTY 1 C-A-T® Holder	 QTY 1 ETD™ 4 in., Flat	 QTY 1 Combat Gauze LE	 QTY 1 Shield, CPR Face	 QTY 1 Nitrile Gloves, Lg
---	---	---	---	--	--





- Individual hemorrhage control kit
- Compact, vacuum-sealed & packaged to fit into a pocket, vest pouch or individual bag
- Rugged, durable pouch features NAR's signature Red Tip Technology® ensuring quality & performance

SPECIFICATIONS

H 4 in. x W 6.5 in. x D 2.75 in.
Weight: 7.4 oz

ORDERING INFORMATION

HEMOSTATIC/GAUZE	ITEM#
S-Rolled Gauze	80-0167
Combat Gauze, LE	80-0168
ChitoGauze	80-0223
Celox Gauze	80-0276

IPOK CONTENTS:



IPOK HOLDER

- Protects your medical kit from the elements
- Enclosure can be secured by Fastex snap buckle, hook-and-loop or both

ORDERING INFO

Color	ITEM#
Black	80-0205
Coyote	80-0207
Digital Camo	80-0208
OD Green	80-0254
Red	80-0358



SPECIFICATIONS

H 6.5 in. x W 6 in. x D 3.5 in.
Weight: 4.1 oz

VISOR MED POUCH

- Securely mounts to a vehicle visor
- Orange webbing handle to grab for easy deployment
- Pouch with easy-open buckle



ORDERING INFO

Color	ITEM#
Black	80-0269
Coyote	80-0270
Red	80-0360

SPECIFICATIONS

H 6.5 in. x W 5.5 in. x D 3 in.
Weight: 5.4 oz

CORE KIT

- Compact design with the essential life-saving equipment needed to address moderate and severe traumatic bleeding
- Contains the #1 CoTCCC Recommended Prehospital Tourniquet and Hemostatic Gauze

CORE CONTENTS:



ORDERING INFO

ITEM #
85-0834

SPECIFICATIONS

H 7 in. x W 3.5 in. x D 2.25 in.
Weight: 6.4 oz





S.T.O.R.M. OPERATOR IFAK

- Provides the necessary equipment to administer immediate self-aid/buddy-aid to casualties in hostile environments
- Attaches directly to the duty belt as a drop down leg rig or to the assault vest/body armor with MOLLE/PALS-style connectors

ORDERING INFO

COLOR	ITEM #
Black	80-0133
Coyote	80-0135
Digital Camo	80-0136
OD Green	80-0134

Bag only option also available at NARescue.com

SPECIFICATIONS

H 7.4 in. x W 6.9 in. x D 5 in.
Weight: 1 lb 12 oz



S.T.O.R.M. OPERATOR IFAK CONTENTS:



HIGH RISK WARRANT CASUALTY

- Designed with the capability of both casualty extraction & medical treatment in a high threat situation
- Provides the capability to address the leading causes of preventable death in the tactical environment.
- Equipped with a Dragon® Handle System to perform rapid assault rescue
- Multiple carriage capabilities (fanny pack, over the shoulder or attachment to packs with MOLLE-style webbing)

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.
Weight: 4 lb 12.8 oz

HIGH RISK WARRANTY CASUALTY KIT CONTENTS:



Bag only option also available at NARescue.com

ORDERING INFO

COLOR	ITEM #
Black	80-0121
OD Green	80-0122



Bag only option
also available at
NARescue.com

- Streamlined, compact IFAK designed to MOLLE to the assault vest or body armor
- Includes essential life-saving equipment to be used at the point of wounding
- For self-aid/buddy-aid to address injuries caused by penetrating or blast trauma

SPECIFICATIONS

H 7 in. x W 4.5 in. x D 2.75 in.
Weight: 1 lb 5 oz

ORDERING INFORMATION

COLOR	LEVEL	Hemostatic	ITEM #
Black	Basic	None	85-0008
Coyote	Basic	None	85-0009
Multi-Cam	Basic	None	85-0012
OD Green	Basic	None	85-0054
Black	Advanced	Combat Gauze	85-0048
Coyote	Advanced	Combat Gauze	85-0049
Digital Camo	Advanced	Combat Gauze	85-0050
Multi-Cam	Advanced	Combat Gauze	85-0052
OD Green	Advanced	Combat Gauze	85-0053
4 options online	Advanced	Chito Gauze**	See NARescue.com
5 options online	Advanced	Celox Gauze	See NARescue.com

BASIC T.O.R.K. CONTENTS:



ADVANCED T.O.R.K. CONTENTS:



NAR QUICK DETACH PANEL

- Quick release accessory compatible with any Mil-spec standard MOLLE/PALS system
- Multiple sizes available at www.NARescue.com
- Available in Black, Coyote, Digital Camo, Multi-Cam and OD Green

ORDERING INFO

COLOR	ITEM #
Coyote	80-0233

See NARescue.com for more options






RAPID RESPONSE KITS

LE EDITION

- Modular, personalized organization of mission essential gear and medical equipment
- Includes 2 C-A-T® tourniquets with holders, double mag pouch and nylon pouch to store additional equipment

RAPID RESPONSE KIT - LE EDITION CONTENTS:



Gun Magazines
are not included
in mag pouch



ORDERING INFORMATION

ITEM #

85-0715

SPECIFICATIONS

H 12 in. x W 12 in. x D 6 in.



Shoulder strap: 28 in. to 50 in.

Weight: 3 lb 12.5 oz

RTF EDITION

- Equipped to provide initial triage of a mass casualty incident and provide treatment for 4 or more patients
- Includes 4 C-A-T® tourniquets with holders and nylon dispenser with 4 rolls of color-coded triage tape

RAPID RESPONSE KIT - RTF EDITION CONTENTS:



Red Rapid
Response
RTF Kit comes
standard with
orange C-A-T®
Tourniquets,
red holders
and a red tape
dispenser

ORDERING INFO

Color	ITEM #
Black	85-0716
Red	85-0717

SPECIFICATIONS

H 12 in. x W 12 in. x D 6 in.

Shoulder strap: 28 in. to 50 in.

Weight: 6 lb 4.5 oz





MEDIC TRAUMA PACK

- 7 removable, color-coded pull bags and 2 removable side compartments provide modular organization
- Over 1,000 cubic inches of gear capacity and MOLLE-style attachment system compatible with removable side pouches

ORDERING INFORMATION

**NATO NSN


COLOR	ITEM #	NSN #
Black	80-0070	6545-01-572-7251
Coyote	80-0072	6545-01-588-1862
Digital Camo	80-0096	----
OD Green	80-0071	6545-01-594-5617**

SPECIFICATIONS

H 19 in. x W 19 in. x D 10 in
Weight: 18 lb 1 oz

Bag only option also available at NARescue.com

MEDIC TRAUMA PACK CONTENTS:



OPERATOR BLS IFAK

- Essential life-saving equipment that can be used for self-aid/ buddy-aid at the point of wounding
- Attaches directly to the duty belt as a drop down leg rig or to the assault vest or body armor with MOLLE/PALS-style connectors

OPERATOR BLS IFAK CONTENTS:



SPECIFICATIONS

H 8.5 in. x W 8.75 in. x D 4.5 in.
Weight: 1 lb 12 oz

ORDERING INFORMATION

COLOR	ITEM #
Black	80-0117
OD Green	80-0118

Bag only option also available at NARescue.com

888-689-6277 • www.NARescue.com

ARMOR

M-A-R-C-H

PABC

KITS

TRAINING

TAC-EVAC

SUPPORT

RESCUE

LIFESTYLE

INDEX



TACTICAL RAPID DEPLOYMENT

- Medium size go-bag with necessary equipment to treat critical casualties
- Kit and its contents are compact, durable and allow reliable and versatile interventions wherever the need arises

ORDERING INFORMATION

COLOR	ITEM #
Black	80-0201
Coyote	80-0203
Digital Camo	80-0204
OD Green	80-0202



TACTICAL RAPID DEPLOYMENT KIT CONTENTS:



SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.
Weight: 4 lb 12.8 oz

Bag only option also available at NARescue.com



EMS RAPID DEPLOYMENT KIT

- Highly visible, compact kit for pre-hospital healthcare providers
- Provides necessary equipment to treat several casualties

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.
Weight: 4 lb 6.5 oz

ORDERING INFORMATION

Color	ITEM #
Orange	80-0151

EMS RAPID DEPLOYMENT KIT CONTENTS:



Bag only option also available at NARescue.com

NAR-4™ AID KIT

- Low-profile, load carriage for Tactical Healthcare Professionals
- Modular system with 5 segregated removable equipment pouches to carry equipment consistent with TCCC Guidelines
- Designed to attach directly to the body armor with four rapid attachment straps or carried using shoulder straps

SPECIFICATIONS

H 19 in. x W 15 in. x D 7 in.

Weight: 12 lb 8.3 oz

NAR-4 AID KIT CONTENTS:

 QTY 2 C-A-T® Tourniquet	 QTY 2 ETD® 6 inch	 QTY 1 ETD Abdom.	 QTY 4 Gauze S-Rolled	 QTY 2 Choice of Hemostatic
 QTY 1 Tactical Comp. Wrap	 QTY 2 Naso Airway & Lube	 QTY 1 King LTS-D	 QTY 1 Cyclone BVM, Pocket	 QTY 1 Tactical Suction Dev.
 QTY 1 CricKit Tactical	 QTY 1 HyFin® Vent, Twin Pk	 QTY 2 ARS® Needle, 14g	 QTY 1 BOA® Standard	 QTY 4 Saline Lock Kit
 QTY 1 Sharps Shuttle	 QTY 1 Splint, Tact. Tract.	 QTY 2 SAM Splint	 QTY 1 Trauma Shears, Lg	 QTY 1 7 Hook Safety Cutter
 QTY 1 Case, Armadillo	 PAIR 8 Nitrile Gloves, Lg	 QTY 6 PES Eye Shield	 QTY 4 Combat Casualty Card	 QTY 1 Casualty Ref. Card
			 QTY 2 Tape, 2 in. Gecko-Grip®	
				 QTY 1 Marker Permanent

ORDERING INFORMATION

COLOR	Hemostatic	ITEM #
Black	Combat Gauze	80-0181
Coyote	Combat Gauze	80-0183
Digital Camo	Combat Gauze	80-0184
Multi-Cam	Combat Gauze	80-0244
OD Green	Combat Gauze	80-0308
5 options online	Chito Gauze	See NARescue.com
5 options online	Celox Gauze	See NARescue.com

Bag only option also available at NARescue.com



RANGE TRAUMA KIT

- Essential equipment for medical and non-medical personnel to address injuries common to firearm ranges and training facilities
- Treatment capabilities include penetrating injuries, chest wounds, burns, eye injuries, minor lacerations, sprains and fractures

ORDERING INFORMATION

*NATO NSN#

Hemostatic	Color	ITEM #	NSN #
Combat Gauze	Orange	80-0299	----
Chito Gauze	Orange	80-0213	6545-01-610-5411*
Celox Gauze	Orange	80-0298	----

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.
Weight: 5 lb 10 oz



RANGE TRAUMA KIT CONTENTS:



Bag only option
also available at
NARescue.com



HARD CASE RANGE TRAUMA AID KIT

- Water-tight case maintains the integrity of your medical products when you need them most
- Kit contents arranged in a modular system featuring rugged lid storage organizer and removable nylon pouches
- Labeled internal pouches for ease of use.

ORDERING INFO

Color	ITEM #
Orange	85-0889

SPECIFICATIONS

H 9.5 in. x W 13.5 in. x D 6 in.
Weight: 8 lb 12 oz



RANGE TRAUMA AID KIT HARD CASE CONTENTS:



K-9 HANDLER IFAK

K-9 HANDLER IFAK CONTENTS:



Bag only option also available at NARescue.com



- Designed for immediate lifesaving treatment to a wounded handler, operator or K-9
- Treatment capabilities range from major trauma typical to operating in tactical or patrol environments to K-9 emergencies

SPECIFICATIONS


H 7 in. x W 4 in. x D 2.75 in

Weight: 1 lb 3 oz

ORDERING INFORMATION

Color	Hemostatic	ITEM #
Black	----	80-0439
Coyote	----	80-0440
Multi-Cam	----	80-0441
OD Green	----	80-0442
Black	Combat Gauze	80-0443
Coyote	Combat Gauze	80-0444
Multi-Cam	Combat Gauze	80-0445
OD Green	Combat Gauze	80-0446

K-9 TACTICAL FIELD KIT



- Designed for the handler or medic to deliver lifesaving treatment to a wounded K-9 at the point of injury
- K-9 treatment capabilities range from major trauma typical to operating in tactical or patrol environments to more common K-9 emergencies

- Provides the essential equipment to render care to the injured K9 when veterinary care is not immediately available



Bag only option also available at NARescue.com

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.
Weight: 5 lb 2 oz

ORDERING INFORMATION

Color	Hemostatic	ITEM #
Black	Combat Gauze	80-0300
Coyote	Combat Gauze	80-0301
OD Green	Combat Gauze	80-0303
3 options online	Chito Gauze	See NARescue.com
3 options online	Celox Gauze	See NARescue.com

K-9 Litter also available - for more information visit www.NARescue.com



CRISIS INCIDENT RESPONSE KIT

- Designed to treat penetrating trauma, major bleeding, rescue breathing and basic first aid on multiple injured people
- Rugged carrying bag constructed of durable nylon with multiple carry options

ORDERING INFO

COLOR	ITEM #
Orange	85-1338
Black	85-0410

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.
Weight: 5 lb

CRISIS INCIDENT RESPONSE KIT CONTENTS:



Bag only option also available at NARescue.com




INDIVIDUAL AID KIT

- Designed for non-medical trained personnel to help control bleeding until medically trained responders arrive on scene
- Vacuum packed & easy-to-open



INDIVIDUAL AID KIT CONTENTS:



ORDERING INFO

ITEM

85-0404

SPECIFICATIONS

H 5.75 in. x W 3.5 in. x D 2 in.
Weight: 6.6 oz



MASS CRISIS INCIDENT KIT


- Designed to store & transport Individual Aid Kits to classrooms, cafeterias, auditoriums or other areas where a large mass of injured persons are gathered
- Contains 40 Individual Aid Kits (see above for Individual Aid Kits contents) packaged in a medium-sized red duffel style bag

ORDERING INFO

ITEM #
85-0411

SPECIFICATIONS

H 10 in. x W 20 in. x D 10 in.
Weight: 16.5 lb



SRO INDIVIDUAL RESPONSE KIT

- Designed to treat penetrating trauma and severe bleeding for immediate care to injuries
- Attaches to the duty belt as a drop down leg rig
- Includes an extraction device for casualty evacuation

ORDERING INFO

ITEM


85-0163

SPECIFICATIONS

H 8.5 in. x W 8.75 in. x D 4.5 in.

Weight: 2 lb 4.8 oz

SRO INDIVIDUAL RESPONSE KIT CONTENTS:



Bag only option also available at NARescue.com

SRO CRISIS RESPONSE KIT

- Equip School Resource Officers with the essential life-saving tools to begin treating injured victims until medical personnel arrive
- Empowers an SRO (School Resource Officer) to administer self-aid in the event he or she is injured or provide immediate care to injured casualties during a school incident
- Multiple carry options

ORDERING INFO

ITEM

85-0417

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.

Weight: 5 lb 16 oz

SRO INDIVIDUAL RESPONSE KIT CONTENTS:



Bag only option also available at NARescue.com

ARMOR

M.A.R.C.H.

PARC

KITS

TRAINING

TAC-EVAC

SUPPORT

REF

FE

INDEX



MEDIC TRAUMA SHEET

- Versatile, rugged, lightweight equipment access platform for managing casualties during TACEVAC operations
- Oversized pockets & large grommet attachment points facilitate rapid set-up and equipment retrieval
- Quick attach and release straps allow for secure attachment to the inside wall of an aircraft or vehicle

ORDERING INFORMATION

Color	ITEM #	NSN #
Black	80-0067	6515-01-540-7612
Coyote	80-0069	6515-01-540-7603
Digital Camo	80-0094	6515-01-540-7579
OD Green	80-0068	6515-01-540-7585

SPECIFICATIONS

Open: H 41.5 in. x W 21.5 in. x D 2.5 in.
Folded: H 20.25 in. x W 21.5 in. x D 5 in.
Weight: 4 lb 12 oz

Bag only; contents are sold separately



OCHO BAGS & ACCESSORIES

- Featuring 2 interchangeable mesh aid pouches*, hook and loop fields and hidden compartments – all within purposely designed exteriors that resemble standard off-the-shelf bags to minimize the threat of theft

OCHO SHOULDER BAG

SPECIFICATIONS

H 12 in. x W 14 in. x D 6 in.
Weight: 1 lb 7 oz

Color	ITEM #
Blue/Gray	80-0257




OCHO BACKPACK

SPECIFICATIONS

H 19 in. x W 14 in. x D 9 in.
Weight: 3 lb 14 oz

Color	ITEM #
Blue/Gray	80-0256
Red/Gray	80-0338
Green/Gray	80-0339
Gray/Black	80-0340



OCHO IFAK/ACCESSORY BAG

Color	ITEM #
Blue/Gray	80-0371
Gray/Black	80-0370

SPECIFICATIONS

H 7 in. x W 2.25 in. x D 3.25 in. Weight: 2.75 oz



**NOTE: OCHO IFAK/Accessory Bag does NOT include Mesh Aid Pouches*

OCHO BRIEFCASE

SPECIFICATIONS

H 12 in. x W 17 in. x D 4.5 in.
Weight: 2 lb 10 oz

Color	ITEM #
Brown	80-0258



OCHO QUICK DETACH GEAR ORGANIZERS

OCHO Shock Cord Panel

SPECIFICATIONS

H 6 in. x W 3 in. x D 0.25 in.
Weight: 116 oz

Color	ITEM #
Black	80-0260



OCHO Small Loop Panel

SPECIFICATIONS

H 7 in. x W 5.5 in. x D 0.25 in.
Weight: 183 oz

Color	ITEM #
Black	80-0261



OCHO Mesh Aid Pouch

SPECIFICATIONS

H 7 in. x W 4.5 in. x D 3 in.
Weight: 187 oz

Color	ITEM #
Black/Foliage	80-0259

Offers additional storage - compatible with any storage solution that has a loop panel attached. Rapid detach capability via red quick-pull tab.

OCHO IFAK



- Compact, lo-vis point-of-wounding IFAK designed for your duty/pistol belt
- Purposely designed to look like a non-tactical, accessory bag
- Equipped with a 3 in. belt loop attachment to accommodate a duty/pistol belt
- Constructed with rip stop nylon, elastic loop organizational storage and easy open zipper

SPECIFICATIONS

L 7 in. x H 2.25 in. x D 3.25 in.
Weight: 13 oz

ORDERING INFORMATION

Color	Level	ITEM #
Blue/Gray	Basic	80-0375
Black/Gray	Basic	80-0374
Blue/Gray	Advanced	80-0373
Black/Gray	Advanced	80-0372

OCHO IFAK BASIC CONTENTS:



OCHO IFAK ADVANCED CONTENTS:



Bag only option also available at NARescue.com

WALK


- Official vehicle life-saving kit of the U.S. Army
- For treatment of multiple trauma casualties & the immediate evacuation of one patient
- Includes full-size Talon® 90C quad-fold litter
- Pocket system holds self or buddy-aid medical supplies
- Back pack straps provide easy transport from the vehicle to the point of need
- Includes Hypothermia Prevention & Management Kit™ to help stabilize patients during transport
- Interfaces with any standard or non-standard TACEVAC asset

SPECIFICATIONS

H 22.75 in. x W 16.5 in. x D 11.5 in.

Weight: 29 lb 15 oz

WARRIOR AID & LITTER KIT CONTENTS:



ORDERING INFORMATION

COLOR	ITEM #	NSN #
Black	80-0097	----
Coyote	80-0089	6545-01-587-1199
Digital Camo	80-0059	6545-01-532-4962
OD Green	80-0060	----

Bag only option also available at NARescue.com

WALK ReSupply Kit (Item# 85-0020):
Contains WALK components for resupply - *excluding nylon bag & litter*



MCI WALK

- All-in-one trauma treatment & evacuation platform for treatment of multiple trauma patients
- Ideal for use as a vehicle trauma kit
- Choose from two different litter configurations: 8 non-rigid, disposable QuikLitters or 1 full-size, quad-fold, lightweight MedEVAC Litter
- Pocket system holds self- or buddy-aid point-of-wounding medical supplies
- Bag features back pack straps for easy transport from the vehicle to the point of need

SPECIFICATIONS

H 22.75 in. x W 16.5 in. x D 11.5 in.

Weight with 8 QuikLitters: 18 lb 12 oz

Weight with MedEVAC Litter: 29 lb 15 oz

Both red configuration options contain orange C-A-T® Tourniquets and the red option in the MedEVAC Edition (#80-0458) contains an orange MedEVAC Litter

MCI WALK - QUIKLITTER EDITION CONTENTS:



MCI WALK - MEDEVAC EDITION CONTENTS:



ORDERING INFORMATION

Color	Litter	ITEM #
Black	8 QuikLitters	80-0564
Red	8 QuikLitters	80-0565
Black	MedEVAC Litter	80-0457
Red	MedEVAC Litter	80-0458

Bag only option also available at NARescue.com



online video available

< scan to view



PATROL VEHICLE TRAUMA KIT

- Provides personnel with an easily accessible, quick to deploy trauma kit from their vehicle
- Modular storage system comes organized with an array of medical products ranging from basic first aid to major trauma care
- Equipment is packed in the kit to offer basic first aid and minor injury equipment on one side and major trauma equipment on the other
- Exterior includes MOLLE-style system to accommodate extra equipment as needed

SPECIFICATIONS

H 7 in. x W 12 in. x D 5 in.

Weight: 6 lb

ORDERING INFORMATION

LEVEL	ITEM #
Basic	80-0542
with Combat Gauze	80-0859

PATROL VEHICLE TRAUMA KIT CONTENTS:



CHOOSE:



OR:





LIFESAVING ► MISSION CRITICAL ► CUSTOM BUILT

North American Rescue® is **the market leader** for designing, producing, integrating and assembling tactical field medical sets and kits **supporting the services of the Department of Defense, federal agencies** as well as **state and local law enforcement**.

Over the past several years, our ability to provide rapid response products and services has aided in the delivery of critical battlefield medical care when and where it's needed most with an **end result of saving more lives**.

A comprehensive list of supply chain services affords a unique opportunity to North American Rescue's customers in the procurement, assembly, storage and management of their **complete medical and material logistics** requirements. NAR maintains a 75,000 square foot warehouse facility

Custom Solutions... for your unique requirements

which consists of 50,000 square feet for receiving and distribution staging and an additional 14,000 square feet solely dedicated to the assembly room where sets, kits and **custom medical solutions are built to order**.

Highly qualified logistics and warehouse staff ensure that orders, inventory and distribution processes **maintain the superior quality standards** associated with the North American Rescue brand. Our integral partnership with the South Carolina Vocational Rehabilitation **NISH Partner** is a pillar in the successful execution of these processes. ●

TRAINING

BLEEDING CONTROL SKILLS TRAINING KITS


- Provides the appropriate equipment to teach up to 20 students the necessary life-saving skills of bleeding control
- Packaged in a self-contained, wheeled, rigid case with internal dividers to keep the equipment organized and easy to inventory
- Skill set training addresses extremity bleeding, wound packing, pressure bandages, chest injuries and patient transport
- Corresponds with NAR's Public Access Bleeding Control Kits
- Available with or without Wound Packing Simulator

SPECIFICATIONS

Hard Storage Case: W 32 in x H 21 in. x D 12.5 in.
Weight without Wound Packing Simulator: 44.8 lb
Weight with Wound Packing Simulator: 50.9 lb

ORDERING INFORMATION

TRAINING LEVEL	ITEM #
Basic + Wound Packing Simulator	80-0880
Intermediate + Wound Packing Simulator	80-0878
Advanced + Wound Packing Simulator	80-0876
Basic	80-0881
Intermediate	80-0879
Advanced	80-0877




BASIC BLEEDING CONTROL TRAINING SKILLS KIT CONTENTS:



INTERMEDIATE BLEEDING CONTROL TRAINING SKILLS KIT CONTENTS:



ADVANCED BLEEDING CONTROL TRAINING SKILLS KIT CONTENTS:



C-A-T® TRAINER

Combat Application Tourniquet

- Fully functioning C-A-T® tourniquet colored blue to identify it as a training device
- Enables instructors and trainees hands-on training & tourniquet application practice
- Providing a distinct visual means to distinguish training devices from personal/operational life-saving devices

SPECIFICATIONS

H 6.5 in. x W 2.4 in. x D 1.5 in.
Weight: 2.7 oz

ORDERING INFO

ITEM #	NSN #
30-0033	6910-01-560-2972

See page 11 for Tactical and Rescue versions of the C-A-T®



HYFIN® TRAINER

Twin Pack

- Unique adhesive gel formulation does not damage mannequins or leave behind excessive residue
- Packaged with two training chest seals in one unique perforated package allowing students to open only one seal at a time

SPECIFICATIONS

Pkg: L 7.5 in. x W 4.5 in. x D 0.25 in.
Chest Seal: H 6 in. x W 6 in. x D 0.1 in.
Weight: 2.5 oz

ORDERING INFO

ITEM #
10-0032



See page 17 for other HyFin® Chest Seal product options



COMBAT GAUZE®

Moulage Trainer

- Allows medical trainers to thoroughly instruct personnel on the proper application of Z-Fold Combat Gauze®
- Completely inert (no hemostatic properties) - training purposes only
- Re-sealable packaging for multiple uses

ORDERING INFO

ITEM #
30-0063

SPECIFICATIONS

Packaged - H 6 in. x W 4.5 in. x D 1.25 in.
Gauze - L 4 yd x W 3 in.
Weight: 0.75 oz



HEMORRHAGE CONTROL TRAINING KIT

with Bio Feedback

- Training tool designed to teach effective wound packing skills with Immediate results for both trainer and student
- Verifies correct or incorrect technique

SPECIFICATIONS

Case: 23 in. x 8 in. x 15 in.
Weight: 17 lb 3.2 oz

ORDERING INFO

ITEM #
93-0019



HEMORRHAGE CONTROL TRAINING KIT

Combat Gauze®

- Comes with 10 Combat Gauze Trainers and 5 Combat Gauze LE Demo Samples
- Includes a "training leg" that features an entry/exit bullet wound and a deep laceration

SPECIFICATIONS

Case:
H 15 in. x W 23 in. x D 8 in.
Weight:
15 lb 8 oz

ORDERING INFO

ITEM #
93-0017



TCCC/TECC MODULAR TRAINING

Complete System - All 3 Modules

- Turnkey solution for teaching the PHTLS TCCC/TECC required tasks for completing the Tactical Combat Casualty Care (TCCC/TECC) course
- All the necessary tactical medical and casualty extraction/evacuation equipment needed for training centers and educational institutions
- Designed to equip up to 20 students
- Complete system consists of three self-contained hard storage case modules: Skill Stations Module, IFAKS Module and Extraction Module

SPECIFICATIONS

Each Case: W 30.75 in. x H 20.5 in. x D 11.6 in.

Module 1 Wgt: 53 lb • Module 2 Wgt: 59 lb • Module 3 Wgt: 45 lb

Complete System Weight: 157 lb

ORDERING INFO

ITEM #

85-0550



TCCC/TECC MODULE 1

SKILL STATIONS MODULE

- Train up to twenty (20) students at individual skill stations in the treatment of: *extremity hemorrhage, wound packing, chest seal application, needle chest decompression, junctional hemorrhage, basic airway adjuncts, surgical airway, splinting, hypothermia prevention/management, IV skills and eye injuries*
- Self-contained, wheeled, rigid case with internal dividers to keep the equipment organized and easy to inventory

ORDERING INFO

ITEM #

85-0519

SPECIFICATIONS

W 30.75 in. x H 20.5 in. x D 11.6 in.

Weight: 53 lb



QTY
20 C-A-T®
Trainer



QTY
20 ETD™
6 inch




QTY
20 Gauze
S-Rolled




QTY
20 Combat
Gauze® Trainer



QTY
1 JETT®
Device



QTY
5 Naso
Airway & Lube



QTY
1 Cyclone
BVM, Pocket



QTY
2 CricKit
Tactical



QTY
10 HyFin®
Trainer, Twin




QTY
5 ARS®
Needle, 14g



QTY
20 Saline
Lock Kit



QTY
5 Sharps
Shuttle



QTY
1 HPMK®
Kit



QTY
5 Cravat
Bandage



QTY
5 SAM
Splint




QTY
6 PES
Eye Shield

TRAINING

TCCC/TECC TRAINING MODULE 2 IFAKS MODULE

- Twenty pre-packed Training IFAKs for use by students during trauma lanes or training scenarios



ORDERING INFO

ITEM #

85-0520

SPECIFICATIONS

W 30.75 in. x H 20.5 in. x D 11.6 in.
Weight: 59 lb

QTY
20

EACH TCCC TRAINING IFAK CONTAINS:



TCCC/TECC TRAINING MODULE 3 EXTRACTION MODULE

- High-threat extraction devices and casualty evacuation platforms packaged in a self-contained, wheeled, rigid case that allows for easy transport and storage

ORDERING INFO

ITEM #

85-0521

SPECIFICATIONS

W 30.75 in. x H 20.5 in. x D 11.6 in.
Weight: 45 lb

TCCC EXTRACTION TRAINING MODULE 3 CONTAINS:



SKILL STATIONS & IFAK MODULE TRAIN 10 STUDENTS on Skill Stations & IFAKS

See Module 2 above for training IFAK contents

ITEM #

85-0740

SKILL STATIONS & IFAK MODULE CONTAINS:



ARMOR

M-A-R-C-H

KITS

PABC

TRAINING

TAC-EVAC

SUPPORT

RESCUE

LIFESTYLE

INDEX

TRAINING

IFAK TRAINER Individual First Aid Kit

- Pre-packaged training kit with all the essential equipment for training initial care
- The same point-of-wounding training kit found in the TCCC/TECC Training Module #2 (IFAKS Module), sold individually
- Easy-to-attach modular carrying bag can be worn as a drop-down leg rig, attached to a belt with a quick-detach buckle system or MOLLE'd to a vest
- ReSupply kit (contents only) also available in a clear plastic pouch



IFAK TRAINER CONTENTS:



ORDERING INFO

Style	ITEM #
Complete Kit	85-0493
ReSupply	85-0518

SPECIFICATIONS

H 8.25 in. x W 71 in. x D 5.5 in.
Weight: 2 lb 10 oz

I/O TRAINER

by EZ-IO®

- Compact design
- Red color indicates training only device

SPECIFICATIONS

H 4.5 in. x W 6.5 in. x D 2.5 in.
Weight: 15 oz

ORDERING INFO

ITEM #
30-0076



TRAINING NEEDLE KIT

T.A.L.O.N. Manual

- Train in manual I/O insertion at all seven optional sites

ORDERING INFO

ITEM #
30-0133



TRAINING NEEDLES

for the G3 Power Driver

- Designed for use with any EZ-IO® Power Driver

ORDERING INFO

Size	ITEM #
15 mm (Pediatric)	30-0082
25 mm (Adult)	30-0084
45 mm (Adult, large)	30-0086



TRAINER NEEDLE SET

EZ-IO® Manual

- Six different extremity sites for easily accessible insertion

ORDERING INFO

ITEM #
30-0079



G3 EDUCATOR KIT

by EZ-IO®

- Contents Include:
 - G3 Training Power Driver
 - EZ-IO® Training Needle Sets: 15mm Pediatric, 25mm Adult, 45mm Adult, Large
 - 4 Wristbands
 - 4 EZ-Connects
 - G3 Vascular Access Pack
 - 4 EZ-Stabilizers
 - Training Bones: Infant Tibia, Adult Proximal, Humerus, Adult Tibia
 - Training CD & softside Nylon Case



ORDERING INFO

ITEM #
30-0074

RANGER MEDIC HANDBOOK Fourth Edition


- Updated official medical reference of the 75th Ranger Regiment's Tactical Healthcare Professionals
- Fourth Edition contains more than 100 added pages of updated content (248 total pages) featuring updated Trauma Protocols, Tactical Medical Emergency Protocols, Pharmacology, Casualty Operations & Planning and Tactical Combat Casualty Care updates
- Ruggedized pocket-sized format made of 12 mil UltraFilm cover, tear-resistant, waterproof 6 mil paper pages and stay-open, flexible coil-spiral binding.
- Larger 8.5 in. x 11 in. textbook format also available with standard paper

ORDERING INFO

Style	ITEM #
Pocket Size	ZZ-0194
Textbook Size	ZZ-0426

SPECIFICATIONS

Pocket: H 7.5 in. x W 6.5 in. x D 1.25 in.
Weight: 1 lb 6 oz
Textbook: H 11 in. x W 8.5 in. x D 1 in.
Weight: 1 lb 2 oz



All net proceeds of the Ranger Medic Handbook are donated to the Sentinels of Freedom to help wounded warrior causes

FLIGHT PARAMEDIC Certification and Review Manual

- 154 page textbook, thoroughly covering ALL required areas of knowledge necessary to be successful on the Flight Paramedic Certification (FP-C) Exam
- Includes the latest CAMTS guidelines (9th Ed.)
- Easy-to-understand, bulleted format and full-color illustrations plus diagrams for the adult learner
- Author: Kyle Faudree

ORDERING INFO

ITEM #
ZZ-0467

SPECIFICATIONS

H 11 in. x W 8.5 in. x D 0.5 in.
Weight: 1 lb 2 oz
154 pages

A portion of the proceeds from the Flight Paramedic Study Guide is donated to wounded warrior foundations supporting Special Operations soldiers and their families



TACTICAL PARAMEDIC Certification & Practical Application

- 182 page training guide, thoroughly covering all materials required for achieving Tactical Paramedic Certification (TP-C), a standardized credential for TEMS providers created by the Board For Critical Care Transport Paramedic Certification (BCCTPC)
- Full-color excellent resource for medics looking to learn the intricacies of TEMS and passing the TP-C Exam
- Author: Kyle Faudree


ORDERING INFO

ITEM #
ZZ-0617

SPECIFICATIONS

H 10.88 in. x W 8.38 in. x D 0.38 in.
Weight: 1 lb 5 oz
182 pages

A portion of the proceeds from the Tactical Paramedic Manual is donated to wounded warrior foundations supporting Special Operations soldiers and their families



TAC·EVAC

MEDEVAC LITTER

- The ultimate, cost efficient casualty evacuation stretcher for First Responders operating in Active Shooter, Mass Casualty and Disaster incidents
- Compact, collapsible, quad-fold litter is easy to deploy and carry
- Two patient securing buckles, six IV-attachment points, ergonomically designed handles & patented joint configuration with locking leg crossbar

SPECIFICATIONS

Folded: H 21.5 in. x W 8.5 in. x D 7.5 in.
Open: H 84.5 in. x W 22.5 in. x D 4 in.
Weight: 13 lb 8 oz
Payload: 1,200 lb (544.3 kg)
Ground Clearance: 1.5 in.

ORDERING INFORMATION

Bed Color	Frame Color	ITEM #
Black	OD Green	60-0051
Orange	Black	60-0052



TALON® II 90C

Collapsible Handle Litter

- Standard 90 in. NATO-compatible evacuation platform
- Developed to U.S. Army requirements
- Advanced, innovative combat evacuation system featuring patented portable quad-fold design with retractable handles
- Slip-resistant, chemical-resistant rip-stop fabric tested and approved by Federal & private labs using live warfare agents

SPECIFICATIONS

Folded, handles collapsed: H 20.5 in. x W 8.5 in. x D 9.5 in.
Open, handles collapsed: L 78 in. x W 22.5 in. x D 6 in.
Open, handles extended: L 90 in. x W 22.5 in. x D 6 in.
Weight: 15 lb 6 oz (6.97 kg)
Payload: 1,200 lb (544.3 kg)
Ground Clearance: 1.5 in.
U.S. Patents: 6,842,923 & 5,598,592, Patent(s) Pending

ORDERING INFORMATION

ITEM#	NSN#
60-0002	6530-01-504-9051

WHEELED LITTER CARRIER

Lightweight with case



- Lightest wheeled litter cart on the market for single person transport to and from MedEVAC vehicles
- Locking clamp system keeps the litter firmly in place and allows for different litter widths
- Compatible with the Talon® II 90C, Raven® 90C, MedEvac & T3, plus all NATO litters
- Quick & easy to deploy or disassemble

SPECIFICATIONS

Deployed: L 56.29 in. x W 27.5 in. x H 36.2 in.
Case: L 48 in. x W 29.5 in. x H 12.5 in.
Weight: Carrier: 42 lb
Carrier with Case: 98 lb
Payload: 550 lb

ORDERING INFORMATION

ITEM#

60-0063



OSL™ LITTER

- Safely evacuate larger patients with 30% larger surface area compared to standard-sized evacuation platforms
- Includes 6 stirrups, 6 padded carrying loops & 3 patient securing straps

SPECIFICATIONS

Folded: H 45 in. x W 9 in. x D 8 in.
Open: L 90 in. x W 31 in. x D 6 in.
Weight: 22 lb 6 oz (10.2 kg) Ground Clearance: 1.5 in.
Payload: Floor: 1,600 lb (725.75 kg) Hanging: 1,200 lb (544.3 kg)

ITEM#	NSN#
60-0023	6530-01-552-7708

RAVEN® 90C

Bifold Litter


- The official litter of the U.S. Navy and certified safe-to-fly by the U.S. Air Force
- Slip-resistant, chemical-resistant, rip-stop fabric

SPECIFICATIONS

Folded: H 45.5 in. x W 7.5 in. x D 8 in.
Open: 90 in. x W 22.5 in. x D 6 in.
Weight: 16 lb 8 oz (7.5 kg)
Payload: 1,200 lb (544.3 kg)
Ground Clearance: 1.5 in.
U.S. Patents: 6,842,923 5,598,592

ORDERING INFORMATION

ITEM#	NSN#
60-0001	6530-01-432-5114



T3™ LITTER

Lightweight, Assault

- Ultra light, ultra strong litter weighing less than 12 lb
- Specifically designed for Special Operations Forces
- Overall length enables floor loading on an UH-60 helicopter

SPECIFICATIONS

Folded: H 21 in. x W 8.75 in. x D 7.25 in.
Open: L 81 in. x W 22.5 in. x D 6 in.
Weight: 11 lb 13 oz (5.34 kg)
Payload: 1,200 lb (544.3 kg)
Ground Clearance: 1.5 in.

ORDERING INFORMATION

ITEM#
60-0022



QUIKLITTER

Available in Standard and XL

- Lowest-cost emergency litter available
- Constructed of water-resistant coated woven polyethylene with ten large hand-holds
- Multi-use disposable evacuation platform rated up to 500 lb for safe, effective casualty transport & transfer
- QuikLitter™ XL is 12 in. wider and 6 in. longer than the original, features more hand-holds plus two built-in patient retention straps

SPECIFICATIONS

STANDARD: Payload 500 lb.
Litter Bed: L 74 in. x W 28 in.
Weight: 13.3 oz
XL: Payload 1,000 lb
Litter Bed: L 80 in. x W 40 in.
Weight: 2.2 lb

ORDERING INFO

Size	ITEM #
Standard	60-0061
XL	60-0065



EEL

Emergency Evacuation Litter

- Cost-effective, non-rigid casualty evacuation platform designed for fast, effective casualty evacuation in the tactical environment
- Rapid deploy handle allows quick removal from storage pouch
- Rugged, lightweight black nylon with 8 carry handles
- Compact, small-cube storage pouch attaches to the assault vest with easy PALS/MOLLE-style connectors

SPECIFICATIONS

Pouch: L 9.75 in. x W 7.5 in. x D 2.5 in.
Litter (Open): L 78 in. x W 30 in.
Weight: 1 lb 14 oz (0.57 kg)
Payload: 700 lb (317.5 kg)

ORDERING INFO

Color	ITEM #
Black	60-0047
Coyote	60-0048
Multi-Cam	60-0050
OD Green	60-0049



STINGRAY®

Poleless Litter

- Ideal for decontamination roller systems & moving wet patients
- Non-skid, flame-retardant monofilament polypropylene fabric that is resistant to mildew, acids, alkalis and penetrating liquids
- Tested by Federal & private labs using live warfare agents
- Includes six heavy-duty carrying handles & two casualty securing straps
- Easy to store and clean

SPECIFICATIONS

Folded: H 20 in. x W 12 in. x D 2 in.
Open: L 78 in. x W 24 in.
Weight: 13 oz (0.37 kg)
Payload: 1,000 lb (453.6 kg)

ORDERING INFO

ITEM #	NSN #
60-0004	6530-01-492-4079



TED

Tactical Extrication Device

- Designed to allow two people to carry a casualty to safety
- Lighter, more compact non-rigid evacuation platform designed specially for use in high-threat environments and limited spaces
- Constructed of heavy duty vinyl with 4 sturdy carrying handles
- Packaged in a sturdy, transport friendly bag with adjustable shoulder strap

SPECIFICATIONS

Folded: H 5 in. x W 10.25 in. x D 3.5 in.
Open: L 43.5 in. x W 19 in.
Weight: 1 lb 8 oz (0.68 kg)
Payload: 600 lb (272 kg)

ORDERING INFORMATION

Color	ITEM #	NSN #
Black	60-0005	6530-01-531-6808
OD Green	60-0006	6530-01-492-4114



LITTER STANDS

33 in. Standard

- Set of 2 lightweight stanchions for standard sized litters
- Constructed of a lightweight aluminum frame that conveniently folds flat for ease in transport & storage
- Quickly locks open to safely secure the litter in place

SPECIFICATIONS

H 34 in. x W 25.5 in. x D 6.92 in.
Extended: 33 in.
Base Width: 23.5 in.
Weight: 8 lb 8 oz

ORDERING INFORMATION

Unit	ITEM#	NSN#
Set of 2	60-0014	6530-01-523-3314



BACKREST

TALON® or RAVEN®

- Allows patients on a Talon® 90C or Raven® 90C litter to be placed at 90° or 45° angles
- Constructed of a lightweight aluminum frame with Nomex/ Kevlar-blend fabric covering



ORDERING INFORMATION

ITEM#	NSN#
60-0027	6530-01-594-7444

SPECIFICATIONS

H 31.5 in. x W 29.25 in. x D 16 in.
Weight: 2 lb 13 oz

OSL™ LITTER STANDS

for the OverSized Litter™

- Set of 3 lightweight aluminum stanchions accommodate our oversized OSL™ Litter
- Quickly opens and securely locks the litter in place
- Easily folds for transport & storage



SPECIFICATIONS

Folded: H 34 in. x W 25 in. x D 2.5 in.
Extended Height: 32 in.
Base Width: 23.5 in.
Weight: 28 lb 8 oz (12.93 kg)

ORDERING INFO

Unit	ITEM#
Set of 3	60-0029

TALON® CARRIER

Assault Litter Carrier

- Transport NAR's quad-fold litters (Talon, MedEvac & T3)
- MOLLE-style attachment points
- Includes thin back plates, drain hole, non-slip, adjustable shoulder straps & top-carrying handle

ORDERING INFORMATION

*NATO NSN

Color	ITEM #	NSN #
Black	60-0009	6530-01-531-6907
Coyote	60-0011	6530-01-531-6854
Digital Camo	60-0025	6530-01-643-9134*
OD Green	60-0010	6530-01-504-9055

SPECIFICATIONS

Collapsed: H 20 in. x W 9 in. x D 2 in.
Expanded: H 20 in. x W 9 in. x D 10 in.
Weight: 1 lb 7 oz



PORTABLE LIGHTING

LED FIELD MEDICAL LIGHT PAXLIGHT®

- Compact, lightweight precision lighting solution with the rugged engineering required for use in field hospitals.
- LED Technology with intense light fifteen times brighter than traditional halogen field lights.
- Patented 3-pod head eliminates shadow interference common with single point light sources.
- Energy efficient, cool white color with an exceptionally long life.
- Stays cool to the touch, projecting no heat on the patient.
- Locking joints make it fully adaptable and secure with adjustable mounting clamp for attachment to hospital bed, litter or litter stand.
- Includes a rugged storage case with wheels for easy transport.

SPECIFICATIONS

Light in case: L 47.2 in. x W 16.5 in. x D 6.7 in.
Light Head: 12.8 in. Diameter
Light Arm/Post: L 78 in. x 2.8 in. Diameter
Weight: Case & Light • 37 lb 8 oz
Light with Arm/Post • 11 lb 8 oz
Light Output: 50,000 lux (4590 foot candles)

ORDERING INFO

ITEM #	NSN #
60-0042	6530-01-549-1089



BATTERY BASE PAXLIGHT®

- Uninterrupted power supply to the PAXLight® for up to 8 hours
- Includes durable hard case

SPECIFICATIONS

Battery: H 12.5 in. x W 19 in. x D 17 in.
Pole: extends to 36 in.
Case: H 20 in. x W 32 in. x D 23.5 in.
Weight: Battery: 86 lb Battery & Case: 132 lb

ORDERING INFO

ITEM #
60-0043



DELUXE BATTERY BASE KIT • PAXLIGHT®

- Includes battery base, stackable transport case, portable charger & spare parts

SPECIFICATIONS

Battery: H 12.5 in. x W 19 in. x D 17 in.
Pole: extends to 36 in.
Case: H 20 in. x W 29.5 in. x D 25 in.
Weight: Battery: 86 lb
Battery & Case: 141 lb
Battery Output: 24VDC, 6.50A

ORDERING INFO

ITEM #	NSN #
60-0044	6530-01-588-8955



WHEELED BASE PAXLIGHT®

- Independent 5 caster wheel base with low center of gravity
- 13.5 lb counter weight for stabilization
- Includes convenient attachment point for power cord & AC adapter

SPECIFICATIONS

Base: W 38 in.
Assembled: H 38 in.
Weight: Base: 13.5 lb
Assembled: 22 lb

ORDERING INFO

ITEM #
60-0045



EZ LIGHT HANDLES PAXLIGHT®


- Sterile, disposable one-piece light handles to eliminate the time and labor of reprocessing
- Conveniently applies to the PAXLight handle via screw-on attachment for a solid feel

ORDERING INFO

Unit	ITEM #
Pack of 10	60-0046

SPECIFICATIONS

L 4.5 in. x W 1 in. x D 1 in.
Shield Diameter: 4.25 in.



FIELD HOSPITAL

MARK IV FIELD HOSPITAL BED

Lightweight solution

- Lightweight solution designed for patient comfort in the combat support or field hospital/treatment tent
- Durable & easy-to-open, supporting up to 850 lb of static weight
- Designed to be rapidly deployed with minimal training in less than 60 seconds by a single person
- Constructed of a bi-fold aluminum frame that is adjustable to 6 positions
- 3 in. polyurethane foam pad with water-resistant decontaminable cover easily cleans with 10% bleach solution or other standard hospital cleaning products
- Includes IV Pole and a removable table platform that rests on its adjustable bed rails

*extremely compact
for easy transport
and storage*

SPECIFICATIONS

Folded: L 42 in. x W 17 in. x D 32 in.
Opened: L 79.25 in. x W 33 in. x H 35.25 in.
Weight: 54 lb (Mattress and Frame)
Weight: 68 lb (Mattress, Frame and Accessories)
U.S. Patents Pending

ORDERING INFO

ITEM #	NSN #
66-0019	6530-01-624-7626

ARMOR

M·A·R·C·H

PABC

KITS

TRAINING

TAC·EVAC

SUPPORT

RESCUE

LIFESTYLE

INDEX

HYPERTHERMIA

POLAR SKIN™

- The latest cooling technology for effective Human Thermo-regulation
- Polar Skin Mini-Mesh,™ a unique smart textile with a hydrofunctional polymer that allows for liquid absorptive capacity 350% greater than base weight
- ArcticAqua,™ a proprietary solution that prevents rigidity at temperatures well below 0°F while providing the requisite moisture to support evaporative, convective and conductive cooling
- Adheres to current U.S. Army and National Athletic Trainers' Association protocols for heat illness treatment
- Protracted direct cooling of high blood flow regions

COMPREHENSIVE COOLING SYSTEM

- Highly-effective portable cooling system for the treatment of heat-related illnesses
- Scientifically shown to reduce core temperature by .09°F/minute in 100°F heat when used according to protocol

ORDERING INFO

ITEM

ZZ-0674

SPECIFICATIONS

L 11 in. x W 10 in. x H 4 in.

Weight: 8 lb 5.3 oz



COMPLETE CARE PACK

- Portable, highly insulated option
- Contains one Polar Skin Comprehensive Cooling System and one heavy-duty, leak & sweat proof cooler that keeps ice from melting for 24 hours in +100°F heat

ORDERING INFO

ITEM

80-0451

SPECIFICATIONS

W 14 in. x L 12 in. x D 7 in.

Weight: 10 lb 3.2 oz



COMPLETE CARE BACKPACK

- Available in either Single Patient (one comprehensive cooling system) or Double Patient (two cooling systems) configurations

ORDERING INFO

Size	ITEM #
Single	80-0543
Double	80-0544

SPECIFICATIONS

W 19.5 in. x L 15 in. x D 7.5 in.

Single Weight: 12 lb 13 oz

Double Weight: 21 lb 13 oz



Available in Standard, Large and Backpack COOLER

- Standard size accommodates one Polar Skin™ Comprehensive Cooling System and Large or Backpack Size accommodates two cooling systems

ORDERING INFORMATION

Color	SIZE	ITEM #
Red	Standard	ZZ-0682
Black	Standard	ZZ-0759
Black	Large	ZZ-0756
Black	Backpack	ZZ-0794

SPECIFICATIONS

Standard: L 14 in. x W 12 in. x H 7 in.
Weight: 1 lb 3.2 oz

Large: L 18.5 in. x W 13 in. x H 9.5 in.
Weight: 2 lb 13 oz

Backpack: L 19.5 in. x W 15 in. x D 7.5 in.
Weight: 3 lb 13 oz



ICE SHEETS Standard

- Principle cooling component of Polar Skin™ cooling systems
- Adheres to current National Athletic Trainers' Assoc and U.S. Army protocols for heat illness treatment
- Sheet Size: W 36 in. x L 30 in.



ORDERING INFORMATION

Unit	ITEM #	NSN #
Single Pack	ZZ-0668	6532-01-653-6054
Double Pack	ZZ-0669	----
Triple Pack	ZZ-0755	----

SPECIFICATIONS

Single: L 4.5 in. x W 7.75 in. x H 0.5 in. Wgt: 9.6 oz
 Double: L 9.5 in. x W 8 in. x H 0.5 in. Wgt: 1 lb 6.4 oz
 Triple: L 12.5 in. x W 7.5 in. x H 0.5 in. Wgt: 2 lb 3.2 oz

ICE SHEETS Large

- Same characteristics as the standard size Ice Sheets, but in a larger 36 in. x 60 in. size

ORDERING INFO

ITEM #	SPECIFICATIONS
ZZ-0758	L 8.5 in. x W 5.5 in. x H 1.5 in. Weight: 1 lb 8 oz

C2E ICE PACK BUNDLE

- Four Comfort Gel Packs
- Facilitates Conductive, Convective, and Evaporative (C2E) cooling of high blood flow (large artery) regions of the body



ORDERING INFO SPECIFICATIONS

ITEM #	SPECIFICATIONS
ZZ-0675	Indiv. Gel Pk: W 6 in. x L 10 in. Pkg: L 8 in. x W 9.5 in. x H 1.25 in. Weight: 4 lb 4 oz

COOLING COLLAR



- Two gel packs encased in integrated pouches
- Designed to align with the carotid artery

ORDERING INFO

ITEM #	SPECIFICATIONS
ZZ-0670	Indiv. Gel Pk: W 6 in. x L 10 in. Pkg: L 8 in. x W 9.5 in. x H 1.3 in. Weight: 4 lb 4 oz

SPLINTING

TC SPLINT

Telescoping Collapsible Splint

- Multi-capable device for splinting a wide variety of major and minor fractures
- Lightweight and compact for easy field carry with 6 easily adjustable interconnected telescoping poles and 4 wide elastic bands with hook and loop securing
- Includes clearly marked directions on the device as well as an easy-to-use, step-by-step instructional insert



online video available

< scan to view

ORDERING INFO

ITEM #
50-0039

SPECIFICATIONS

Pkgd: L 8 in. x W 5 in. x D 3 in.
 Open: W 9.5 in. x L 34 in.
 Weight: 1 lb 4 oz



COMBAT READY™

Tactical Body Tape

- Easy-to-use solution for reduced pain and swelling
- Precut strips designed for specific body area, can be worn for multiple days
- Allows the user to stay in the fight and complete the mission



MultiPack: One each of all 5 Combat Ready™ Tape

Average of **60% pain reduction** observed in over 1400+ patient encounters

ORDERING INFORMATION

Style	ITEM #	NSN #
Alpha	ZZ-0630	6510-01-659-7999
Bravo	ZZ-0631	6510-01-659-7998
Charlie	ZZ-0632	6510-01-659-8001
Delta	ZZ-0633	6510-01-659-8353
Echo	ZZ-0634	6510-01-659-8349
MultiPack	ZZ-0635	6510-01-659-3942

SPECIFICATIONS

Alpha, Bravo, Charlie, Delta, Echo Packaged: H 8 in. x W 5 in. x D 0.2 in. Weight: 1.5 oz

MultiPack Packaged: H 8 in. x W 5 in. x D 1 in. Weight: 7.5 oz

SPLINTING & ACCESSORIES

TACTICAL TRACTION SPLINT (TTS) Compact mechanical traction device

SPECIFICATIONS

Packaged: L 11 in. x W 5.5 in. x D 3 in.
Weight: 1 lb 5 oz

ORDERING INFO

Color	ITEM #	NSN #
Black	50-0015	6515-01-346-9186
Coyote	50-0033	----
Digital Camo	50-0034	----
OD Green	50-0032	----

- Relieve pain, increase patient comfort and reduce the risk of vascular injury or nerve damage in femur fractures.
- Lightweight and portable kit with clearly marked elastic straps for the proper fit to the thigh, knee & shin.



SAM® SPLINTS

Soft Shell 12 or 15 in.

SPECIFICATIONS

12 in: L 12 in. x W 4.25 in. x D 0.5 in.
15 in: L 15 in. x W 4.25 in. x D 0.5 in.

Weight: 2.7 oz
Weight: 3.2 oz

ORDERING INFO

Size Unit	ITEM #	NSN #
12 in. each	50-0006	6515-01-523-5403
15 in. each	50-0008	6515-01-523-5440

Also available in boxes of 20, please see NARescue.com

SAM® PELVIC SLING

for Pelvic Ring Fractures

- Force-controlled, ideal compression every time without the guesswork
- Rapid, simple, three-step application, usually in less than one minute

SPECIFICATIONS

H 725 in. x W 7.75 in. x D 4 in.
Hip Circum: 32 in. to 50 in.
Weight: 9.8 oz

ORDERING INFO

ITEM #	NSN #
50-0004	6515-01-509-6866

SAM® SPLINT II

- Versatile, malleable splint easily molds to fit most fractures
- Lightweight, comfortable and strong
- Waterproof, Reusable
- Easy to clean, Compact
- Illustrated step-by-step instructions included

SPECIFICATIONS

Pkgd: H 9.25 in. x W 4.25 in. x D 1 in.
Unfolded: L 36 in. x W 4.25 in.
Weight: 4.5 oz

ORDERING INFO

ITEM #	NSN #
50-1005	6515-01-494-1951



TRAUMA SHEARS

SPECIFICATIONS

SMALL: L 6.25 in. x W 2.75 in. x D 0.3 in.
Weight: 15 oz
LARGE: L 7.25 in. x W 3.7 in. x D 0.4 in.
Weight: 2 oz

ORDERING INFO

*NATO NSN

Size	ITEM #	NSN #
S 6.25 in.	ZZ-0064	6510-01-622-1624*
L 7.25 in.	ZZ-0063	6515-01-538-9276



SCISSOR LEASH™

- Salt water resistant
- Self cleaning
- Helps keep track of your shears

ORDERING INFO

ITEM #	NSN #
70-0005	6515-01-540-7226

SPECIFICATIONS

H 6 in. x W 16 in. x D 0.6 in.
Cable Length - 36 in.
Weight: 16 oz



PERSONAL PROTECTION

NITRILE GLOVES

BEAR CLAW®



- Form-fitting, ambidextrous design with textured fingertips
- Thick palm ensures durability
- Textured for improved wet and dry gripping
- Longer and thicker than standard nitrile gloves

Nitrile gloves available in Black or Tan colors and come either as pre-packaged pairs (kits) or in bulk (box)



BLACK TALON®



SPECIFICATIONS

Kit: L 4.75 in. x W 6 in. x D 4 in.
Weight: 1 lb 1 oz
Box: L 11 in. x W 5 in. x D 3.75 in.
Weight: 2 lb 4 oz

ORDERING INFORMATION

*NATO NSN

Size	Pkg Qty	ITEM #	NSN #
S	Kit 25 pair	ZZ-0246	6515-01-619-5350
S	Box 50 pair	70-0009	6515-01-600-6059*
M	Kit 25 pair	ZZ-0247	6515-01-619-5689
M	Box 50 pair	70-0010	6515-01-600-6057*
L	Kit 25 pair	ZZ-0248	6515-01-600-0262
L	Box 50 pair	70-0011	6515-01-600-6039*
XL	Kit 25 pair	ZZ-0249	6515-01-619-5367
XL	Box 45 pair	70-0012	6515-01-600-6052*

ORDERING INFORMATION

Size	Pkg Qty	ITEM #	NSN #
S	Kit 25 pair	ZZ-0218	----
S	Box 50 pair	70-0001	6515-01-515-0195
M	Kit 25 pair	ZZ-0217	----
M	Box 50 pair	70-0002	6515-01-515-0200
L	Kit 25 pair	ZZ-0216	6515-01-632-0059
L	Box 50 pair	70-0003	6515-01-515-0197
XL	Kit 25 pair	ZZ-0215	----
XL	Box 45 pair	70-0004	6515-01-515-0190

PES Polycarbonate Eye Shield

- Meets CoTCCC treatment protocols for penetrating eye trauma
- Convex eye shield accommodating either eye
- Translucent orange shell features perforations to allow air circulation

SPECIFICATIONS

3 in. x 2.38 in.
Weight: 0.2 oz



ORDERING INFO

Unit	ITEM #	NSN #
Single	30-0142	6515-01-647-5456
Box of 100	30-0143	6515-01-630-5551

FOX Eye Shield

- Convex rigid eye shield made of lightweight aluminum with perforations allowing air circulation
- Meets CoTCCC treatment protocols for penetrating eye trauma
- Available with or without garter



SPECIFICATIONS

3 in. x 2.38 in.
Weight: 0.2 oz

ORDERING INFO

Style	Unit	ITEM #	NSN #
No Garter	Box of 12	30-0144	6515-01-449-1016
Garter	Single	30-0055	----
Garter	Box of 50	30-5550	6515-01-598-1877

DIAGNOSTIC EQUIPMENT

FIELD CORPSMAN KITS

Portable ENT diagnostics

- Lightweight, durable ENT diagnostic sets for portable field care of the eyes, ears & nose
- A wide selection of high intensity lights & magnifiers for enhancing both day and night examinations



Items to the right included in Deluxe Kit only



SPECIFICATIONS

H 8 in. x W 4.25 in. x D 2.25 in.
Weight: Basic 6.2 oz
Deluxe 11.6 oz

ORDERING INFORMATION

Kit	ITEM #	NSN #
Basic	20-0001	6515-01-458-6057
Deluxe	20-0002	6515-01-458-6178

T2® TAG NAR

- Triage & treatment tag utilizing the S.A.L.T. Mass Casualty Triage Standard
- To aid in documenting critical patient information throughout the entire mass casualty response process

ORDERING INFORMATION


Unit	ITEM #	NSN #
20 Pack	20-0031	6515-01-624-7148
1000 Pack	20-0032	6545-01-592-3261

SPECIFICATIONS

L 5.69 x W 3.88 in.
Weight: 0.7 oz



INCLUDED WITH THE T2 TAG:



COMBAT CASUALTY CARD

- Triage card that meets all TCCC standards for casualty management

SPECIFICATIONS

H 6 in. x W 4.5 in. each
Weight: 0.22 oz ea | 5.5 oz 25-Pack

ORDERING INFORMATION

Unit	ITEM #	NSN #
25 Pack	ZZ-0262	6515-01-587-1700



COMBAT CASUALTY REFERENCE CARD

- Reference card for emergency management featuring TCCC Guidelines

SPECIFICATIONS

H 6 in. x W 5 in. x D 0.1 in.
Weight: 2.3 oz

ORDERING INFORMATION

ITEM#	NSN #
ZZ-0057	6515-01-537-4161



BURN DRESSINGS

BURNTEC® Hydrogel Dressings

- Sterile non-stick hydrogel dressings for treating 1st, 2nd & 3rd-degree burns (thermal, chemical, electric)
- Also treats other skin injuries, such as sunburns, abrasions, skin irritations, bruises and swelling (insect bites)
- Applied dressings provide a cool, soothing effect for up to 24 hours
- Easy to handle without spilling or damaging the hydrogel
- Permeable to water vapor and oxygen but impermeable to bacteria
- Biocompatible and hypo-allergenic; Reduces risk of hypertrophic scars

ORDERING INFO

Size	ITEM #	NSN #
2 in. x 2 in.	30-0114	6510-01-631-1913
2.5 in. x 5 in.	30-0115	6510-01-631-2080
4 in. x 4 in.	30-0110	6510-01-631-1883
5 in. x 5 in.	30-0116	6510-01-631-2138
5 in. x 10 in.	30-0117	6510-01-631-1907
8 in. x 8 in.	30-0118	6510-01-631-1897
8.5 in. x 11 in.	30-0119	6510-01-631-1890
4 in. x 16 in.	30-0111	6510-01-631-2124
8 in. x 16 in.	30-0112	6510-01-626-2990
16 in. x 24 in.	30-0120	6510-01-631-1886
10 in. x 10 in. Face	30-0121	6510-01-627-1834
12 in. x 16 in. Face	30-0113	6510-01-627-1855
16 in. x 24 in. Face	30-0122	----

BURN DRESSING KIT

- Convenient array of BurnTec® dressings in a handy carrying case.
- Kit contains (1) each of the following dressings:

2 in. x 2 in.	5 in. x 10 in.	8 in. x 16 in.
2.5 in. x 5 in.	8 in. x 8 in.	10 in. x 10 in. (face)
4 in. x 4 in.	4 in. x 16 in.	12 in. x 16 in. (face)
5 in. x 5 in.	8.5 in. x 11 in.	

ORDERING INFORMATION

ITEM #	NSN #
80-0361	6515-01-630-9177

SPECIFICATIONS

L 12 in. x W 15.5 in. x D 3.5 in.
Weight: 6 lb 4 oz



online video
available

< scan to view

MINOR BURN DRESSING KIT

- Variety of BurnTec® dressings in a handy carrying case.
- Kit contains the following:

Two 4 in. x 4 in. dressings
Two 4 in. x 16 in. dressings
One 8.5 in. x 11 in. dressing

ORDERING INFORMATION

ITEM #	NSN #
80-0415	6510-01-644-5427

SPECIFICATIONS

L 14 in. x W 8.5 in. x D 2.5 in.
Weight: 2 lb 13 oz

DRY STERILE BURN DRESSING

- Dry sterile, non adherent burn dressing & bandage
- Includes two No. 3 closed safety pins and rubber band individually packaged

SPECIFICATIONS

Bandage: 45 in. x 45 in. x 63 in.
Pkg: W 4.5 in. x H 4 in. x D 0.5 in.
Weight: 1 oz

ITEM #	NSN #
30-0162	6510-01-587-6579



FIRE ACCESSORIES

NITROSTRIKE™ Nitrogen Foam Systems

- Patented CAFS (Compressed Air Foam System) technology that brings increased capabilities and performance in a cost effective solution
- Unique blending system generates a more durable bubble matrix, allowing better adhesion to surfaces and more effective fire suppression
- Maximizes foam expansion rates and bubble durations for optimized fire fighting, meeting NFPA standards
- System allows the use of virtually all foams, including decontamination solutions, without the need to modify the base system*
- Reduces mechanical complexity by blending foaming agents at the manifold without moving parts
- Features the added durability required by extreme operational environments and reduces long-term maintenance costs
- Well suited for mobile and austere operations and ideal for tactical situations
- Single device configuration allows the use of AFFF (Aqueous Film Forming Foam) concentrate and decontamination solutions
- Offers a narrow profile for ease of storage and maximum maneuverability
- Compatible with Nitrogen or Compressed Air**

*AFFF freeze protect solutions should be utilized when conditions warrant

**Conversion fitting required for compressed air use
Fitting is sold separately



NitroStrike™ Backpack also available:

ITEM# 91-0014
See NARescue.com



NITROSTRIKE™ 2.5 GALLON

- Portable Fire Suppression System
- Allows 1 minute and 45 seconds of sustained fire fighting/decontamination up to a 25 feet
- Tested to UL standards and rated for 2A and 10B fires
- Simple to refill and ideal for mobile, austere and tactical operations

SPECIFICATIONS

H 23.5 in. x L 14.5 in. x W 8.5 in.
Weight: ~20 lb (dry) | 40 lb (full)

ORDERING INFORMATION

Case Color	ITEM #
Black	91-0018
Red	91-0000

COMPRESSED AIR ADAPTER for NitroStrike™

- Quickly converts Nitrogen CGA 580 cylinder to compressed air CGA 346 for filling NitroStrike™ units with compressed air
- Kit includes conversion adapter and universal pouch
- Optional mounting bracket available for permanent installs



SPECIFICATIONS

L 3.38 in. x Diameter 1.13 in.
Weight: 7.58 oz

ORDERING INFO

ITEM #
91-0016

FIRE ACCESSORIES

NITROSTRIKE™ 30 GALLON

- Allows 5 min of sustained fire fighting/decontamination at a distance of up to 75 ft.
- Tested to UL standards and rated for 10A and 40B fires
- Allows adjustments to the water content of the foam solution to better tailor the foam consistency (wet or dry) to your requirements
- Narrow profile for doorway access or limited storage, large wheels to navigate uneven surfaces and a protective frame for tank security

ORDERING INFORMATION

Hose Color	ITEM #
Black Hose	91-0017
White Hose	91-0001

SPECIFICATIONS

H 32 in. x L 40 in. x W 27.5 in.
Weight: ~330 lb (empty) | 560 lb (full)



NITROSTRIKE™ 60 GALLON

- Twice the duration (10 minutes) with the same shooting distance (up to 75 ft), discharge rate of volume, operation and use as the 30 gallon unit
- Equipped with caster wheels in front and flat free tires in rear to be pushed on flight lines, at airports, on warehouse floors, etc.
- Geared toward being mounted to the back of a truck, ATV, HUMVEE, etc.

SPECIFICATIONS

H 37 in. x L 55 in. x W 34 in.
Weight: ~500 lb (empty) | 980 lb (full)

ORDERING INFORMATION

Hose Color	ITEM #
Black Hose	91-0019
White Hose	91-0013



NAR FIRE BLANKET

- Made of inherently fire retardant 100% non-woven PCA fiber
- Resistant to flame, heat, acid, alkali, solvents and UV
- Rated to 350°C - will not carry flame, melt or drip
- Non-toxic, non-abrasive and displays anti-microbial properties

SPECIFICATIONS

Packaged: 20 in. x 15 in. x 4 in.
Opened: 78 in. x 70 in.
Weight: 4 lb 14 oz

ORDERING INFO

Color	ITEM #
Gray	70-0093
Red	70-0092



online video
available

< scan to view

TRANSPORT & STORAGE

ARMADILLO CASE Medication storage



Case only,
medication not included

SPECIFICATIONS

Ext: L 7.95 in. x W 5.1 in. x D 2.14 in.
Int: L 6.45 in. x W 3.87 in. x D 1.75 in.
Weight: 10.5 oz

- Compact solution for keeping critical medications safe, secure and easily accessible
- Hook & loop attachment system allows storage compartment to be tailored to hundreds of custom configurations
- Features Automatic Pressure Equalization Valve
- Hardware constructed of Marine grade 316 SS
- Case is water-resistant and crush-proof

ORDERING INFORMATION

* NATO NSN

Lid Style	ITEM #	NSN #
Black	ZZ-0081	8145-01-573-2533
Clear	ZZ-0259	6545-01-594-9399*

NOTE: Clear lid option does not include hook and loop system

Also Available: Half Inch Hook Discs with pressure sensitive adhesive
Set of 250 • Item# ZZ-0258 • NSN# 8145-01-617-9075*

M.U.L.E.

Multi-Use Luggage & Equipment Carrier

- Heavy-duty, water-resistant nylon platform with large compartment and mesh pockets on the inside lid for maximum storage of critical items
- Outside end pockets feature open channels for expanded storage capability
- Suitcase-style carrying handles on top and side; Telescoping handle and wheels for easy transport
- Hidden rucksack shoulder pads


ORDERING INFO

ITEM

ZZ-0212

SPECIFICATIONS

H 30.5 in. x W 18.5 in. x D 11 in.
Weight: 11 lb 10 oz Payload: 70 lb



I&A CASE

Injury & Ailment Case

- Superbly engineered medical equipment organization system
- Lid's clear pockets and hook & loop fasteners allow secure storage as well as superior inventory management
- Polymer O-ring provides dust and water-proof seal
- Automatic Gore-TEX® pressure equalization valve stops moisture from entering the case – also prevents vacuum-lock
- Double-throw latches engineered to tolerate up to 400 lb of pressure
- All-weather, sure-grip handle and reinforced padlock protectors



SPECIFICATIONS

H 8.5 in. x W 20.5 in. x D 17 in.
Weight: 13 lb 12 oz

ORDERING INFO

ITEM

ZZ-0222

TACTICAL LIGHTS

TOMAHAWK NIGHT VISION

by First Light®

- Right angle tactical light with rotating head and 3 button control panel for one-handed thumb actuation
- Provides White, IR, and Red+Green LED mix for highlighting blood during casualty treatment
- Featuring Integrated high intensity IR signaling beacon
- Includes momentary and three intensity levels of constant-on (all output colors)
- Can be used handheld or hands-free



ORDERING INFORMATION

Case Color	Task Light	ITEM #
Coyote	UV Task	70-0015
OD Green	UV Task	70-0112
Coyote	Blue Task	70-0111
OD Green	Blue Task	70-0110

SPECIFICATIONS

Output Wavelengths: White, IR, & Red+Green (Blood Highlighting), IR Signaling Beacon
 Maximum White Light Output: 425 lumens
 Max Output Runtime: 2 hours (ANSI Tested)
 Power Source: CR123 Batteries (2)
 Material: Aerospace-Grade Aluminum
 With mount: H 3.4 in. x L 2.4 in. x W 1.5 in. | Weight: (with mount & batteries) 8 oz

TOMAHAWK MULTI-COLOR

by First Light®

- Right angle tactical light with rotating head and 3 button control panel
- Provides White, Red, Green, and Red+Green LED mix for highlighting blood during casualty treatment
- Includes momentary and three intensity levels of constant-on



SPECIFICATIONS

Output Wavelengths: White, Red, Green, & Red+Green (Blood Highlighting)
 Maximum White Light: 425 lumens
 Max Output Runtime: 2 hours (ANSI Tested)
 Power Source: CR123 Batteries (2)
 Material: Aerospace-Grade Aluminum
 With mount: H 3.4 in. x L 2.4 in. x W 1.5 in.
 Weight (with mount & batteries): 7.2 oz

ORDERING INFORMATION

Case Color	Task Light	ITEM #
Coyote	None	70-0016
OD Green	None	70-0117
Coyote	UV Task	70-0116
OD Green	UV Task	70-0115
Coyote	Blue Task	70-0114
OD Green	Blue Task	70-0113

TORQ by First Light®

- Right angle tactical light, rotating head and 3 button control
- White, Red, Green, and Red+Green LED mix for highlighting blood
- Includes momentary and three intensity levels of constant-on (all output colors)
- Features integrated high intensity white tactical strobe and red/white/green flashing safety beacon
- Provides 155 lumens of maximum white light
- 10 hour Runtime on high output

ORDERING INFO

ITEM #	NSN #
70-0118	6515-01-646-2636

SPECIFICATIONS

Wavelengths: White, Red, Green, & Red+Green (Blood Highlighting)
 Max Runtime: 10 hours (ANSI Tested)
 With mount: H 3.9 in. x L 2.6 in. x W 1.3 in.
 Weight (with mount & batteries): 7.1 oz



PEGASUS ROPE®

WARHORSE™ SERIES

100% Innegra™ fiber braided core


- High quality, multipurpose rope
- Static rope provides reliability, high performance & great handling
- Excellent knot-ability with high strength-to-weight ratio
- Low stretch – wet, dry or in any other climatic condition
- Ultra-lightweight, strong Innegra™ core for reliable strength, durability & flotation
- Properties remain unchanged even when wet, dry or contaminated with chemicals
- Low dielectric and UV stable

SPECIFICATIONS

Diameter	9.1 mm	9.5 mm
Weight:	4.88 lb/100 ft	3.34 lb/100 ft
Tensile Strength:	5,670 lbf (22.7 kN)	5,100 lbf (22.7 kN)
Fibers (core/sheath):	Innegra™/Polyester	Innegra™/Polyester
Sheath Ratio:	63% Polyester	28% Polyester
Ratings:	NFPA 1983/2012 Throwline (Pending)	NFPA 1983/2012 Throwline (Pending)

ORDERING INFORMATION

Diameter	Color	Runner	ITEM #
9.1 mm	Black	Gray	90-0207
9.5 mm	Black	Gray Reflective	90-0206
9.5 mm	Yellow	Blue Reflective	90-0205
9.5 mm	White	Blue	90-0208



PEGASUS ROPE®

OLYMPUS™ SERIES

Engineered for rough terrain



ORDERING INFORMATION

Diameter	Color	Runner	ITEM #
9.5 mm	Gold	Black/Red	90-0401
11 mm	Gold	Black/Blue	90-0400

- Ideal for use in challenging and abusive environments
- Durable, lightweight and great handling offers reliability in all weather conditions
- Technora sheath adds strength and abrasion resistance
- Buoyant even with a 62% Technora sheath (9.5mm)

SPECIFICATIONS

Diameter	9.5 mm	11 mm
Weight:	3.82 lb/100 ft	4.99 lb/100 ft
Tensile Strength:	5,495 lbf (24.4 kN)	6,834 lbf (30.4 kN)
Fibers (core/sheath):	Innegra™/Technora	Innegra™/Technora
Sheath Ratio:	62% Technora	60% Technora

STALLION™ SERIES

Low visibility - high performance



ORDERING INFORMATION

Diameter	Color	Runner	ITEM #
9.5 mm	Black	Blue/Reflective	90-0103
9.5 mm	Black	Blue	90-0104

- 100 ft of this 9.5 mm rope offers 5,100 lb of strength under 3 lb – wet or dry
- Chemically resistant to acids, bases & common industrial chemicals
- Excellent strength-to-weight ratio and hydrophobic

SPECIFICATIONS

Diameter	9.5 mm
Weight:	2.9 lb/100 ft
Tensile Strength:	5,100 lbf (22.7 kN)
Fibers (core/sheath):	Innegra™/Innegra™
Sheath Ratio:	24% Innegra™

TRITON™ SERIES

Specially designed for water rescue

- 10% lighter & 250% stronger than polypropylene alone
- 100% Innegra™ fiber braided core, the lightest structural fiber available, provides additional buoyancy, low elongation and consistent strength

SPECIFICATIONS

Diameter	7.5 mm
Weight:	2 lb/100 ft
Tensile Strength:	3,050 lbf (13.6 kN)
Fibers (core/sheath):	Innegra™/Polypropylene
Sheath Ratio:	46% Polypropylene
Ratings:	NFPA 1983/2012 Throwline (Pending)

ORDERING INFORMATION

Diameter	Color	Runner	ITEM #
7.5 mm	Yellow	Blue Reflective	90-0303
7.5 mm	Red	Yellow	90-0305

ROPE RESCUE

R.A.C.E. SYSTEM Rapid Access Casualty Extraction

- Lightweight, multiple use casualty evacuation system enabling first responders/medical providers to ascend, descend, hoist & lower both rescuer & casualty to safety
- Combines 150 ft of lightweight durable 9.5 mm Olympus Pegasus Rope® with the latest in lightweight multipurpose hardware
- Capable of rigging 2:1, 3:1, 4:1, 5:1, 6:1, and 9:1 systems in different configurations for a wide variety of rescuer options with fully illustrated instruction booklet
- Front pocket allows for optional storage
- Carriage options include 4 Point, D-ring configurable modified side and/or backpack carry-style adjustable strap

R.A.C.E. SYSTEM CONTENTS:



ORDERING INFORMATION

Color	ITEM #
Coyote	90-0077
Multi-Cam	90-0063

SPECIFICATIONS

H 15 in. x W 10 in. x D 7.5 in.
Weight: 10 lb 14 oz

WATER RESCUE BELT THROW BAG

- Multipurpose water rescue solution designed for easy deployment and retrieval for rescuers
- Quick release reverse tapered throw bag for rapid packing of Triton rope and a slim profile for users
- 2 in. adjustable belt with quick release strap, 1.5 in. plastic cam buckle, and red toggle ball for easy ID



ORDERING INFO

ITEM #
90-0065

SPECIFICATIONS

Rope Bag: L 11 in. x W 4.5 in. x D 4 in.
Gear Pouch: L 5.5 in. x W 5 in. x D 2 in.
Total Weight: 2 lb 2 oz

50 FOOT THROW BAG

- Must-have necessity for water rescue
- Durable, easy-to-use rescue device with highly visible orange nylon bag and mesh panels that drain quickly for use in multiple rescue operations
- Packed with 50 ft of 7.5 mm Triton Pegasus Rope® rated at 3,050 lbf
- Triton rope is strong enough for multiple uses: as a throw line, tension diagonal, tow line and for pinned kayak recovery



ORDERING INFO

ITEM #
90-0064

SPECIFICATIONS

L 8.5 in. x W 4.5 in. x D 4.5 in.
Weight: 1 lb 5 oz

HELO TAGLINE KIT

Help control litter movement during hoist

- Help control litter movement during hoist operations – a requirement on all Army MEDEVAC helicopters
- Contains 250 ft of 7.5 mm Triton Pegasus Rope® featuring a highly visible polypropylene sheath with a super strong Innegra™ core
- Includes one locking carabiner, ten weak links (designed to break away if the litter is snagged during the hoist) and one stainless steel quick link


SPECIFICATIONS

L 11 in. x W 7 in. x D 7 in.
Weight: 5 lb 13 oz

ORDERING INFO

ITEM

90-0069



HIGH THREAT EXTRACTION KIT


Essential tools for rescue ops in high threat environments

- Facilitates immediate casualty extraction so that Care Under Fire casualty treatment may be performed as soon as possible
- Self-contained kit with all of the essential tools for tactical operators and rescue personnel
- Includes a C-A-T® tourniquet to immediately address extremity hemorrhage at the point of wounding
- Rugged nylon bag attaches to the duty belt as a drop down leg rig or to the assault vest/body armor with PALS/MOLLE-style connectors

SPECIFICATIONS

L 8.5 in. x W 8.75 in. x D 4.5 in.
Weight: 3 lb 6 oz

HIGH THREAT EXTRACTION CONTENTS:



ORDERING INFORMATION

Color	ITEM #
Black	80-0139
Coyote	80-0141
Digital Camo	80-0142
OD Green	80-0140



RESCUE EQUIPMENT

RAT STRAP

Rescue Assault Tether


- For head/feet first linear drags during shield-assisted, open-field/indoor-corridor rescues or single-rescuer extractions
- Innovative Daisy Chain configuration offers multiple reinforced attachment points

SPECIFICATIONS

Pkg: H 1 in. x W 4 in. x D 5.5 in.
Extended: L 62 in.
Minimum Break Strength: 2,702 lb
Weight: 12.7 oz

ORDERING INFO

ITEM #	NSN #
90-0002	6530-01-632-7872



EX-TRACT STRAP

Point of Wounding Extraction


- Point of wounding extraction
- Back pack carry
- Hasty harness for traditional and confined space lowering and hauling

SPECIFICATIONS

Pkg: H 18 in. x W 3.5 in. x D 1 in.
Extended: L 58 in.
Weight: 9.3 oz

ORDERING INFO

ITEM #
90-0082



DRAGON HARNESS

Under the vest harness

- Built-in handle system for immediate high threat extraction
- Operator harness to aid in the rescue of downed personnel without the need to holster or sling weapons
- Prevents the "ride up" of the casualty's ballistic vest, reducing exposed, unprotected body regions

SPECIFICATIONS

Height: 20 in. from waist to non-deployed handle
Length: 14 in. from harness attachment point to top of handle
Weight: 1 lb 1.6 oz.
Waist Size: 30 in. to 45 in.

ORDERING INFO

ITEM #
90-0003



DRAGON HANDLE

Tool offering head/feet first drag options


- Maintain an anatomically correct posture for rapidly extracting greater weight
- Decreases manpower needs in "zeroed in" locations
- Black Omega, non-locking "D" carabiner rated for 31kN (6,968.8 lb).

SPECIFICATIONS

H 6.3 in. x W 4.1 in. x D 2.25 in.
Weight: 11.3 oz
Device Length: L 42 in.
US Patent(s): Pending

ORDERING INFO

ITEM #	NSN #
90-0004	6530-01-598-9595



TACTICAL BAIL-OUT KIT

Rapid extraction of individuals or teams


- Integral for multiple hasty rescues including rapid casualty lowering & hoisting, rescuer throw bag and rappel line
- Facilitates a variety of anchor locations & rescue options

ITEM

90-0020

SPECIFICATIONS

L 10 in. x W 6 in. x D 4 in.
Weight: 2 lb 14 oz



HASTY HARNESS

Multi-purpose harness

- Resourceful, multi-purpose harness for unlimited applications & maximum versatility
- Indispensable gear for every webbing equipment toolbox
- Constructed of 22 ft loop of 1 in. Tech Tape (MBS: 4,400 lb) with opposing color thread to aid in visual safety inspection

ORDERING INFO

Color	ITEM #
Black	90-0001
Red	90-0044

SPECIFICATIONS

Loop: L 10 ft 11 in. x W 1 in.
Circumference - 21 ft 10 in.
Min. Break Str: Tech Tape (4,400 lb)
Sewn Sling (6,000 lb)



RESCUE EQUIPMENT

rockD CARABINER

by Rock Exotica®

Narrow asymmetric "D" shape machined frame with snag resistant keylock nose and large, ergonomic angled gate opening

SPECS

Gate Opening: 24 mm
Weight: 2.6 oz (73 g)

ITEM

ZZ-0543



BASIC ASCENDER

by Petzl®

Compact ascender with excellent grip for rope ascent on a fixed line or as progress capture in a hauling system

SPECS

Rope Size: 8-11 mm
Weight: 3 oz (85 g)

ITEM

ZZ-0554



GRIGRI 2 BELAY DEVICE

by Petzl®

Compact, ultra-light belay device allowing excellent control during the descent with assisted braking capability

SPECS

Rope Size: 8-11 mm
Weight: 6 oz (170 g)

ITEM

ZZ-0545



Mini Machined PULLEY

by Rock Exotica®

Mini pulley with recessed axle screw, machined sideplates and integrated axle for a stronger, compact, snag-free profile

SPECS

Rope Size: up to 11 mm
Weight: 2.9 oz (83 g)

ITEM

ZZ-0548



TriRig RIGGING PLATE

by Rock Exotica®

Lightweight, high performance 3 hole rigplate with numbered holes to assist with detailed rigging

SPECS

Break Strg: 36 kN
Weight: 1.8 oz (51 g)

ITEM

ZZ-0553



REVERSO 4

Belay/Rappel device by Petzl®

Regulates the amount of braking friction on the rope according to diameter and condition

SPECS

Weight: 2.08 oz (59 g)

ITEM

ZZ-0546



TIBLOC

Ascender/Rope Grab by Petzl®

Ultra-light Emergency Ascender/Rope Grab that can be used to build hauling systems or as a Prusik-replacement in self-rescue situations

SPECS

Rope Size: 8-11 mm
Weight: 3 oz (85 g)

ITEM

ZZ-0547



FOOTAPE

Adjustable Foot Loop by Petzl®

Adjustable nylon foot loop for attachment to the BASIC rope clamp for rope ascents

SPECS

Weight: 2.3 oz (65 g)

ITEM

ZZ-0551



ANNEAU SLINGS

by Petzl®

Polyester sling to set up a belay or an anchor

SPECS

W 19 mm
L 60 or 120 cm
Wt: 60 cm - 2 oz
120 cm - 3.5 oz

LENGTH

60 cm

ITEM

ZZ-0549

120 cm

ZZ-0550



SHIELD CARABINER

Hot forged "I-Beam" construction with patented anti-snag twisted wirelock gate

SPECS

Length: 100 mm
Weight: 1.1 oz

ITEM

90-0005



AERO SCREWGATE

Offset locking D carabiner with keylock nose that is the new standard for entry level lockers

SPECS

L 4 in x 2.46 in
Weight: 2.1 oz

ITEM

90-0006



REVOLVER CARABINER

Unique roller insert performs as an emergency pulley

SPECS

Gate Clear: 24 mm
Weight: 1.6 oz

ITEM

90-0007



RESCUE ACCESSORIES

RAPTOR™ MULTI-TOOL by Leatherman®

ARMOR

M-A-R-C-H

PABC

KITS

TRAINING

TAC-EVAC

SUPPORT

RESCUE

LIFESTYLE

INDEX



- High-performance multi-tool equipped with a variety of mission critical tools essential to emergency medical practitioners in the field

1. 420HC Stainless Steel Folding Medical Shears
2. Strap Cutter
3. Oxygen Tank Wrench
4. Ring Cutter
5. Carbide Glass Breaker
6. 5 cm Ruler
7. Holster with Clip

SPECIFICATIONS

Shear Length: 1.9 in.
Closed Length: 5 in.
Weight: 5.8 oz

ORDERING INFO

ITEM

70-0014

See our Treatment Support Section
on page 87 for the compatible Scissor Leash

SUPERTOOL 300 by Leatherman®

- Heavy-duty, all purpose multi-tool featuring 19 devices in one easy to carry configuration

SPECIFICATIONS

Closed: L 4.5 in.
Primary Blade: L 3.2 in..
Weight: 9.6 oz

ORDERING INFO

ITEM

ZZ-0542



1. Needlenose Pliers
2. Regular Pliers
3. Replaceable Wire Cutters
4. Hard-Wire Replaceable Cutters
5. Stranded-Wire Cutters
6. Electrical Crimper
7. Wire Stripper
8. 420HC Knife
9. 420HC Serrated Knife

10. Saw
11. Awl with Thread loop
12. 9 in. Ruler
13. Can Opener
14. Bottle Opener

15. Wood/Metal File
16. Phillips Screwdriver
17. Large Screwdriver
18. Medium Screwdriver
19. Small Screwdriver



MUT® by Leatherman®

- Tactical and practical multi-tool for military, law enforcement and civilian shooters - features 16 devices

1. Needlenose Pliers
2. Regular Pliers
3. Replaceable Wire Cutters
4. Hard-Wire Replac. Cutters
5. Stranded-Wire Cutters
6. Electrical Crimper
7. 420HC Knife
8. Saw
9. Replaceable Cutting Hook
10. Hammer
11. Bolt Override Tool
12. Bronze Carbon Scraper
13. #8-32 Cleaning Rod/Brush Adapter
14. Firearm Disassembly Punch
15. Carabiner/Bottle Opener Punch
16. Large Bit Drive

SPECIFICATIONS

Closed: L 5 in.
Primary Blade: L 3 in.
Weight: 11.2 oz

ORDERING INFO

ITEM

ZZ-0539

RESCUE ACCESSORIES

BENCHMADE® KNIVES & CUTTERS

8.50 in.
8.00 in.
7.50 in.
7.00 in.
6.50 in.
6.00 in.
5.50 in.
5.00 in.
4.50 in.
4.00 in.
3.50 in.
3.00 in.
2.50 in.
2.00 in.
1.50 in.
1.00 in.
0.05 in.



1 7 HOOK SAFETY CUTTER by Benchmade®

- Hook cutters designed to cut through a wide array of materials quickly and reliably
- More than 1.5 minutes faster than traditional trauma shears when cutting two articles of clothing

ORDERING INFO

ITEM #	NSN #
ZZ-0652	4240-01-543-9618

SPECIFICATIONS

Open: 4.3 in.
Weight: 1.5 oz

2 8 HOOK SAFETY CUTTER by Benchmade®

- 440C Cutting Surface
- Extensively tested and widely used by the US Military for medical and egress purposes

ORDERING INFO

ITEM #	NSN #
ZZ-0250	4240-01-568-3219

SPECIFICATIONS

Open: 6.5 in.
Weight: 1.5 oz

3 NAVY SAR SOML 112 SBK by Benchmade®

- Approved for Navy SAR with opposing bevel blade style with blunt tip and integrated rescue hook
- Opposing fixed bevel blade style with blunt tip and integrated rescue hook

ORDERING INFO

ITEM #
90-0078

SPECIFICATIONS

Overall: 8.02 in. Blade: 3.5 in.
Weight: 3.2 oz

4 TRIAGE 916 SBK by Benchmade®

- Dynamic triple utility tool featuring three essential rescue tools in a single compact design
- N680 high corrosion resistant, partially serrated, opposing bevel blade

ORDERING INFO

ITEM #	NSN #
90-0060	4240-01-598-4341

SPECIFICATIONS

Overall: 8.2 in. Blade: 3.5 in.
Weight: 5.1 oz

5 AUTOTRIAGE 9160 SBK by Benchmade®

- Packaged in high-strength aluminum with the speed & efficiency of AXIS® auto pull release and integrated safety
- Auto opening N680 highly corrosion resistant blade with modified drop-point

ORDERING INFO

ITEM #
90-0061

SPECIFICATIONS

Overall: 8.31 in. Blade: 3.58 in.
Weight: 5.75 oz



OPERATORS CAP NAR

- 100% cotton
- Adjustable strap
- 3.5 in. by 2.5 in. cut loop on the front designed for attaching patches (sold separately)
- 1 in. square of cut loop on the top for attaching an IR square

SPECIFICATIONS

One size fits most
Weight: 3.2 oz

ORDERING INFO

ITEM#

ZZ-0092



EMBROIDERED LOGO CAP NAR

- All American Cap
- Made in the USA
- Comfortable, attractive Fit-4-You™ design fits most wearers
- American flag on the front, NAR logo on back

SPECIFICATIONS

Fit-4-You™ S/M
Weight: 3.2 oz

ORDERING INFO

ITEM#

ZZ-0176



MEDIC CROSS PATCHES PVC

- Soft flexible PVC patch
- Backed with hook material for easy attachment and removal
- Ideal for use on medical kits such as the M-FAK and T.O.R.K.

ORDERING INFORMATION

Colors	ITEM #
Coyote & Black	ZZ-0797
Black & Gray	ZZ-0798
OD Green & Black	ZZ-0800
Black & Orange	ZZ-0799



SPECIFICATIONS

L 2 in. x W 2 in. x D 0.125 in.
Weight: 0.6 oz

IR FLAGS NAR

- Clearly visible under IR light
- Black hook fasteners on the back for easy hook-and-loop attachment and removal
- 2 flags per pack

SPECIFICATIONS

L 2 in. x W 2 in. x D 0.125 in.
Weight: 0.6 oz

ORDERING INFORMATION

Color	Stars	ITEM #
Full Color	L	ZZ-0093
Coyote	L	ZZ-0094
Coyote	R	ZZ-0235
OD Green	L	ZZ-0107



MED ID PATCH LUMINOUS

Glow in the Dark

- Soft flexible custom PVC with hoop and loop backing
- Ultimate glow-in-the-dark identification technology
- 1 in. block letters
- 2 med patches per package

SPECIFICATIONS

Each: L 3 in. x W 2 in. x D 0.125 in.
Weight: 0.4 oz

ORDERING INFO

ITEM#

ZZ-0121



IR SQUARES NAR

- Infrared patch printed on honeycomb reflective material
- Black hook fasteners on the back for easy attachment
- Clearly visible under IR light illumination
- Six squares per pack

SPECIFICATIONS

Each: 1 in. x 1 in.
Weight: 0.2 oz.

ORDERING INFO

ITEM#

ZZ-0108



LIFESTYLE/APPAREL

NAR T-SHIRTS *from* RANGER UP



VINTAGE LOGO

ORDERING INFO

ITEM#

TS-1



GRUNGE EAGLE

ORDERING INFO

ITEM#

TS-2



AMERICAN FLAG & EAGLE

ORDERING INFO

ITEM#

TS-3



LOGO FLAG

ORDERING INFO

ITEM#

TS-4



EKG LOGO

ORDERING INFO

ITEM#

TS-5



BEFORE YOU ASK

ORDERING INFO

ITEM#

TS-6



COMBAT MEDIC

ORDERING INFO

ITEM#

TS-7



CALL ME DOC

ORDERING INFO

ITEM#

TS-8

FEATURES

- 100% Soft Cotton
- Pre-Shrunk Shirts
- Set-in Ribbed Collar
- Double Needle Hems
- Available in S-XXL
- Cutting Edge Design

For individual colors/logo combinations please see NARescue.com

CUSTOMER SUPPORT

- **COMMITMENT:** We pride ourselves in providing an exceptional customer experience beginning with your initial inquiry up to order placement and beyond. Our customer service staff members are trained by certified instructors on all NAR products and are standing by to help with any questions you may have
- **SUPPORT:** Our in-house technical support professionals consist of an elite group of Nationally Registered Paramedics and Special Operations Medics from military, law enforcement, tactical law enforcement, fire-rescue and EMS backgrounds

ORDERING

- **ONLINE:** www.NARescue.com
GOVERNMENT PRICING: automatically available ONLINE for all military customers with a .mil e-mail address! Register your account in just minutes by creating an online account using your .mil e-mail address
ESTABLISHED CUSTOMERS:
Want to order online with special pricing? Contact Customer Service at 888.689.6277 to have your online account set up
NEW CUSTOMERS: Simply create an online account and you will have access to customer specific pricing – anytime, anywhere you are logged in to NARescue.com
QUOTE REQUESTS: Online accounts have no waiting to receive your custom pricing to place an order or, we will gladly accept your submittal for a pricing quote request
- **CONTRACT:**
Tax ID: 27-1024029
Cage Code: 06ST7
DUNS: 83-24-26782
CCR Registered: Yes
ORCA Reg: Yes
ECAT: dmmonline.dscp.dla.mil
Prime Vendor: dmmonline.dscp.dla.mil/primevendor/primehome
- **EMAIL:** info@narescue.com
- **PHONE:** 888.689.6277
- **ADDRESS:** 35 Tedwall Ct | Greer, SC
- **FAX:** 864.675.9880

SHIPPING

- In-stock merchandise ships within two business days
- Need it ASAP? Ask Customer Service about overnight or second day delivery

STAY informed

and KEEP IN TOUCH



Keep up with the latest North American Rescue® news and happenings!

Join us on Facebook:
facebook.com/NARescueLLC



Providing medical solutions and celebrating the life-savers who use them!

Join us on Instagram:
instagram.com/northamericanrescue



Keep up with North American Rescue® with 140 characters or less!

Join us on Twitter:
twitter.com/NARescue



View North American Rescue® instructional and informational videos!

Join us on You Tube:
youtube.com/NARescue01



Visit North American Rescue® online to sign up for our newsletter!

Sign up at:
www.NARescue.com



Office Hours: M-F 8:00am - 5:00pm EST



900 535 295 | 688 66 089 | 944 008 847 | comercial@iesmedical.es | iesmedical.es

PRODUCT index by name

A

Agent Armor Carrier System	6
Agent RTF Kit	6

Airway Management

AMBU® RDIC	4, 16
AMBU® SPUR II BVM	15
Cyclone® Improved Pocket BVM	15
Laryngoscope Blades, NAR	16
King LTS-D™	16
Laryngoscope Set, Special Ops	16
Nasopharyngeal Airway	15
Tactical Crickit®	15
Tactical Suction Device, NAR	15
Tracheal Hook, NAR	15

AMBU® RDIC	16
AMBU® Spur BVM	4, 16
Ankle Trauma Holster	27
Ankle Trauma Kit	27

Apparel

Embroidered Logo Cap	78
Flags	78
IR Squares	78
Medic Cross Patch	78
Operators Cap	78
Armadillo Case	68
ARS® Needle Decompression Kit	18

B

Bail-Out Kit, Tactical	74
------------------------------	----

Ballistic Accessories

Gear Bags	4, 9
Identification Safety Banners	5, 9
Rescue Task Force Placards	9
Soft Litter Holder	9

Ballistic Protection

Agent Armor Carrier System	6
Agent RTF Kit	6
Batlskin® Viper A1	10
Batlskin® Viper A3	10
Batlskin® Viper P2	10
DKX Level III Rifle Plate	9
IPRO® by NAR	10
OTW Goggles	10
Rescue Responder Side Armor	8
Rescue Task Force Vest Kit	7
Tactical Responder Vest	8
Batlskin® Viper Helmets	10
Bear Claw® Gloves	63
Belt Trauma Kit	30

Benchmade® Knives

7 Hook Safety Cutter	77
8 Hook Safety Cutter	77
Autotriage 9160 SBK	77
Navy SAR SOML 112 SBK	77
Triage 916 SBK	77

Black Talon® Gloves	63
Bleeding Control Training Kits	5, 48
BOA® IV Constricting Band	19

Books/Manuals

Flight Paramedic	53
Ranger Medic Handbook	53
Tactical Paramedic	53
Burn Dressing, Dry Sterile	65
BURNTec® Hydrogel Dressings	65

C

C3 - Concealed C•A•T® Carrier	12
-------------------------------------	----

Canine (K-9)

K-9 Handler IFAK	39
K-9 Tactical Field Kit	39

C•A•T® Tourniquet

Ankle Trauma Holster	27
C3 - Concealed C•A•T® Carrier	12
C•A•T® Holder	12
Rigid C•A•T® TQ Case	5, 12
Rigid C•A•T® TQ Case & Cover	12
TQ 911 C•A•T® Holder	12

C•A•T® Holder	12
---------------------	----

Celox™	13
--------------	----

Chitogauze®	13
-------------------	----

Circulation

BOA® IV Constricting Band	19
G3 Power Driver	20
Manual Extremity Needle Set	20
Needleless Saline Lock Kit, NAR	19
Power Driver Needle Set	20
Raptor® IV Securing Device	20
Saline Lock Kit	19
Sharps Shuttle	20
T.A.L.O.N.™ Manual Needle Set	20

Combat Casualty Card	64
----------------------------	----

Combat Casualty Reference Card	64
--------------------------------------	----

Combat Gauze®	13
---------------------	----

Combat Gauze® LE	13
------------------------	----

Combat Gauze® Trainer	49
-----------------------------	----

Combat Ready™ Tactical Tape	61
-----------------------------------	----

Compressed Gauze, NAR	14
-----------------------------	----

CORE	31
------------	----

Cravat	14
--------------	----

Crickit®	15
----------------	----

Crisis Incident Response Kit	40
------------------------------------	----

Custom Kit Solutions	47
----------------------------	----

Cyclone® BVM	15
--------------------	----

D

DKX LEVEL III Rifle Plate	9
---------------------------------	---

Dragon Handle	74
---------------------	----

Dragon Harness	74
----------------------	----

E

Eagle IFAK	28
------------------	----

EEL, Emergency Evacuation Litter	56
--	----

Emergency Survival Wrap	5, 21
-------------------------------	-------

EMS Rapid Deployment Kit	36
--------------------------------	----

ETD™

Emergency Trauma Dressings	4, 14
----------------------------------	-------

Responder ETD™	14
----------------------	----

Eye Protection

Fox Eye Shield	63
----------------------	----

IPRO® Tactical Goggles	10
------------------------------	----

Outside the Wire Goggles	10
--------------------------------	----

Polycarbonate Eye Shield, PES	63
-------------------------------------	----

EZ-IO® G3 Educator Kit	52
------------------------------	----

EZ-IO® G3 Power Driver	20
------------------------------	----

EZ-IO® I/O Trainer	52
--------------------------	----

EZ-IO® Extremity Needle Set	20
-----------------------------------	----

EZ-IO® Power Driver Needle Set	20
--------------------------------------	----

EZ-IO® Training Needles	52
-------------------------------	----

F

Field Corpsman Kit	64
--------------------------	----

Field Hospital

Battery Base PaxLight®	58
------------------------------	----

Deluxe Battery Base PaxLight®	58
-------------------------------------	----

EZ Light Handles PaxLight®	58
----------------------------------	----

LED Field Medical Light	58
-------------------------------	----

Mark IV Field Hospital Bed	59
----------------------------------	----

Wheeled Base PaxLight®	58
------------------------------	----

Fire Blanket	67
--------------------	----

First Light®

Tomahawk	69
----------------	----

Torq	69
------------	----

Fox Eye Shield	63
----------------------	----

G

Gauze

Celox™	13
Chitogauze®	13
Combat Gauze®	13
Combat Gauze® LE	13
Combat Gauze® Trainer	49
Compressed Gauze, NAR	14
QuikClot Bleeding Control	13
S-Rolled, NAR	14

Gear Bags	9
-----------	---

Gloves

Bear Claw® Gloves	63
Black Talon® Gloves	63

Goggles

IPro® Tactical Goggles	10
Outside the Wire Goggles	10

H

Hasty Harness	74
Helmets, Batlskin® Viper	10
Helo Tagline Kit	73
Hemorrhage Control Training	49
High Risk Warrant Casualty Kit	32
HPMK®	21
HRS® Heat Reflective Shell	21

HyFin® Chest Seals

HYFIN® Trainer	49
HyFin® Vent Compact	17
HyFin® Vent Individual	17
HyFin® Vent Occlusive	17
HyFin® Vent Twin Pack	17

Hyperthermia	60
--------------	----

Hypothermia Prevention

Emergency Survival Wrap	21
HPMK®	21
HRS® Heat Reflective Shell	21
Ready Heat™ Blanket	21

I

Identification Banners	9
IFAK Trainer	52
Individual Aid Kit	40
Injury & Ailment Case	68
Innegra™	71
IPOK	31
IPOK Holder	31
IPro® Tactical Goggles	10

J

JETT® Junctional Tourniquet	11
JETT® Pack	11

K

K-9 Handler IFAK	39
K-9 Tactical Field Kit	39
King LTS-D™	16

Kits

Agent RTF Kit	6
Ankle Trauma Holster	27
Ankle Trauma Kit	27
Belt Trauma Kit	30
Bleeding Control Skills Training	48
Burn Dressing Kit	65
C.O.R.E. Kit	31
Crisis Incident Response Kit	40
Eagle IFAK	28
EMS Rapid Deployment Kit	36
Field Corpsman Kits	64
First Responder Go Bags	34
Helo Tagline Kit	73
High Risk Casualty Warrant Kit	32
High Threat Extraction Kit	73
Individual Aid Kit	40
IPOK	31
IPOK Holder	31
K-9 Handler IFAK	39
K-9 Tactical Field Kit	39
Mass Crisis Incident Kit	40
MCI WALK	45
Medic Trauma Pack	35
Medic Trauma Sheet	42
M-FAK	26
Minor Burn Dressing Kit	65
NAR-4™ Aid Kit	37
NAR-4™ Chest Pouch	30
OCHO IFAK	43
Operator BLS IFAK	35
Patrol Vehicle Trauma Kit	5, 46
Polar Skin	60
Public Access Individual Kits	22
Public Access Packs	4, 5, 23-25
R.A.C.E. System	72
Range Trauma Aid Kit, Hard Case	38
Range Trauma Kit	38
Rapid Response Kit, LE Edition	34
Rapid Response Kit, RTF Edition	34
Rescue Task Force Vest Kit	7
RIG Series Eagle IFAK	29
SRO Crisis Response Kit	41
SRO Individual Response Kit	41
S.T.O.R.M. Operator IFAK	32
Tactical Bail-Out Kit	74
Tactical Rapid Deployment	36
T.O.R.K.	33
Unconventional Medic	47
Visor Med Pouch	31
WALK ReSupply Kit	44
W.A.L.K. Warrior Aidd & Litter Kit	44

L

Laryngoscope Blades, Fiber Optic	16
Laryngoscope Set, Special Ops	16

Leatherman®

MUT®	76
Raptor™ Multi-Tool	76
Supertool 300	76
Light, Field Medical LED	58
Lights, Tactical	69

Litter Accessories

Backrest, Talon or Raven	57
--------------------------	----

Litter Stands	57
OSL™ Litter Stands	57
TALON® Carrier	57

Litters

EEL, Emergency Evacuation Litter	56
QuikLitter	56
QuikLitter XL	56
Soft Litter Carrier	4
STINGRAY®	56
TED, Tactical Extrication Device	56
Litter Stands	57

M

Mark IV Field Hospital Bed	59
Mass Casualty Incident WALK	45
Mass Crisis Incident Kit	40

Massive Hemorrhage

C-A-T®	11
Celox™	13
Chitogauze®	13
COMBAT GAUZE®	13
Combat Gauze® LE	13
Compressed Gauze	14
Cravat (Triangular Bandage)	14
ETD™	4, 14
JETT®	11
QuikClot® Bleeding Control	13
Responder ETD™	14
S-Rolled Gauze	14
Tactical Compression Wrap	14
MCI WALK	45
MedEVAC Litter	54
Medic Cross Patch	26, 78
Medic Trauma Pack	35
Medic Trauma Sheet	42
M-FAK	26
MKII Tactical Responder Vest	8
M.U.L.E.	68

N

NAR-4™ Aid Kit	37
NAR-4™ Chest Pouch	30
Nasopharyngeal Airway	15
Needless Saline Lock Kit, NAR	19
Nitrostrike™	66

O

Ocho Bags & Accessories

Ocho Backpack	42
Ocho Briefcase	42
Ocho IFAK	43
Ocho IFAK/Accessory Bag	42
Ocho Mesh Aid Pouch	42
OCHO Shock Cord Panel	42
Ocho Shoulder Bag	42
Ocho Small Loop Panel	42
Ocho IFAK	43
Operator BLS IFAK	35
OSL™ Litter	55
OSL™ Litter Stands	57
Outside the Wire Goggles	10

P

Patches

Flags	78
IR Squares	78
Medic Cross Patch	26, 78
Med ID Patch	78
Patrol Vehicle Trauma Kit	5, 46
PaxLight®	58
Pegasus Rope®	70-71
Placards, Rescue Task Force	9

Polar Skin™

C2E Ice Pack Bundle	61
Complete Care Backpack	60
Complete Care Pack	60
Comprehensive Cooling System	60
Cooler	60
Cooling Collar	61
Ice Sheets, Large	61
Ice Sheets, Standard	61
Polycarbonate Eye Shield, PES	63

Public Access

5-Pack Bleeding Control	23
8-Pack Bleeding Control	23
Bleeding Control Skills Training	48
Individual Kits	22
Nylon Bag Individual Kits	22
Vacuum Sealed Individual Kits	22
Wall Station, Clear Polycarbonate	24
Wall Station, Semi-Recessed	5, 25

Q

Quick Detach Panel, NAR	33
QuikClot® Bleeding Control	13
QuikLitter	56

R

R.A.C.E. System	72
Range Trauma Aid Kit, Hard Case	38
Range Trauma Kit	38
Rapid Response Kit, LE Edition	34
Rapid Response Kit, RTF Edition	34
Raptor® IV Securing Device	20
RAT Strap	74
RAVEN® 90c	55
Raven® Backrest	57
Ready Heat™ Blanket	21

Reference

Combat Casualty Reference Card	64
Flight Paramedic	53
Ranger Medic Handbook	53
Tactical Paramedic	53

Rescue Accessories

Aero Screwgate	75
Anneau Slings	75
Basic Ascender	75
Benchmade® Knives/Cutters	77
Dragon® Harness	74
FOOTAPE	75
GRIGRI 2 Belay Device	75

Hasty Harness	74
Helo Tagline Kit	73
H.I.D.E. Suspenders	74
High Threat Extraction Kit	73
Leatherman® Tools	76
Mini Machined Pulley	75
R.A.C.E.	72
RAT Strap	74
Reverso 4	75
Revolver Carabiner	75
rockD Carabiner	75
Shield Carabiner	75
Tactical Bail-Out Kit	74
Throw Bag, 50 foot	72
Tibloc Ascender	75
TriRig Rigging Plate	75
Water Rescue Belt Throw Bag	72
Rescue Responder Side Armor	8
Rescue Task Force Placards	9
Rescue Task Force Vest Kit	7

Respiratory Support

ARS® Needle Decompression	18
HyFin® Occlusive Individual	17
HyFin® Twin Pack Trainer	17, 49
HyFin® Vent Compact	17
HyFin® Vent Individual	17
HyFin® Vent Twin Pack	17
Rigid C-A-T® TQ Case	5, 12
Rigid C-A-T® TQ Case & Cover	12
RIG Series Eagle IFAK	29

Ropes

Pegasus Rope®	70
Olympus™ Series	71
Stallion™ Series	71
Triton™ Series	71
Warhorse™ Series	70

S

Saline Lock Kit, NAR	19
Saline Lock Kit, Needleless	19
SAM® Splints	62
Scissor Leash™	62
Sharps Shuttle	20

Shears

Scissor Leash™	62
Trauma Shears	62
Side Armor, Rescue Responder	8
Social Media	80
Soft Litter Holder	9

Splinting

SAM® Pelvic Sling	62
SAM® Splint II	62
SAM® Splints	62
Tactical Traction Splint	62
Telescoping Collapsible Splint	61
SRO Crisis Response Kit	41
SRO Individual Response Kit	41
S-Rolled Gauze, NAR	14
Stingray® Poleless Litter	56
S.T.O.R.M. Operator IFAK	32
Survival Wrap, Emergency	21

T

T2® Tag	64
T3™ Lightweight Assault Litter	55
Tactical Bail-Out Kit	74
Tactical Body Tape, Combat Ready	61
Tactical Compression Wrap	14
Tactical Crickit®	15
Tactical Lights	69
Tactical Rapid Deployment Kit	36
Tactical Responder Vest	8
Tactical Suction Device, NAR	15
Tactical Traction Splint	62
TALON® Backrest	57
TALON® Carrier	57
TALON® II 90C	54
T.A.L.O.N. Manual Training	52
T.A.L.O.N. Manual Needle Set	20

Training

Bleeding Control Skills Training	48
Extraction Module, Module 3	51
IFAKS Module, Module 2	51
IFAK Trainer	52
Skill Stations & IFAKS, 10 Students	51
Skill Stations Module, Module 1	50
TCCC Complete Training System	50
TC Splint	61
TED, Tactical Extrication Device	56
Throw Bag, 50 foot	72
T.O.R.K.	33
TQ 911 C-A-T® Holder	12
Tracheal Hook, NAR	15
Trauma Shears, NAR	62
T-Shirts	79

U

Unconventional Medic	47
----------------------	----

V

Vehicle Trauma Kit, Patrol	46
----------------------------	----

Videos, Online

ARS Needle Decompression Kit	18
BOA IV Constricting Band	19
Burntec	65
Combat Application Tourniquet®	11
Fire Blanket, NAR	67
HyFin Chest Seals	17
JETT	11
MCI WALK	45
Public Access	22
Telescoping Collapsible Splint	61
Viper Helmets, Battskin	10
Visor Med Pouch	31

W

WALK®	44
WALK®, Mass Casualty Incident	45
Water Rescue Belt Throw Bag	72
Wheeled Litter Carrier	55

PRODUCT

index

by item number

10

10-0001.....	16
10-0015.....	17
10-0017.....	15
10-0018.....	15
10-0019.....	15
10-0020.....	15
10-0025.....	16
10-0026.....	16
10-0027.....	16
10-0029.....	17
10-0032.....	49
10-0037.....	17
10-0042.....	17
10-0043.....	15
10-0047.....	16

20

20-0001.....	64
20-0002.....	64
20-0031.....	64
20-0032.....	64

30

30-0001.....	11
30-0003.....	14
30-0006.....	20
30-0009.....	19
30-0012.....	14
30-0013.....	19
30-0014.....	20
30-0016.....	14
30-0023.....	11
30-0031.....	14
30-0032.....	14
30-0033.....	11,49
30-0044.....	12
30-0046.....	19

30-0052.....	14
30-0053.....	12
30-0055.....	63
30-0057.....	12
30-0059.....	12
30-0063.....	49
30-0071.....	19
30-0074.....	52
30-0075.....	20
30-0076.....	52
30-0077.....	20
30-0078.....	20
30-0079.....	52
30-0081.....	20
30-0082.....	52
30-0083.....	20
30-0084.....	52
30-0085.....	20
30-0086.....	52
30-0088.....	11
30-0089.....	14
30-0093.....	14
30-0094.....	14
30-0110.....	65

30-0111.....	65
30-0112.....	65
30-0113.....	65
30-0114.....	65
30-0115.....	65
30-0116.....	65
30-0117.....	65
30-0118.....	65
30-0119.....	65
30-0120.....	65
30-0121.....	65
30-0122.....	65
30-0128.....	12
30-0132.....	20
30-0133.....	52
30-0137.....	12

30-0142.....	63
30-0143.....	63
30-0144.....	63
30-00146.....	20
30-0160.....	12
30-0162.....	65
30-0165.....	14
30-0168.....	14
30-0171.....	14
30-5550.....	63

50

50-0004.....	62
50-0006.....	62
50-0008.....	62
50-0015.....	62
50-0032.....	62
50-0033.....	62
50-0034.....	62
50-0039.....	61
50-1005.....	62

60

60-0001.....	55
60-0002.....	54
60-0004.....	56
60-0005.....	56
60-0006.....	56
60-0009.....	57
60-0010.....	57
60-0011.....	57
60-0014.....	57
60-0022.....	55
60-0023.....	55
60-0025.....	57
60-0027.....	57
60-0029.....	57
60-0042.....	58
60-0043.....	58

60-0044.....	58
60-0045.....	58
60-0046.....	58
60-0047.....	56
60-0048.....	56
60-0049.....	56
60-0050.....	56
60-0051.....	54
60-0052.....	54
60-0061.....	56
60-0063.....	55
60-0065.....	56
66-0019.....	59

70

70-0001.....	63
70-0002.....	63
70-0003.....	63
70-0004.....	63
70-0005.....	62
70-0009.....	63
70-0010.....	63
70-0011.....	63
70-0012.....	63
70-0014.....	76
70-0015.....	69
70-0016.....	69
70-0017.....	8
70-0018.....	8
70-0019.....	8
70-0020.....	8
70-0021.....	8
70-0022.....	8
70-0024.....	9
70-0025.....	10
70-0026.....	10
70-0027.....	10
70-0028.....	10
70-0033.....	10
70-0034.....	10

70-0035.....	10	80-0096.....	35	80-0301.....	39	80-0467.....	22
70-0036.....	10	80-0097.....	44	80-0303.....	39	80-0468.....	23
70-0041.....	8	80-0117.....	35	80-0308.....	37	80-0469.....	23
70-0042.....	8	80-0118.....	35	80-0338.....	42	80-0470.....	23
70-0043.....	8	80-0121.....	32	80-0339.....	42	80-0471.....	24
70-0044.....	8	80-0122.....	32	80-0340.....	42	80-0472.....	24
70-0045.....	8	80-0133.....	32	80-0355.....	11	80-0473.....	24
70-0046.....	8	80-0134.....	32	80-0358.....	31	80-0474.....	24
70-0048.....	10	80-0135.....	32	80-0360.....	31	80-0475.....	24
70-0049.....	10	80-0136.....	32	80-0361.....	65	80-0476.....	24
70-0050.....	10	80-0139.....	73	80-0363.....	7	80-0477.....	9
70-0051.....	10	80-0140.....	73	80-0364.....	7	80-0479.....	12
70-0070.....	6	80-0141.....	73	80-0365.....	7	80-0480.....	12
70-0071.....	6	80-0142.....	73	80-0366.....	7	80-0481.....	12
70-0072.....	6	80-0151.....	36	80-0370.....	42	80-0482.....	12
70-0073.....	6	80-0167.....	31	80-0371.....	42	80-0483.....	12
70-0081.....	10	80-0168.....	31	80-0372.....	43	80-0484.....	12
70-0082.....	10	80-0173.....	30	80-0373.....	43	80-0485.....	12
70-0083.....	10	80-0175.....	30	80-0374.....	43	80-0486.....	12
70-0084.....	10	80-0176.....	30	80-0375.....	43	80-0487.....	12
70-0092.....	67	80-0181.....	37	80-0383.....	7	80-0494.....	26
70-0093.....	67	80-0183.....	37	80-0384.....	7	80-0495.....	26
70-0094.....	8	80-0184.....	37	80-0385.....	7	80-0496.....	26
70-0095.....	8	80-0201.....	36	80-0386.....	7	80-0497.....	26
70-0097.....	8	80-0202.....	36	80-0410.....	6	80-0498.....	26
70-0110.....	69	80-0203.....	36	80-0411.....	6	80-0499.....	26
70-0111.....	69	80-0204.....	36	80-0415.....	65	80-0500.....	26
70-0112.....	69	80-0205.....	31	80-0426.....	30	80-0501.....	26
70-0113.....	69	80-0207.....	31	80-0427.....	30	80-0502.....	26
70-0114.....	69	80-0208.....	31	80-0439.....	39	80-0503.....	26
70-0115.....	69	80-0213.....	38	80-0440.....	39	80-0504.....	26
70-0116.....	69	80-0223.....	31	80-0441.....	39	80-0505.....	26
70-0117.....	69	80-0244.....	37	80-0442.....	39	80-0509.....	29
70-0118.....	69	80-0247.....	30	80-0443.....	39	80-0510.....	29
		80-0254.....	31	80-0444.....	39	80-0511.....	29
		80-0256.....	42	80-0445.....	39	80-0512.....	29
		80-0257.....	42	80-0446.....	39	80-0513.....	29
		80-0258.....	42	80-0447.....	7	80-0514.....	29
		80-0259.....	42	80-0448.....	7	80-0515.....	29
		80-0260.....	42	80-0452.....	22	80-0516.....	29
		80-0261.....	42	80-0453.....	22	80-0517.....	29
		80-0269.....	31	80-0454.....	22	80-0518.....	29
		80-0270.....	31	80-0457.....	45	80-0519.....	29
		80-0276.....	31	80-0458.....	45	80-0520.....	29
		80-0289.....	30	80-0460.....	23	80-0522.....	22
		80-0298.....	38	80-0462.....	23	80-0523.....	22
		80-0299.....	38	80-0465.....	22	80-0524.....	23
		80-0300.....	39	80-0466.....	22	80-0525.....	23

80

80-0027.....	21
80-0059.....	44
80-0060.....	44
80-0067.....	42
80-0068.....	42
80-0069.....	42
80-0070.....	35
80-0071.....	35
80-0072.....	35
80-0089.....	44
80-0094.....	42

80-0526.....	24	85-0012.....	33	90-0077.....	72	ZZ-0247.....	63
80-0527.....	24	85-0020.....	44	90-0079.....	74	ZZ-0248.....	63
80-0542.....	46	85-0048.....	33	90-0081.....	74	ZZ-0249.....	63
80-0543.....	60	85-0049.....	33	90-0103.....	71	ZZ-0250.....	77
80-0544.....	60	85-0050.....	33	90-0104.....	71	ZZ-0258.....	68
80-0545.....	9	85-0052.....	33	90-0205.....	70	ZZ-0259.....	68
80-0547.....	9	85-0053.....	33	90-0206.....	70	ZZ-0262.....	64
80-0556.....	9	85-0054.....	33	90-0207.....	70	ZZ-0298.....	18
80-0557.....	9	85-0163.....	41	90-0208.....	70	ZZ-0426.....	53
80-0558.....	9	85-0177.....	28	90-0303.....	71	ZZ-0467.....	53
80-0560.....	12	85-0178.....	28	90-0305.....	71	ZZ-0509.....	9
80-0561.....	12	85-0179.....	28	90-0400.....	71	ZZ-0510.....	9
80-0562.....	12	85-0180.....	28	90-0401.....	71	ZZ-0537.....	9
80-0564.....	45	85-0189.....	28	91-0000.....	66	ZZ-0539.....	76
80-0565.....	45	85-0190.....	28	91-0001.....	67	ZZ-0542.....	76
80-0573.....	23	85-0191.....	28	91-0013.....	67	ZZ-0543.....	75
80-0574.....	23	85-0192.....	28	91-0016.....	66	ZZ-0545.....	75
80-0575.....	23	85-0404.....	40	91-0017.....	67	ZZ-0546.....	75
80-0576.....	23	85-0410.....	40	91-0018.....	66	ZZ-0547.....	75
80-0577.....	25	85-0411.....	40	91-0019.....	67	ZZ-0548.....	75
80-0578.....	25	85-0417.....	41	93-0017.....	49	ZZ-0549.....	75
80-0579.....	25	85-0493.....	52	93-0019.....	49	ZZ-0550.....	75
80-0580.....	25	85-0518.....	52			ZZ-0551.....	75
80-0581.....	24	85-0519.....	50			ZZ-0553.....	75
80-0582.....	24	85-0520.....	51			ZZ-0554.....	75
80-0583.....	24	85-0521.....	51			ZZ-0617.....	53
80-0584.....	24	85-0550.....	50			ZZ-0630.....	61
80-0859.....	46	85-0715.....	34			ZZ-0631.....	61
80-0861.....	25	85-0716.....	34			ZZ-0632.....	61
80-0862.....	25	85-0717.....	34			ZZ-0633.....	61
80-0863.....	25	85-0740.....	51			ZZ-0634.....	61
80-0864.....	25	85-0834.....	31			ZZ-0635.....	61
80-0865.....	25	85-0889.....	38			ZZ-0651.....	9
80-0866.....	25					ZZ-0652.....	77
80-0867.....	25					ZZ-0668.....	61
80-0868.....	25					ZZ-0669.....	61
80-0872.....	26					ZZ-0670.....	61
80-0873.....	26					ZZ-0674.....	60
80-0874.....	26					ZZ-0675.....	61
80-0876.....	48					ZZ-0682.....	60
80-0877.....	48					ZZ-0755.....	61
80-0878.....	48					ZZ-0756.....	60
80-0879.....	48					ZZ-0758.....	61
80-0880.....	48					ZZ-0759.....	60
80-0881.....	48					ZZ-0794.....	60
						ZZ-0797.....	78
						ZZ-0798.....	78
						ZZ-0799.....	78
						ZZ-0800.....	78
						ZZ-0834.....	21

85

85-0008.....	33
85-0009.....	33

90

90-0001.....	74
90-0002.....	74
90-0003.....	74
90-0004.....	74
90-0005.....	75
90-0006.....	75
90-0007.....	75
90-0020.....	74
90-0044.....	74
90-0063.....	72
90-0064.....	72
90-0065.....	72
90-0066.....	10
90-0069.....	73

ZZ

ZZ-0034.....	15
ZZ-0051.....	21
ZZ-0052.....	21
ZZ-0056.....	18
ZZ-0057.....	64
ZZ-0063.....	62
ZZ-0064.....	62
ZZ-0081.....	68
ZZ-0092.....	78
ZZ-0093.....	78
ZZ-0094.....	78
ZZ-0107.....	78
ZZ-0108.....	78
ZZ-0109.....	10
ZZ-0121.....	78
ZZ-0176.....	78
ZZ-0194.....	53
ZZ-0212.....	68
ZZ-0215.....	63
ZZ-0216.....	63
ZZ-0217.....	63
ZZ-0218.....	63
ZZ-0222.....	68
ZZ-0235.....	78
ZZ-0246.....	63

PRODUCT index by NSN Number

4240-01-543-9618	Benchmade 7 Hook Safety Cutter	77	6515-01-540-7612	Medic Trauma Sheet, Black	42
4240-01-568-3219	Benchmade 8 Hook Safety Cutter	77	6515-01-541-0635	ARS Needle Decompression, 14 gauge	18
4240-01-598-4341	Triage 916 SBK	77	6515-01-562-9734	Tracheal Hook, NAR	15
6510-00-201-1755	Cravat Bandage	14	6515-01-568-0193	Cyclone Improved Pocket BVM	15
6510-01-503-2117	Compressed Gauze	14	6515-01-587-1700	Combat Casualty Card, 25 Pack	64
6510-01-532-8019	HyFin Occlusive	17	6515-01-591-0397	King LTS-D, Purple, Size 5	16
6510-01-541-8121	ETD, Abdominal/Stump	14	6515-01-591-1578	King LTS-D, Red, Size 4	16
6510-01-558-4108	ETD, 4 in. Standard	14	6515-01-598-1877	Fox Eye Shield, With Garter, Box of 50	63
6510-01-558-4114	ETD, 6 in. Standard	14	6515-01-600-0262	Bear Claw Gloves, Kit, 25 pr, Large	63
6510-01-562-3325	Combat Gauze, Standard	13	6515-01-600-6039*	Bear Claw Gloves, Box, 50 pr, Large	63
6510-01-591-7740	ChitoGauze	13	6515-01-600-6052*	Bear Claw Gloves, Box, Extra Large	63
6510-01-600-4705	Combat Gauze, XL	13	6515-01-600-6057*	Bear Claw Gloves, Box, 50 pr, Medium	63
6510-01-605-2234	S-Rolled Gauze	14	6515-01-600-6059*	Bear Claw Gloves, Box, 50 pr, Small	63
6510-01-614-0938	Tactical Compression Wrap	14	6515-01-607-1943	Needleless Saline Lock Kit, NAR	19
6510-01-619-8636	Celox Gauze, 10 ft	13	6515-01-612-7129	CAT, Rescue Orange	11
6510-01-622-1624*	Trauma Shears, 6.25 inch, Small	62	6515-01-616-5841	JETT	11
6510-01-623-9910	Celox Gauze, 5 ft	13	6515-01-619-5350	Bear Claw Gloves, Kit, 25 pr, Small	63
6510-01-624-0840	HyFin Vent Individual	17	6515-01-619-5367	Bear Claw Gloves, Kit, Extra Large	63
6510-01-642-6212	HyFin Vent Twin Pack	17	6515-01-619-5671	Sharps Shuttle, 6 pack	20
6510-01-654-6801	ETD, 8x10 inch	14	6515-01-619-5689	Bear Claw Gloves, Kit, Medium	63
6510-01-659-3942	Combat Ready, Multi-Pack	61	6515-01-624-7148	T2 Tag, 20 Pack	64
6510-01-659-7998	Combat Ready, Bravo	61	6515-01-630-5551	Polycarbonate Eye Shield, Box of 100	63
6510-01-659-7999	Combat Ready, Alpha	61	6515-01-632-0059	Black Talon Gloves, Kit, 25 pr, Large	63
6510-01-659-8001	Combat Ready, Charlie	61	6515-01-646-2636	Torq First Light	69
6510-01-659-8349	Combat Ready, Echo	61	6515-01-647-5456	Polycarbonate Eye Shield, Single	63
6510-01-659-8353	Combat Ready, Delta	61	6530-01-492-4079	Stingray Poleless Litter	56
6515-01-338-6602	AMBU RDIC Mark III Resuscitator	16	6530-01-492-4114	Tactical Extrication Device, OD Green	56
6515-01-346-9186	Tactical Traction Splint, Black	62	6530-01-504-9055	TALON Carrier, OD Green	57
6515-01-449-1016	Fox Eye Shield, No Garter, Box of 12	63	6530-01-523-3314	33 in. Litter Stands	57
6515-01-450-9790	Laryngoscope Set, Special Operations	16	6530-01-531-6808	TED. Tactical Extrication Device, Black	56
6515-01-458-6057	Field Corpsman Kit, Basic	64	6530-01-531-6854	TALON Assault Litter Carrier, Coyote	57
6515-01-458-6178	Field Corpsman Kit, Deluxe	64	6530-01-531-6907	TALON Assault Litter Carrier, Black	57
6515-01-494-1951	Sam Splint II	62	6530-01-549-1089	PaxLight LED Field Medical Light	58
6515-01-509-6866	SAM Pelvic Sling	62	6530-01-588-8955	PaxLight Deluxe Battery Base Kit	58
6515-01-515-0190	Black Talon Gloves, Box, Extra Large	63	6530-01-594-7444	Backrest, Talon or Raven Litter	57
6515-01-515-0195	Black Talon Gloves, Box, 50 pr, Small	63	6530-01-598-9595	Dragon Handle	74
6515-01-515-0197	Black Talon Gloves, Box, 50 pr, Large	63	6530-01-605-9432*	CAT Holder, OD Green	12
6515-01-515-0200	Black Talon Gloves, Box, 50 pr, Medium	63	6530-01-624-7626	Mark IV Field Hospital Bed	59
6515-01-518-8551	Sharps Shuttle, Single Pack	20	6530-01-632-7872	RAT Strap	74
6515-01-519-3612	Ambu Spur II BVM	15	6530-01-643-9134*	TALON Carrier, Digital Camo	57
6515-01-521-7976	CAT, Tactical Black	11	6532-01-627-7888	HRS Heat Reflective Shield	21
6515-01-523-3317	Raptor IV Securing Device, 6 Pack	20	6532-01-653-6054	Polar Skin Ice Sheets, Standard	61
6515-01-523-5403	Sam Splints, 12 inch	62	6545-01-532-4962	W.A.L.K. Digital Camo	44
6515-01-523-5440	Sam Splint, 15 inch	62	6545-01-572-7251	Medic Trauma Pack, Black	35
6515-01-529-1187	Nasopharyngeal Airway & Lube, 28F	15	6545-01-587-1199	W.A.L.K. Coyote	44
6515-01-532-8056	HPMK	21	6545-01-588-1862	Medic Trauma Pack, Coyote	35
6515-01-537-2611	BOA IV Constricting Band, Standard	19	6545-01-592-3261	T2 Tag, 1000 Pack	64
6515-01-537-4094	Saline Lock Kit, NAR	19	6545-01-594-5617*	Medic Trauma Pack, OD Green	35
6515-01-537-4161	Combat Casualty Reference Card	64	6545-01-594-9399*	Armadillo Case, Clear Lid	68
6515-01-538-9276	Trauma Shears, 7.25 in., Large	62	6545-01-610-5411*	Range Trauma Kit	38
6515-01-540-7206	Tactical Suction Device, NAR	15	6910-01-560-2972	CAT, Trainer Blue	11, 49
6515-01-540-7568	Tactical CricKit	15	8145-01-573-2533	Armadillo Case, Black Lid	68
6515-01-540-7579	Medic Trauma Sheet, Digital Camo	42	8145-01-617-9075*	Armadillo Hook Discs, Set of 250	68
6515-01-540-7585	Medic Trauma Sheet, OD Green	42			
6515-01-540-7603	Medic Trauma Sheet, Coyote	42			

*NATO NSN



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

☎ Teléfono gratuito

900 535 295

www.iesmedical.es

✉ email contacto

comercial@iesmedical.es

☎ Atención comercial

688 666 089

☎ Oficinas centrales

944 008 847

IESMEDICAL
Innovation for better care