

Our reference: REP179
Subject: Evone - airway adapter use

Eindhoven, 20 Apr 2020

To whom it may concern,

In the current situation, it may be required to use alternative disposables for the Evone system because of supply shortages. In this letter we provide guidance for selecting these alternatives. However, please be aware that when using non-approved parts, it remains the responsibility of the hospital and the involved medical professionals to ensure correct functioning of the complete Evone system.

In response to questions from our users regarding the use of alternative airway adapters with the Evone® system we like to state the following:

- The Evone Airway adapter, ref 6125, is the approved and verified adapter for use with the Evone system.
- We have understood that some of our customers use the Philips Respironics M2516A (Infant/Neonatal) Airway Adapter and obtain accurate measurements. As we have not verified this adapter, using this adapter remains the responsibility of the hospital.
- Other Philips Respironics Airway Adapters may work when used with the Evone system but it remains the responsibility of the hospital to ensure correct functioning. When using non-approved airway adapters, please take into account the following:
 - Use the zeroing function in the Start-up menu when switching adapters.
 - The Adult/Pediatric Adapters will add approximate 4cc of dead space.
 - Cross check measured CO2 values with e.g. blood gas analysis.
 - Single-Patient Disposable Adapters are expected to wear out faster due to different pressures at the Evone Cartridge side.

With kind regards,



Rienco Verschelling
COO
Ventinova Medical B.V.

Ventinova Medical B.V.

A Meerenakkerplein 7
5652 BJ Eindhoven
The Netherlands

T +31 (0)40 751 60 20
E info@ventinova.nl
I www.ventinova.nl

Chamber of Commerce 58307834
VAT NL8529.77.694.B01
IBAN NL71 RABO 0314885706
BIC RABONL2U